

An Overview of the Capillary Electrophoresis Process at the DOE Joint Genome Institute Production Genomics Facility: The Dual Operation of the AB 3730xl & GE MegaBACE 4500 DNA Sequence Analyzers

Christopher Daum, Lena Philip, Danielle Mihalkanin, Cailyn Spurrell, Don Miller, Susan M. Lucas, Alex Copeland, Damon Tighe, Eric Abbott, Marlon Arcaina, Nicholas Eattock, Melanie Lafrades, Albert Linkowski, Adrienne Loero-Pequignot, Andy Yuen, Mathew Zane



ABSTRACT

At the center of the Department of Energy's (DOE) Joint Genome Institute (JGI) Production Genomics Facility (PGF), lies a highly efficient and automated production line devoted to the generation of high-quality genomic DNA sequence. The JGI utilizes a dual platform of DNA sequence analyzers within its Production Capillaries group: seventy Applied Biosystems 3730xl and thirty-six GE Healthcare MegaBACE 4500 instruments. The Capillaries group is comprised of eleven employees that are responsible for operating and maintaining both platforms; the group is also involved in the daily monitoring of performance stats and the troubleshooting of DNA sample and instrument related issues. The operation of these high-throughput fluorescence-based DNA sequence analyzers will be assessed on the strengths and benefits of each platform, including instrument overview of operational parameters and mechanical/component specifications. In addition, instrument setups for production operation, operation schedules, loading, and maintenance strategies as well as the various sequencing strategies for each platform will be compared. Throughput numbers and sequencing quality results will be presented.

Introduction

The MegaBACE DNA Analysis platform has been an indispensable component of the JGI production sequencing line for several years; with the inception of this capillary based platform in 1999, significant advances were made over the industry standard slab gel instruments of the time – in areas of automation, workflow, and operation. The current MB4500 platform offers much improved sample throughput with the ability to process 384 sample lanes in a single run and improved sequence quality over its predecessors with updated components.

In 2002 the ABI3730xl was released - this highly automated DNA analyzer was poised to shift the paradigm of high-throughput sequencing as the JGI and genomics facilities world-wide adopted this new platform; a new industry standard had been set.

The JGI has chosen to utilize both the MB4500 - 36 instruments, and the ABI3730xl - 70 instruments. The highly automated ABI3730xl allows for 24hr-7day/week processing of sample plates, with only limited interaction by a small staff of Technicians who load/unload samples and change out reagents a few times per week. Alternatively, the 384-capillary array system of the MB4500 allows for high sample throughput on each sequencing run performed, however it requires manual interaction with a Technician to load each run.

The MB4500 uses an efficient high power solid-state laser, which matched to the scanning confocal optics system provides enhanced detection sensitivity for long reads; attaining on average 100bases more per lane (based upon JGI run parameters). Additionally, the MB capillary arrays and solid-state laser achieve much longer operational lifetimes when compared to the arrays and argon laser of the ABI. Although in JGI production, the MB instruments have an operational uptime efficiency of 97.9% compared to 98.8% for the ABI; downtime on both platforms is primarily due to instrument errors.

AB 3730xl DNA Analyzer



MegaBACE 4500



Instrument Error Response

Capillary operators are alerted of instrument errors in a variety of ways. When a MegaBACE error, it usually occurs while the operator is attempting to follow the command prompts while loading the machine. The MegaBACE software initializes a pop-up that notifies the user that an error has occurred, where the operator then has the option to respond to the error and correct it. When the error can not be corrected the service engineer is notified.

ABI 3730xl errors are found in a different way. Because the ABIs operate autonomously, errors usually occur when the operator is not close to the instrument. Therefore each instrument is connected to a system monitoring program called BioMonitor and when an error occurs notifications are sent out via email and to pagers that the operators carry. When BioMonitor is not working correctly, and notifications are not sent, then down instruments are located during routine walk-throughs. Operators are trained to fix many minor errors and some major errors themselves; however, repeating errors or serious system errors and failures (i.e.: laser, camera, electronics) are reported to the service engineer.




Donald Miller: Donald Miller is the service engineer for the MegaBACE 4500s. He has been a JGI employee since Jan. 2002, and has completed training offered through Amersham (GE Healthcare) in order to be certified to work on the MegaBACE platform. He performs all of the scheduled preventative maintenance and troubleshoots all instrument errors that occur. When errors occur during hours that Don is not available, the instruments are stored until he can look at them.




Service Engineers

Dennis Schuster: Dennis Schuster is the service engineer for the ABI3730xl. He has been an Applied Biosystems employee for 11 years, and has worked on-site at JGI 5 days / week since Sept. 2005. He performs all of the scheduled ABI preventative maintenance as well as troubleshoots 99% of the instrument errors that occur. ABI operators fix minor errors during hours when Dennis is not on-site.


JGI Production Process Summary



JGI Throughput Comparison of 36 MB4500s vs. 70 ABI3730s



Performance Comparison of ABI 3730xl & MegaBACE 4500 Platforms



JGI Quick Reference – Operational Parameters & Specifications

Operational Parameter	ABI 3730xl Specification	MegaBACE 4500 Specification
Run Parameters:	8.9kV for 70min	8.9kV for 120min
Sequencing Parameters:	1.8kV for 15min	1.8kV for 15min
Sequencing Temperature:	40°C	40°C
Sequencing Method:	PGM™ polymerase	3'2'2' LPA
Capillary Array Type:	70 cm - 70 capillaries 150 meters w.d. x 70 meters i.d.	70 cm - 384 capillaries 200 meters w.d. x 70 meters i.d.
Sample Lane Throughput (based upon JGI Run parameters):	1,440 samples/day	2,304 samples/day
Laser power type:	27Watt argon-ion	100Watt solid state
Laser accessories (readings):	External RFLAC system with argon-ion gas leading to four 400 nm lines	None - use external Power Supply for module in RFLAC system output
Optical system - Excitation & Emission:	Scanning confocal optics for the capillary laser lenses & detectors	Scanning confocal optics for the capillary laser lenses & detectors
Instrument accessories:	SPS for loading during power outages 10" Monitor, Barcode Scanner	SPS for loading during power outages 10" Monitor, Barcode Scanner Step & Scan power supply
Reagents/Materials Handling:	Automated reagent delivery pump & on-board reagent for 2 day maintenance operation. Reagent and waste sample & sample plate collection (16)	Reagent and waste sample & sample plate collection (16) Reagent and waste sample & sample plate collection (16)
Sequencing Software:	Sequencing Analysis v2.2 with 3'2'2' LPA	Sequencing Analysis v2.2 with 3'2'2' LPA
Sequencing Software:	Sequencing Analysis v2.2 with 3'2'2' LPA	Sequencing Analysis v2.2 with 3'2'2' LPA



Capillary Operators Left to Right: Andy Yuen, Nicholas Eattock, Donald Miller, Mathew Zane, Adrienne Loero-Pequignot, Christopher Daum, Marlon Arcaina, Melanie Lafrades, Lena Philip, Albert Linkowski, Don Miller, Cailyn Spurrell, Dennis Schuster, Andrew Yuen

Group History

The Production Capillaries group was established in 1999 with the operation of 84 MegaBACE 1000 instruments. The group acquired its first 5 ABI3730xl instruments in May 2002, and by 2004 70 instruments total had been brought online. The MegaBACE 1000s were replaced with 21 MegaBACE 4000s in early 2002 and an additional 15 were brought online in 2004. In 2003 two MegaBACE 4000s were upgraded to the developmental MegaBACE 4500 platform. After extensive testing and development, all of the MegaBACE 4000s were upgraded to MegaBACE 4500s; a process that was completed in early 2005.

The capillaries group was absorbed into the QC group in December 2002 and capillary operators functioned as both instrument operators and QC members. Operators not only ran the instruments, but ran electrophoresis gels, engaged in the troubleshooting of statistics and production line problems, and ran experiments. In January 2006 after some facility-wide restructuring, the capillaries group once again became an autonomous group within the production line under the supervision of Chris Daum. While the capillaries group is autonomous, members still partially function as a QC group and continue to run gels and troubleshoot. Operators also continue to run experiments, but they are tailored more towards capillary improvement as opposed to general QC experiments.

Conclusion

The JGI has chosen to utilize both the MegaBACE 4500 and the ABI3730xl DNA sequencer analyzer instruments, and has found the running of these competing platforms to be beneficial to the production sequencing process. The highly automated ABI3730xl allows for 24hr-7day/week processing of sample plates, with only limited interaction. Alternatively, the 384-capillary array system of the MB4500 allows for high sample throughput during each sequencing run performed by a Technician, and its solid-state laser mated with the scanning confocal optics system provides enhanced detection sensitivity for longer reads. The MB platform also offers slightly reduced per lane cost when comparing the associated operational and sample preparation costs for each platform at the JGI; this reduced operational cost is primarily due to the long lifetimes achieved by the MB capillary arrays. Moreover, these two DNA sequencing platforms complement each other in the sequence they produce - thereby providing for better coverage when sequence derived from both platforms are used in genomic project assemblies at the JGI.

Weekly ABI Tasks	MB4500 Run Schedule
ABI Tasks: M T W T F S S Load Plates: M T W T F S S Unload Plates: M T W T F S S Array QC: M T W T F S S Array Cleanup: M T W T F S S	MB4500 Run Schedule: M T W T F S S Run 1: M T W T F S S Run 2: M T W T F S S Run 3: M T W T F S S Run 4: M T W T F S S Run 5: M T W T F S S Run 6: M T W T F S S Run 7: M T W T F S S Run 8: M T W T F S S Run 9: M T W T F S S Run 10: M T W T F S S Run 11: M T W T F S S Run 12: M T W T F S S Run 13: M T W T F S S Run 14: M T W T F S S Run 15: M T W T F S S Run 16: M T W T F S S Run 17: M T W T F S S Run 18: M T W T F S S Run 19: M T W T F S S Run 20: M T W T F S S Run 21: M T W T F S S Run 22: M T W T F S S Run 23: M T W T F S S Run 24: M T W T F S S Run 25: M T W T F S S Run 26: M T W T F S S Run 27: M T W T F S S Run 28: M T W T F S S Run 29: M T W T F S S Run 30: M T W T F S S Run 31: M T W T F S S Run 32: M T W T F S S Run 33: M T W T F S S Run 34: M T W T F S S Run 35: M T W T F S S Run 36: M T W T F S S Run 37: M T W T F S S Run 38: M T W T F S S Run 39: M T W T F S S Run 40: M T W T F S S Run 41: M T W T F S S Run 42: M T W T F S S Run 43: M T W T F S S Run 44: M T W T F S S Run 45: M T W T F S S Run 46: M T W T F S S Run 47: M T W T F S S Run 48: M T W T F S S Run 49: M T W T F S S Run 50: M T W T F S S Run 51: M T W T F S S Run 52: M T W T F S S Run 53: M T W T F S S Run 54: M T W T F S S Run 55: M T W T F S S Run 56: M T W T F S S Run 57: M T W T F S S Run 58: M T W T F S S Run 59: M T W T F S S Run 60: M T W T F S S Run 61: M T W T F S S Run 62: M T W T F S S Run 63: M T W T F S S Run 64: M T W T F S S Run 65: M T W T F S S Run 66: M T W T F S S Run 67: M T W T F S S Run 68: M T W T F S S Run 69: M T W T F S S Run 70: M T W T F S S Run 71: M T W T F S S Run 72: M T W T F S S Run 73: M T W T F S S Run 74: M T W T F S S Run 75: M T W T F S S Run 76: M T W T F S S Run 77: M T W T F S S Run 78: M T W T F S S Run 79: M T W T F S S Run 80: M T W T F S S Run 81: M T W T F S S Run 82: M T W T F S S Run 83: M T W T F S S Run 84: M T W T F S S Run 85: M T W T F S S Run 86: M T W T F S S Run 87: M T W T F S S Run 88: M T W T F S S Run 89: M T W T F S S Run 90: M T W T F S S Run 91: M T W T F S S Run 92: M T W T F S S Run 93: M T W T F S S Run 94: M T W T F S S Run 95: M T W T F S S Run 96: M T W T F S S Run 97: M T W T F S S Run 98: M T W T F S S Run 99: M T W T F S S Run 100: M T W T F S S

Sample MegaBACE run schedule: Four runs / day split between 2 primary operators and 4 secondary operators. Each run is comprised of 5 sets and each set has 7 instruments; except for the "off-set", which has 8 instruments. This allows for a total of 144 plates to be processed/day. In between runs, operators prepare reagents for the following shift.

Five walk-throughs are performed each day (every 4 hours), which consist of checking the status of the ABI 3730xls and unloading plates that have already been processed. Group members also update the arriving shift of problems that may have occurred during the day.