



## Introduction

The ability to sequence genomes far outpaces our ability to systematically determine gene function. Subsequently, systems-level analyses of less studied bacteria are limited by the presence of numerous uncharacterized genes and an over reliance on annotations from well studied bacteria such as *E. coli*. To meet this challenge, we are developing a flexible mutagenesis technique and applying it to the environmental bacteria *Shewanella oneidensis* MR1, *Shewanella loihica* PV4, and *Desulfovibrio desulfuricans* G20. The hallmarks of our approach are the sequencing and archiving of thousands of transposon mutants and the use of molecular "barcodes" (tags) for the parallel phenotypic analysis of defined mutant pools. The successful completion of this project will enable the quantitative phenotypic analysis of thousands of mutants across a wide range of conditions. In addition, our genetic resources provide a framework for the systematic genetic interrogation of individual pathways.


## Cloning of Tag Modules



## Transfer of Tag Modules via Gateway Reaction




## Our Focus Organisms




We have converted pRL27 (Tn5), PHIMAR-RB1 (mariner), and an *in vitro* Tn5 system (epicentre) into Gateway destination vectors. The tagged transposons are functional in the indicated organisms.

## High-throughput Transposon Insertion Mapping




- After random mutagenesis with the collection of tagged transposons, we determine the insertion location of the transposon in the genome and the tag identify using a 2-step degenerate PCR approach. Sequence analysis of the second round PCR product is used to infer both the tag identify and the genome insertion site.
- We anticipate sequencing ~30,000 colonies per bacterial genome in order to get adequate coverage. Mutants with unique insertion sites will be selected from these libraries and combined into pools of ~5000 mutants (each with a unique tag module).
- To facilitate this process, we have optimized the use of a liquid handling robot (Biomek FX) to process (PCR and sequencing) 1536 mutants per day.

## Parallel Analysis of Mutant Pools



## Application to *Shewanella oneidensis* MR1



## Future Directions

- Generate ~2000 additional sequence-defined tag modules.
- Go into production phase and rapidly generate thousands of sequence-defined transposon mutants in multiple bacterial genomes.
- Explore the potential for double deletion construction for global genetic interaction studies.
- Automate and miniaturize mutant pool experiments using procedures developed for the yeast deletion collection. Such an effort will enable to profile comprehensive mutant libraries across hundreds of diverse conditions.

## ACKNOWLEDGEMENT

ESPP is part of the Virtual Institute for Microbial Stress and Survival supported by the U. S. Department of Energy, Office of Science, Office of Biological and Environmental Research, Genomics Program-GTL through contract DE-AC02-05CH11231 between Lawrence Berkeley National Laboratory and the U. S. Department of Energy.