



Lawrence Berkeley Laboratory

UNIVERSITY OF CALIFORNIA

CHEMICAL BIODYNAMICS DIVISION

I. CONFORMATIONAL ANALYSIS OF CYCLIC PEPTIDES
II. MULTINUCLEAR NMR SPECTROMETER FOR THE STUDY OF BIOLOGICAL SYSTEMS

Willy Chao-Wei Shih
(Ph.D. thesis)

November 1979


RECEIVED
LAWRENCE
BERKELEY LABORATORY

JAN 14 1980

LIBRARY AND
DOCUMENTS SECTION

TWO-WEEK LOAN COPY

*This is a Library Circulating Copy
which may be borrowed for two weeks.
For a personal retention copy, call
Tech. Info. Division, Ext. 6782.*



DISCLAIMER

This document was prepared as an account of work sponsored by the United States Government. While this document is believed to contain correct information, neither the United States Government nor any agency thereof, nor the Regents of the University of California, nor any of their employees, makes any warranty, express or implied, or assumes any legal responsibility for the accuracy, completeness, or usefulness of any information, apparatus, product, or process disclosed, or represents that its use would not infringe privately owned rights. Reference herein to any specific commercial product, process, or service by its trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or favoring by the United States Government or any agency thereof, or the Regents of the University of California. The views and opinions of authors expressed herein do not necessarily state or reflect those of the United States Government or any agency thereof or the Regents of the University of California.

I. Conformational Analysis of Cyclic Peptides

II. Multinuclear NMR Spectrometer for the Study of Biological Systems

The United States Department of Energy has
the right to use this thesis for any purpose
whatsoever including the right to reproduce
all or any part thereof.

(c) Copyright 1979

by

Willy Chao-Wei Shih

- I. Conformational Analysis of Cyclic Peptides
- II. Multinuclear NMR Spectrometer for the Study of Biological Systems

Willy Chao-Wei Shih
Lawrence Berkeley Laboratory, University of California
Berkeley, California 94720


ABSTRACT

A systematic approach to the employment of ^1H NMR data to the analysis of the solution conformations of small polypeptides is outlined. Two-dimensional homonuclear J-spectra along with the corresponding 45° projections and contour maps simplify the task of homonuclear decoupling experiments, and the assignments of lines is made very straightforward. The vicinal couplings that are likely to vary as the conformation changes are then examined. Spectrum simulations coupled with x-ray data for model systems allow the formulation of Karplus relationships for the vicinal couplings, and solution conformations may then be inferred.

This approach to conformational analysis was applied to an examination of the metal affinities of *cyclo*-[3-(4- β -aminoethyl)-phenyloxypropanoyl-L-prolyl] and *cyclo*-[3-(4- β -N-methylaminoethyl)-

phenyloxypropanoyl-L-prolyl]. The NMR data revealed that the metal affinity of the former is probably a result of the orientation of both carbonyl groups towards the same face of the ring; the sharply lower affinity of the N-methyl cyclopeptide for divalent cations is likely a result of the carbonyl facing opposite sides of the ring.

Part II describes the design approach, major construction details, and the test and characterization procedures associated with the construction of a broadband multinuclear NMR spectrometer system. The design emphasis was on incorporating features that would be useful for biological experiments. Summaries of original designs are included with schematic documentation.

A handwritten signature in cursive script, appearing to read "William P. Klein", is written in dark ink on the right side of the page.

ACKNOWLEDGMENTS

It is nearly impossible to acknowledge the assistance of all the people who contributed to the completion of this work. I would particularly like to thank Prof. Melvin Calvin for his support, and for providing a unique educational opportunity. I am especially grateful to Dr. Melvin P. Klein for his wisdom, counsel, and guidance. I shall always respect him as a scientist, an educator, and as a friend.

Prof. Henry Rapoport and Dr. J. Clark Lagarias made major contributions to the cyclopeptide work, and Prof. F. Bordner kindly supplied the x-ray data on the N-methyl cyclopeptide. John Griffin, Ferenc Kovac, Dick O'Brien, Cees DeGroot, and Mike Press provided technical assistance; Sheldon Wong and Ken Wiley assisted in the VAX data transmission and contour plotting work.

Special thanks are due Gary Lee, for filling the liquid nitrogen and liquid helium on the magnets countless times, and for helping on innumerable other tasks. Karen Cornwell and Maria Fassett provided expert assistance in the nitrogenase isolation work; Ben Gordon and Wally Erwin assisted with GCs, plumbing, and many other problems. I thank Howard Wood for teaching me how to use machine tools, and for helping me salvage many projects.

Special thanks are due members of the Electronics Research and Development Group at the Lawrence Berkeley Laboratory, headed by Dr. Branko Leskovar. C.C. Lo and Eric Young provided much valuable advice (and the loan of much valuable test equipment), and I deeply appreciate their help. Thanks are also due Dr. J.D. Swalen for discussions on the iterative fit procedures, and for providing program tapes.

I am especially grateful to Dr. David Dalrymple of Nicolet Technology Corp. Even though we had a customer relationship, he displayed extraordinary patience during many hours on the phone.

I thank my friends around LCB -- Beth Klingel and Lois Soule in particular -- for helping me to maintain a sense of humor, and I particularly thank Craig Hodges and Mary McLean for many valuable consultations.

A special thanks is reserved for the people without whom this work would not have been completed. Gary Smith, in addition to being a good friend, taught me an approach to doing design work and troubleshooting. Without his guidance, the spectrometer system would not approach its current performance level.

I thank Janice DelMar for her friendship, understanding, and counsel. Her wisdom and guidance during many difficult times will always be remembered and appreciated. I thank Barbara Plowman for her friendship and understanding, and for her patience during the preparation of this work. Finally I thank my parents, who have stood by me from the start, hoping to see this day.

This work was supported, in part, by the Biomedical and Environmental Research Division of the U.S. Department of Energy under contract W-7405-eng-48.

TABLE OF CONTENTS

Abstract	1
Acknowledgments	i
Table of Contents	iii
Figure List	vi
List of Tables	xvi
I. Conformational Analysis of Cyclic Peptides	1
1.1 Introduction	3
1.2 Two-Dimensional J-Resolved Spectroscopy	7
1.3 Cyclopeptide Alkaloids	15
A. Occurrence, Structure, Model Systems	15
B. Metal Binding Properties	19
1.4 NMR Spectral Data for Cyclopeptide Alkaloid Model Compounds	27
A. Assignment of Lines	27
B. Comparison of Assignments	34
1.5 Analysis of NMR Data	73
A. Summary of Approach	73
B. C ₈ - C ₉ Four-Spin System	74
C. C ₁ - C ₂ Four-Spin System	76
D. Proline Seven-Spin Systems	77
E. Aromatic Four-Spin Systems	77
F. Analysis of Spectral Data of <i>Cyclo</i> -[3-(4- β - aminoethyl)phenoxy-4-methylpentanoyl- L-prolyl]	78
G. Metal Affinity of N-Methyl Cyclopeptide	78
1.6 Summary	127

A. Cyclopeptide Alkaloid Model Systems	
Conformations	127
B. Methodology	127
1.8 References	131
II. Multinuclear NMR Spectrometer for the Study of	
Biological Systems	133
2.1 Introduction	135
2.2 System Design Considerations	137
2.3 Major Subsystems	145
2.3.1 Magnet	145
2.3.2 Observe Channel Transmitter	145
2.3.3 Observe Channel Receiver	150
2.3.4 Decouple Channel Transmitter	155
2.3.5 Lock Channel	160
2.3.6 Probes	165
2.3.7 Data Processing and Control	186
2.4 Test and Characterization Procedures	192
2.4.1 Complex Impedance, Reflection Coefficients,	
and VSWR	192
2.4.2 Measurement Techniques	193
2.4.3 Evaluation of System Noise Performance	199
2.4.4 Noise Sources	200
2.4.5 Noise Figures	201
2.4.6 Noise Figure Measurement Procedure	204
2.5 System Components	207
2.5.1 16X892 System Control Logic	208
2.5.2 16X894 Lock Channel Control Unit	222

2.5.3	16X902 Remote Controlled Low-Pass Filter . . .	236
2.5.4	16X903 Low Frequency Rf Preamplifier	243
2.5.5	16X904 Nicolet 1180 I/O Bus Intercept	247
2.5.6	16X919 Lock Channel Single Sideband Generator	253
2.5.7	16X935 63.42 KGauss Magnet Room Temperature Shim Control	261
2.5.8	16X950 Nicolet 1180 Driven Probe Temperature Control	265
2.5.9	16X955 270 MHz Preamplifier	280
2.5.10	16X956 Nicolet 1180/CalComp 565 Incremental Plotter Interface	283
2.5.11	16X964 41.45 MHz Receiver - T/R Switch . . .	294
2.5.12	16X965 27.36 MHz / 270 MHz Dual Frequency Receiver - T/R Switch	298
2.5.13	16X966 67.89 MHz / 109.30 MHz Dual Frequency Receiver - T/R Switch	302
2.5.14	Diablo 31 Disk Memory Subsystem -- Power Supply Interconnect	306
2.5.15	Rf Power Amplifiers -- Power Supply Interconnect	306
2.5.16	16X973 Rf Switch	306
2.5.17	16X974 Preamplifier Power Supply	307
2.5.18	16X978 41.45 MHz Amplifier	307
2.5.19	16X980 Rf Switch	307
2.5.20	Data Processor Power Fail Interlock	317

2.6 References

Appendix 1	Nicolet 1180 /VAX 11/780 Data Transmission Protocol	321
Appendix 2	Complex Impedance Calculation	351
Appendix 3	Supplementary Schematic Documentation	353

FIGURE LIST

Figure 1.1	90°-τ-180°-τ Spin-Echo Experiment	10
Figure 1.2	Phase-Time Diagram for a Coupled Two-spin System, IS	13
Figure 1.3	J-Resolved Two-dimensional Experiment	14
Figure 1.4	Metal Binding Behavior of Ceanothine-B as Shown by Circular Dichroism	22
Figure 1.5	Circular Dichroism Spectra of Cyclopeptide Model Systems, 250 - 300 nm Region	23
Figure 1.6	Circular Dichroism Spectra of Cyclopeptide Model Systems, 215 - 250 nm Region	24
Figure 1.7	Metal Binding Behavior of <i>Cyclo</i> -[3-(4-β- aminoethyl)phenyloxypropanoyl-L-prolyl] as shown by Circular Dichroism	25
Figure 1.8	Stereo X-ray Crystal Structure of <i>Cyclo</i> - [3-(4-β-N-methylaminoethyl)phenyloxy- propanoyl-L-prolyl].	26
Figure 1.9	270 MHz ¹ H NMR Spectrum of <i>Cyclo</i> -[3-(4-β-N- methylaminoethyl)phenyloxypropanoyl- L-prolyl], 6a , with tentative line assignments	29
Figure 1.10	270 MHz ¹ H NMR Spectrum of <i>Cyclo</i> -[3-(4-β- aminoethyl)phenyloxypropanoyl-L-prolyl], 6b , with tentative line assignments	30
Figure 1.11	270 MHz ¹ H Two-dimensional Homonuclear J-spectrum of <i>Cyclo</i> -[3-(4-β-N-methyl- aminoethyl)phenyloxypropanoyl-L-prolyl].	31
Figure 1.12	270 MHz ¹ H Two-dimensional Homonuclear J-spectrum of <i>Cyclo</i> -[3-(4-β-aminoethyl)- phenyloxypropanoyl-L-prolyl]	32
Figure 1.13	270 MHz ¹ H NMR Spectrum of <i>Cyclo</i> -[3-(4-β- aminoethyl)phenyloxypropanoyl-L-prolyl] in the 2.60 - 3.10 δ Spectral Region	38
Figure 1.14	Contour Plot of the 2.62 - 3.08 δ Spectral Region of the Two-dimensional Homo- nuclear J-Spectrum of the N-H Cyclo- peptide	39

Figure 1.15	Contour Plot of the 7.34 - 6.88 δ Spectral Region of the Two-dimen- sional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	40
Figure 1.16	Contour Plot of the 6.98 - 6.52 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	41
Figure 1.17	Contour Plot of the 5.00 - 4.54 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	42
Figure 1.18	Contour Plot of the 4.46 - 4.00 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	43
Figure 1.19	Contour Plot of the 3.71 - 3.25 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	45
Figure 1.20	Contour Plot of the 3.14 - 2.68 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	46
Figure 1.21	Contour Plot of the 2.89 - 2.43 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	47
Figure 1.22	Contour Plot of the 2.46 - 2.00 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	48
Figure 1.23	Contour Plot of the 2.06 - 1.60 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide	49
Figure 1.24	Contour Plot of the 7.38 - 6.92 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	50

Figure 1.25	Contour Plot of the 7.06 - 6.60 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	51
Figure 1.26	Contour Plot of the 4.87 - 4.41 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	53
Figure 1.27	Contour Plot of the 4.50 - 4.04 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	54
Figure 1.28	Contour Plot of the 4.03 - 3.57 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	55
Figure 1.29	Contour Plot of the 3.62 - 3.16 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	56
Figure 1.30	Contour Plot of the 2.50 - 2.05 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	57
Figure 1.31	Contour Plot of the 2.28 - 1.82 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	58
Figure 1.32	Contour Plot of the 1.75 - 1.29 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide	59
Figure 1.33	270 MHz ^1H NMR Spectrum of <i>Cyclo</i> -[3-(4- β -N- methylaminoethyl)phenoxypropanoyl- L-prolyl] -- 0 $^\circ$ Projection Sum of the 2-D Homonuclear J-Spectrum	61
Figure 1.34	270 MHz ^1H NMR Spectrum of <i>Cyclo</i> -[3-(4- β -N- methylaminoethyl)phenoxypropanoyl- L-prolyl] -- 45 $^\circ$ Projection Sum of the 2-D Homonuclear J-Spectrum	62

Figure 1.35	270 MHz ^1H NMR Spectrum of <i>Cyclo</i> -[3-(4- β -aminoethyl)phenoxypropanoyl-L-prolyl] -- 0° Projection Sum of the 2-D Homonuclear J-Spectrum	63
Figure 1.36	270 MHz ^1H NMR Spectrum of <i>Cyclo</i> -[3-(4- β -aminoethyl)phenoxypropanoyl-L-prolyl] -- 45° Projection Sum of the 2-D Homonuclear J-Spectrum	64
Figure 1.37	Approximate Chemical Shift and Coupling Constant Pattern of the $\text{C}_8\text{-C}_8'\text{-C}_9\text{-C}_9'$ Four-Spin Pattern, N-Methyl Cyclopeptide .	82
Figure 1.38	Calculated Chemical Shift and Coupling Constant Pattern of the $\text{C}_8\text{-C}_8'\text{-C}_9\text{-C}_9'$ Four-Spin Pattern, N-Methyl Cyclopeptide .	85
Figure 1.39	Calculated Spectrum for the $\text{C}_8\text{-C}_8'\text{-C}_9\text{-C}_9'$ Four-Spin Pattern, N-Methyl Cyclopeptide .	85
Figure 1.40	Approximate Chemical Shift and Coupling Constant Pattern of the $\text{C}_8\text{-C}_8'\text{-C}_9\text{-C}_9'$ Four-Spin Pattern, N-H Cyclopeptide . . .	88
Figure 1.41	Calculated Chemical Shift and Coupling Constant Pattern of the $\text{C}_8\text{-C}_8'\text{-C}_9\text{-C}_9'$ Four-Spin Pattern, N-H Cyclopeptide . . .	91
Figure 1.42	Calculated Spectrum for the $\text{C}_8\text{-C}_8'\text{-C}_9\text{-C}_9'$ Four-Spin Pattern, N-H Cyclopeptide . . .	92
Figure 1.43	Approximate Chemical Shift and Coupling Constant Pattern of the $\text{C}_1\text{-C}_1'\text{-C}_2\text{-C}_2'$ Four-Spin Pattern, N-Methyl Cyclopeptide .	94
Figure 1.44	Calculated Chemical Shift and Coupling Constant Pattern of the $\text{C}_1\text{-C}_1'\text{-C}_2\text{-C}_2'$ Four-Spin Pattern, N-Methyl Cyclopeptide .	97
Figure 1.45	Calculated Spectrum for the $\text{C}_1\text{-C}_1'\text{-C}_2\text{-C}_2'$ Four-Spin Pattern, N-Methyl Cyclopeptide .	98
Figure 1.46	Approximate Chemical Shift and Coupling Constant Pattern of the $\text{C}_1\text{-C}_1'\text{-C}_2\text{-C}_2'$ Four-Spin Pattern, N-H Cyclopeptide . . .	100
Figure 1.47	Calculated Chemical Shift and Coupling Constant Pattern of the $\text{C}_1\text{-C}_1'\text{-C}_2\text{-C}_2'$ Four-Spin Pattern, N-H Cyclopeptide . . .	103

Figure 1.48	Calculated Spectrum for the $C_1-C_1'-C_2-C_2'$ Four-Spin Pattern, N-H Cyclopeptide . . .	104
Figure 1.49	Approximate Chemical Shift and Coupling Constant Pattern of the $C_5-C_{17}-C_{17}'-$ $C_{18}-C_{18}'-C_{19}-C_{19}'$ Seven-Spin Prolyl Pattern, N-Methyl Cyclopeptide	106
Figure 1.50	Approximate Calculated Spectrum for the $C_5-C_{17}-C_{17}'-C_{18}-C_{18}'-C_{19}-C_{19}'$ Seven-Spin Prolyl Pattern, N-Methyl Cyclopeptide	107
Figure 1.51	Approximate Chemical Shift and Coupling Constant Pattern of the $C_{12}-C_{13}-C_{15}-C_{16}$ Four-Spin Aromatic Pattern, N-Methyl Cyclopeptide	109
Figure 1.52	Calculated Chemical Shift and Coupling Constant Pattern of the $C_{12}-C_{13}-C_{15}-C_{16}$ Four-Spin Aromatic Pattern, N-Methyl Cyclopeptide	112
Figure 1.53	Calculated Spectrum for the $C_{12}-C_{13}-C_{15}-C_{16}$ Four-Spin Aromatic Pattern, N-Methyl Cyclopeptide	113
Figure 1.54	Approximate Chemical Shift and Coupling Constant Pattern of the $C_{12}-C_{13}-C_{15}-C_{16}$ Four-Spin Aromatic Pattern, N-H Cyclopeptide	115
Figure 1.55	Calculated Chemical Shift and Coupling Constant Pattern of the $C_{12}-C_{13}-C_{15}-C_{16}$ Four-Spin Aromatic Pattern, N-H Cyclopeptide	118
Figure 1.56	Calculated Spectrum for the $C_{12}-C_{13}-C_{15}-C_{16}$ Four-Spin Aromatic Pattern, N-H Cyclopeptide	119
Figure 1.57	270 MHz 1H NMR Spectrum of <i>Cyclo</i> -[3-(4- β - aminoethyl)phenoxy-4-methylpentanoyl- L-prolyl], Showing Line Assignments . . .	121
Figure 1.58	270 MHz 1H Two-dimensional Homonuclear J- Spectrum of <i>Cyclo</i> -[3-(4- β -aminoethyl)- phenoxy-4-methylpentanoyl-L-prolyl] . .	122

Figure 1.59	270 MHz ^1H NMR Spectrum of <i>Cyelo</i> -[3-(4- β -aminoethyl)phenoxy-4-methylpentanoyl-L-prolyl] -- 0° Projection Sum of the 2-D Homonuclear J-Spectrum	123
Figure 1.60	270 MHz ^1H NMR Spectrum of <i>Cyelo</i> -[3-(4- β -aminoethyl)phenoxy-4-methylpentanoyl-L-prolyl] -- 45° Projection Sum of the 2-D Homonuclear J-Spectrum	124
Figure 1.61	Metal Binding Behavior of <i>Cyelo</i> -[3-(4- β -N-methylaminoethyl)phenoxypropanoyl-L-prolyl] as shown by Circular Dichroism .	126
Figure 2.1	8 - 270 MHz Spectrometer System Block Diagram	140
Figure 2.2	Phase Shift Modulator Scheme for Single Sideband Frequency Conversion	141
Figure 2.3	Superheterodyne Receiver Structure	144
Figure 2.4	8 - 270 MHz Spectrometer System, Observe Channel Transmitter Interconnect [16X9673-B3]	148
Figure 2.5	8 - 270 Mhz Spectrometer System, Observe Channel Receiver Interconnect [16X9673-B4]	152
Figure 2.6	8 - 270 MHz Spectrometer System, Decouple Channel Transmitter Interconnect for Homonuclear Decoupling (^1H) [16X9673-B7]	157
Figure 2.7	8 - 270 MHz Spectrometer System, Decouple Channel Transmitter Interconnect for Heteronuclear Proton Decoupling [16X9673-B8]	158
Figure 2.8	8 - 270 MHz Spectrometer System, ^2H Lock System Interconnect [16X9673-B6]	163
Figure 2.9	Matching of a Parallel Resonant Circuit with a Passive Lossless Two-Port	170
Figure 2.10	5 mm $^1\text{H}/^2\text{H}$ Probe: 270 MHz Port and 41.45 MHz Port Circuit Details	172

Figure 2.11	10 mm ^{13}C Probe: 270 MHz Port Circuit Detail	173
Figure 2.12	10 mm ^{13}C Probe: 67.90 MHz Port and 41.45 MHz Port Circuit Details	174
Figure 2.13	10 mm ^{31}P Probe: 270 MHz Port Circuit Detail	175
Figure 2.14	10 mm ^{31}P Probe: 109.30 MHz Port and 41.45 MHz Port Circuit Details	176
Figure 2.15	15 mm ^{15}N Probe: 270 MHz Port Circuit Detail	177
Figure 2.16	15 mm ^{15}N Probe: 27.36 MHz Port and 41.45 MHz Port Circuit Details	178
Figure 2.17	10 mm $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ Probe: 270 - 286 MHz Port Circuit Detail	179
Figure 2.18	10 mm $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ Probe: 67.90 MHz Port and 41.45 MHz Port Circuit Details	180
Figure 2.19	$^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ 10 mm Probe, Insert Detail and 67.89 MHz/41.45 MHz Matching Network Detail	182
Figure 2.20	$^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ 10 mm Probe, Insert Detail and 270 MHz Matching Network Detail	183
Figure 2.21	Typical High Frequency Coil Cutting Pattern, Dimensions Shown for Placement on the Inside of a 13 mm O.D. Insert.	185
Figure 2.22	8 - 270 MHz Spectrometer System, Digital Interconnect [16X9673-B1]	189
Figure 2.23	Angle θ and Magnitude E_r/E_i of Incident and Reflected Waves as Determined by Load Impedance	196
Figure 2.24	Complex Impedance Measurement Using HP 8405A Vector Voltmeter and HP 778D Dual Directional Coupler	197
Figure 2.25	Swept Frequency VSWR Measurement Using Wiltron 610B Frequency Sweeper and VSWR Bridge	198
Figure 2.26	Cascade of Subsystems Making Up a System, with Signal-to-Noise Ratios Defined at Each Point	203

Figure 2.27	270 MHz Spectrometer System/Nicolet Instrument 293A Interface	213
Figure 2.28	270 MHz Spectrometer System/Nicolet Instrument 293A Interface, State Diagram and Transition Table	214
Figure 2.29	J1 Karnaugh Map	217
Figure 2.30	K1 Karnaugh Map	218
Figure 2.31	J0 Karnaugh Map	219
Figure 2.32	K0 Karnaugh Map	220
Figure 2.33	270 MHz Spectrometer System, Lock Channel 5 KHz Four Phase Generator (Board 1), and Gate Generator (Board 2) [16X8943-S1]	224
Figure 2.34	270 MHz Spectrometer System, Lock Channel 5 KHz Bandpass Filter [16X8943-S2] . . .	227
Figure 2.35	270 MHz Spectrometer System, Lock Channel Dispersion Signal Phase Detector and Lock Control Signal Processing (Board 4) [16X8943-S3]	228
Figure 2.36	270 MHz Spectrometer System, Lock Channel Absorption Signal Phase Detector (Board 5) [16X8943-S4]	231
Figure 2.37	270 MHz Spectrometer System, Lock Channel Sweep Control Logic (Board 7) [16X8943-S5]	232
Figure 2.38	270 MHz Spectrometer System, Lock Channel Analog Buffer (Board 6) and Sweep Generator (Board 9) [16X8943-S6]	233
Figure 2.39	Remote Controlled Low-Pass Filter, Analog Signal Processing [16X9023-S1]	238
Figure 2.40	Remote Controlled Low-Pass Filter, Remote/ Local Function Selection [16X9023-S2A] .	239
Figure 2.41	Remote Controlled Low-Pass Filter, Device Select and Remote Controller [16X9023-S3A]	240
Figure 2.42	Remote Controlled Low-Pass Filter, Differential Amplifier Buffer [16X9023-S4]	241

Figure 2.43	Remote Controlled Low-Pass Filter, LED Remote Frequency and Mode Indicator [16X9023-S5]	242
Figure 2.44	Low Frequency Rf Preamplifier [16X9033] . . .	245
Figure 2.45	Construction Details, Low Frequency Rf Preamplifier	246
Figure 2.46	Lock Channel Single Sideband Generator [16X9193]	257
Figure 2.47	16X919 Lock Channel SSB Generator, Internal Construction Detail (1)	258
Figure 2.48	16X919 Lock Channel SSB Generator, Internal Construction Detail (2)	259
Figure 2.49	16X919 Lock Channel SSB Generator, Internal Construction Detail (3)	260
Figure 2.50	63 KG Magnet System -- Room Temperature Shim Control, Front Panel Wiring [16X9353-S1A]	263
Figure 2.51	63 KG Magnet System -- Room Temperature Shim Control, Z/H ₀ Summing Amplifier Board [16X9353-S2]	264
Figure 2.52	8 - 270 MHz Spectrometer System, Probe Temperature Control System Interconnect [16X9503-B1]	271
Figure 2.53	Temperature Control Device Select and Interrupt Generator [16X9503-S1]	272
Figure 2.54	Temperature Control 20 Bit D-Register Board [16X9503-S2]	273
Figure 2.55	Temperature Control Analog Buffer and Status Indicator Board [16X9503-S3] . . .	274
Figure 2.56	Assembled Listing of Temperature Control Interrupt Service Overlay to NTCFT-1180 .	276
Figure 2.57	270 MHz Preamplifier [16X9553]	282
Figure 2.58	Nicolet 1180/CalComp 565 Interface, Plot Control Logic [16X9563-S1]	287
Figure 2.59	Nicolet 1180/CalComp 565 Interface, Line Driver Detail [16X9563-S2]	288

Figure 2.60	Assembled Listing of CalComp 565 Modifications to NTCFT-1180	291
Figure 2.61	41.45 MHz Receiver - T/R Switch [16X9643]	297
Figure 2.62	27.36 MHz / 270 MHz Dual Frequency Receiver - T/R Switch Assembly [16X9653]	301
Figure 2.63	67.89 MHz / 109.30 MHz Dual Frequency Receiver - T/R Switch Assembly [16X9663]	305
Figure 2.64	Nicolet 1180 Data System -- Disk Memory Subsystems DC Supplies and Signal Interconnections [16X9673-B2]	309
Figure 2.65	8 - 270 MHz Spectrometer System, Rf Power Amplifier Connections [16X9673-B5]	310
Figure 2.66	DPDT Rf Switch [16X9733-S1]	311
Figure 2.67	16X973 Rf Switch PC Layout (Signal Side) . . .	312
Figure 2.68	16X973 Rf Switch PC Layout (Ground Plane Side)	313
Figure 2.69	Preamplifier Power Supply [16X9743]	314
Figure 2.70	41.45 MHz Amplifier [16X9783]	315
Figure 2.71	Rf Switch [16X9803]	316
Figure 2.72	Data Processor Power Fail Interlock [16X9813]	319

LIST OF TABLES

Table 1.1	NMR Line Assignments for <i>Cyclo</i> -[3-(4- β -N-methylaminoethyl)phenoxypropanoyl-L-prolyl]	65
Table 1.2	NMR Line Assignments for <i>Cyclo</i> -[3-(4- β -aminoethyl)phenoxypropanoyl-L-prolyl]	68
Table 1.3	Cyclopeptide Prolyl Seven-Spin System: Summary of Chemical Shift Differences Between N-Methyl and N-H Cases	72
Table 1.4	N-Methyl Cyclopeptide -- Dihedral Angles from X-ray Data	80
Table 1.5	Calculated Transition Frequencies and Assignment of Observed Lines for the Simulation of the C ₈ -C ₈ '-C ₉ -C ₉ ' Four-Spin System, N-Methyl Cyclopeptide	83
Table 1.6	Calculated Transition Frequencies and Assignment of Observed Lines for the Simulation of the C ₈ -C ₈ '-C ₉ -C ₉ ' Four-Spin System, N-H Cyclopeptide	89
Table 1.7	Calculated Transition Frequencies and Assignment of Observed Lines for the Simulation of the C ₁ -C ₁ '-C ₂ -C ₂ ' Four-Spin System, N-Methyl Cyclopeptide	95
Table 1.8	Calculated Transition Frequencies and Assignment of Observed Lines for the Simulation of the C ₁ -C ₁ '-C ₂ -C ₂ ' Four-Spin System, N-H Cyclopeptide	101
Table 1.9	Calculated Transition Frequencies and Assignment of Observed Lines for the C ₁₂ -C ₁₃ -C ₁₅ -C ₁₆ Four-Spin Aromatic System, N-Methyl Cyclopeptide	110
Table 1.10	Calculated Transition Frequencies and Assignment of Observed Lines for the C ₁₂ -C ₁₃ -C ₁₅ -C ₁₆ Four-Spin Aromatic System, N-H Cyclopeptide	116
Table 1.11	Comparison of Bond Angles -- N-Methyl Cyclo- peptide (X-ray Data) and N-H Cyclopeptide (Calculated)	129

Table 2.1	8 - 270 MHz Spectrometer System, Principal Operating Frequencies	146
Table 2.2	Observe Channel Receiver Characteristics	153
Table 2.3	8 - 270 MHz Spectrometer System: Sample Probe Specifications	168
Table 2.4	Observe Channel System Noise Performance	206
Table 2.5	Summary of 1180/293A I/O Functions	210
Table 2.6	16X892 Transition Table	216
Table 2.7	Nicolet 1180 I/O Bus Intercept Wire-Wrap Cross Reference	250
Table 2.8	16X956 Nicolet 1180/CalComp 565 Plot Command Summary	284
Table A1.1	1180 Data System -- Resident Monitor Reserved Memory Locations	327
Table A1.2	Interrupt Vector Addresses and Device Priority Level Assignment -- Nicolet 1180 Data System	328

I. Conformational Analysis of Cyclic Peptides

"Nothing will ever be attempted if all possible objections must first be overcome."

- Samuel Johnson

1.1 Introduction

The physical properties and biological activity of proteins is primarily dependent on levels of structure higher than just the primary amino acid sequence. The analysis of peptide backbone conformations therefore is an important route to the understanding of the three-dimensional arrangement of amino acids in proteins. In recent years, a great deal of attention has been focused on conformational studies of small cyclic polypeptides, both natural and synthetic. These substances, in addition to being interesting because of their own biological activity, have been used as models for larger proteins because of the simplified analysis that is possible.

Of the spectroscopic techniques that are applicable to conformational analysis, nuclear magnetic resonance (NMR) has found increasing use.¹⁻⁵ In part this is due to the wealth of information that is provided by the NMR spectrum, *i.e.* chemical shifts and spin-spin couplings. Proton NMR has been used for the majority of these studies, primarily because of the ease of acquiring and analyzing data. Proton spectral parameters which can provide the most information with regard to conformation are:

1. Vicinal coupling constants, which can be interpreted in terms of Karplus type relationships:^{6,7}

$$^3J(\theta) = A \cos^2\theta + B \cos\theta + C \quad (1)$$

where θ is the dihedral angle between the planes formed by H-C₁-C₂ and C₁-C₂-H, and A = 4.22, B = -0.5, C = 4.5.

Strictly speaking, (1) is only valid for saturated tetrahedral carbons. The presence of electronegative substituents alters this relationship; many groups have assumed that the general form of the equation is correct in such instances, and they have empirically added or modified constants. The simple relationship (1) is often restated in the form

$$^3J(\theta) = K_1 \cos^2\theta - 0.28 \quad 0^\circ \leq \theta \leq 90^\circ \quad (2)$$

$$^3J(\theta) = K_2 \cos^2\theta - 0.28 \quad 90^\circ \leq \theta \leq 180^\circ \quad (3)$$

with $K_1 = 8.5$ and $K_2 = 9.5$.

Another variation on (1) has been proposed for the dependence of $^3J_{\text{NH-C}^\alpha\text{H}}$ on dihedral angle:^{8,9}

$$^3J_{\text{NH-C}^\alpha\text{H}} = A \cos^2\theta + B \cos \theta + C \sin^2\theta \quad (4)$$


The values $A = 7.9$, $B = -1.55$, and $C = 1.35$ gave a reasonable fit for cyclic oligopeptides where both NMR and other types of structural data were available.

It is readily apparent that the application of coupling constant data to the analysis of bond angles and solution conformations depends on the use of good model systems for the evaluation of constants.

2. Chemical shift assignments, which are of more limited utility in systems of higher complexity because of the difficulty in making assignments. The advent of high field superconducting magnet spectrometers has helped to spread out

resonances, making this task easier. An important aspect of making assignments is that line positions together with coupling constant estimates make theoretical simulation of spectra possible.

3. Temperature dependence of NH chemical shifts -- in polar protic solvents, amide protons that are intramolecularly hydrogen bonded tend to experience less of an upfield shift than amide protons that are exposed to the solvent.
4. *Cis-trans* peptide bond isomerism -- the planer nature of the peptide bond, owing to resonance stabilization



means that *cis* and *trans* peptide bond configurations are possible. NMR can often reveal the presence of both forms due to the generally high barrier to internal rotation.

Natural abundance ^{13}C NMR has also come into wider use for solution spectra of small polypeptides. The principal advantage is the greater range of chemical shifts (> 200 PPM) due to the increased sensitivity of shifts to the local environment. Most ^{13}C work has employed broadband decoupling of the coupled protons to improve the signal to noise, as a result, ^{13}C - ^1H J couplings have not been revealed. The low natural abundance (1.1%) of ^{13}C implies that ^{13}C - ^{13}C couplings are not

generally seen. The primary utility of ^{13}C NMR has been looking at chemical shift variation resulting from conformational effects.

In the text that follows, a different approach to conformational analysis using ^1H NMR will be delineated. One of the great difficulties in using ^1H NMR for conformational studies has been the extraction of chemical shift and coupling information, because of the multitude of lines and the incidence of non-first-order coupled systems. The approach that has been taken employs the acquisition of two-dimensional NMR J spectra and the generation of suitable projections to extract shifts and couplings, followed by computer simulation of spin subsystems to generate accurate coupling constants for bond angle calculations.

1.2 Two-Dimensional J-Resolved Spectroscopy

In most NMR experiments, the relative intensities of the observed lines are subject to a variety of external influences. Observed intensities may show a time dependence due to relaxation processes, or observed linewidths (and relaxation rates) may show variation with field strength or other parameters. In such instances, the logical extension for the presentation of spectral data has been the utilization of a second dimension to represent this other functional dependence.

The idea of applying a two-dimensional Fourier transform to experiments where the external variable is time was first suggested by Jeener in 1971.¹⁰ Since that time, an increasing volume of work has appeared, along with several reviews.¹¹⁻¹⁶ The concept is quite simple; normal one-dimensional spectra are acquired as a time domain decay of a transient free precession signal, which may be called $S(t_2)$ where t_2 is the time variable associated with the data acquisition process. Fourier transformation of this signal yields the normal frequency domain spectrum, $S(f_2)$. The appearance of $S(f_2)$, however, depends on the behavior of the nuclear magnetization prior to time $t_2 = 0$. If as a result of some pulse sequence, the nuclear spins execute some motion during this time (called the evolution time t_1), this history will be reflected in the relative phases and amplitudes of the components of $S(t_2)$. The systematic variation of the evolution time t_1 over an appropriate range will then yield a signal matrix $S(t_1, t_2)$ which upon double Fourier transformation yields a two-dimensional spectrum $S(f_1, f_2)$ in two orthogonal frequency dimensions that correlates the behavior of spins during the evolution time t_1 to their behavior during detection





(t_2) .

The purpose of acquiring data in two dimensions is to try to spread out complex spectra that are not amenable to analysis in the conventional fashion. The feature of 2-D that makes this possible is that it is often possible to restrict different nuclear parameters to distinct dimensions. This separation is exemplified by the case of echo modulation by scalar coupling that was first described by Hahn and Maxwell in the basic 90° - τ - 180° - τ spin-echo experiment (Figure 1.1).¹⁷ For a coupled homonuclear spin system, the 180° pulse produces a refocusing of magnetization components that are dephased by chemical shift and field inhomogeneity effects, but the J splitting is not refocused because the 180° pulse also interchanges the spin state labels of the multiplet components. This is conveniently illustrated by the phase-time diagram for a coupled two-spin system IS of Figure 1.2 (after Freeman).¹⁵ The phase shift of the two vectors representing the I spins in the rotating frame are shown as a function of time during the defocusing and refocusing intervals. If the 180° pulse is applied only to the I spins (top), the signs of the accumulated phase errors are changed, resulting in reconvergence at the end of the refocusing interval with zero phase error. If on the other hand the 180° pulse is simultaneously applied to the S spins (bottom), the two vectors are also interchanged. At the end of the refocusing period, the accumulated phase errors are $\pm 1/2 J_{IS}t_1$.



The J-resolved two-dimensional experiment may be summarized by Figure 1.3. There is an evolution period of length t_1 followed by a detection period with running time variable t_2 . During these two times, different effective Hamiltonians apply. During t_2 , the

Figure 1.1 90° - τ - 180° - τ Spin-Echo Experiment

CARR-PURCELL SPIN ECHO SEQUENCE

1. 90° PULSE ALONG X AT TIME 02. MACROSCOPIC MAGNETIZATIONS
DEPHASE3. 180° PULSE ALONG X AT TIME τ 

4. REPHASING

5. ECHO AT TIME 2τ 

6. DEPHASING

Hamiltonian is just that for a spin system under high resolution conditions:

$$\mathcal{H}^{(2)} = \sum_i h\gamma_i(1 - \sigma_i) H \cdot \mathbf{I}_i + \sum_{i < j} \sum J_{ij} \mathbf{I}_i \cdot \mathbf{I}_j \quad (6)$$

which consists of chemical shielding terms and spin-spin coupling terms. During the evolution period t_1 , chemical shift effects are eliminated by the refocusing process, so the effective Hamiltonian is only the spin-spin coupling term:

$$\mathcal{H}^{(1)} = \sum_{i < j} \sum J_{ij} \mathbf{I}_i \cdot \mathbf{I}_j \quad (7)$$



Thus the f_1 dimension of the transformation should reflect only spin coupling effects, while the f_2 dimension should show both spin couplings and chemical shift effects.


Figure 1.2 Phase-Time Diagram for a Coupled Two-spin System, IS

Top: The two components of the I spin multiplet are exactly refocused at the end of the evolution period.

Lower: The 180° pulse applied to the S spins leads to phase modulation at the end of the evolution period.

Figure 1.3 J-Resolved Two-dimensional Experiment





1.3 Cyclopeptide Alkaloids


A. Occurrence, Structure, Model Systems

The term *peptide alkaloids*, first suggested by Goutarel, et al,¹⁸ has been used to describe a large class of polyamide plant bases that are particularly common in plants of the Rhamnaceae family.¹⁹ These compounds often represent *ca.* 0.02 to 0.9% of the dry weight of the plant, and they are usually distributed throughout the plant although often they are more abundant in the root bark and seeds. Isolation procedures produce complex alkaloidal mixtures, the composition of which may depend on the relative maturity and the region of growth of the plant.


Most peptide alkaloids that have been isolated are cyclic, and it has been suggested that the term *cyclopeptide alkaloids* is a more apt description of the class.²⁰

The largest group of these compounds are based on 14 membered rings; they are characterized by a β -hydroxyamino acid unit, and an aryl ether linkage to a *p*-hydroxystyrylamine unit. The β -hydroxyamino acid may be β -hydroxyleucine (as in frangulanine, 1), β -hydroxyphenylalanine (as in integerrine, 2), *trans*-3-hydroxyproline (amphibine-B type, 3), β -hydroxyisoleucine (*e.g.* ceanothine-B, 4), among the more prevalent types. Functionalization of the benzylic position may also vary (5).


NMR has been one of the important techniques used for the structure elucidation work on the cyclopeptides. Its use has been far more limited, though, in the analysis of solution conformations




1
~





2
~



3
~



4
~




5
~

and conformational changes in these compounds. The latter question is of interest because of the suggestion that cyclopeptides may function as ionophores in plants,^{22,23} and that they undergo metal induced conformational changes.²⁴

The analysis of metal induced conformational changes in cyclopeptide alkaloids has been addressed via the study of models for the

natural products: *cyclo*-[3-(4- β -N-methylaminoethyl)phenoxypropanoyl-L-prolyl], 6a, *cyclo*-[3-(4- β -aminoethyl)phenoxypropanoyl-L-prolyl], 6b, and *cyclo*-[3-(4- β -aminoethyl)phenoxy-4-methylpentanoyl-L-prolyl], 6c. These compounds differ from the ring backbone of most of the natural products in the saturated C₁-C₂ linkage (which results in fewer conformational possibilities) and the substitution of an L-proline residue for the β -hydroxy- α -amino acid residue (which eliminates the need to separate diastereomers during synthesis). All three compounds were prepared via high dilution active ester cyclizations from the *p*-nitrophenyl esters of the corresponding *tert*-butoxycarbonyl protected L-prolyl-(4- β -aminoethyl)phenoxypropanoates.²² The *N-tert*-butoxycarbonyl group was removed by dissolving in anhydrous trifluoroacetic acid; the resulting amine salt was dissolved in N,N-dimethylacetamide and added slowly to pyridine maintained at 90° C, giving cyclic monomers (*ca.* 25%) contaminated by cyclic dimers and oligomers. Pure cyclopeptide monomers were isolated by chromatography, then sublimed prior to NMR experiments.



6a: R₁ = CH₃, R₂ = H

6b: R₁ = H, R₂ = H

6c: R₁ = H, R₂ = CH(CH₃)₂

B. Metal Binding Properties

Circular dichroism (CD) studies of metal binding by certain natural and synthetic polypeptides have shown that in aprotic solvents like acetonitrile, large shifts in the dichroic absorption are observable upon the addition of certain metals.^{24,25} The natural cyclopeptide alkaloid Ceanothine-B²⁶ exhibits such behavior; Figure 1.4 shows the change in the near UV region (aromatic) observed upon the addition of Ca^{2+} or Mg^{2+} . No change was observed upon addition of Na^+ however.

The CD spectra of the two principal models, the N-methyl cyclopeptide 6a, and the N-H cyclopeptide 6b, are shown in Figures 1.5 and 1.6 (near and far UV regions). The metal binding properties of the N-H case 6b have been reported;²² the observed shift in dichroic absorption and metal selectivity mimicked the natural product surprisingly well (Figure 1.7).


The metal binding properties of the N-H cyclopeptide 6b were consistent with a conformation in which both carbonyl groups were directed towards the same face of the ring. The inability to crystallize the compound, however, made it impossible to obtain an x-ray crystal structure. An x-ray crystal structure was obtained for the N-methyl cyclopeptide 6a, though, and it showed the carbonyl groups at C_4 and C_7 pointing towards opposite sides of the ring (Figure 1.8).

Infrared absorption spectra of 6a and 6b indicated that the two rings might have quite different conformations. The N-methyl case shows a single amide carbonyl stretch at 1630 cm^{-1} , which is consistent with a *trans* configuration of both amide bonds as shown by the x-ray.

The N-H case $\underline{6b}$ on the other hand shows two amide carbonyl absorptions, one at 1675 cm^{-1} and the other at 1615 cm^{-1} . A lower frequency amide absorption is characteristic of a *cis* amide bond, thus the infrared data suggests one *cis* and one *trans* linkage.


Cumulatively, this data suggested that ring conformation might be a key factor in the metal binding properties of the cyclopeptide model systems. If it was necessary for both carbonyls to be on the same side of the ring to participate in metal complexation, the N-methyl cyclopeptide would likely exhibit vastly different or perhaps no metal affinity. Before testing this hypothesis on the small sample of $\underline{6a}$ that had been prepared, a series of NMR experiments was performed to see if detailed solution conformation data could be obtained. If the solution conformations of $\underline{6a}$ and $\underline{6b}$ could be ascertained by NMR, an approach that is very complementary to x-ray crystallography could be demonstrated for small polypeptides.

- Figure 1.4 Metal Binding Behavior of Ceanothine-B as Shown by
Circular Dichroism
- Figure 1.5 Circular Dichroism Spectra of Cyclopeptide Model Systems,
250 - 300 nm Region
- Figure 1.6 Circular Dichroism Spectra of Cyclopeptide Model Systems,
215 - 250 nm Region
- Figure 1.7 Metal Binding Behavior of *Cyclo*-[3-(4- β -aminoethyl)
phenyloxypropanoyl-L-prolyl] as shown by Circular
Dichroism
- Figure 1.8 Stereo X-ray Crystal Structure of *Cyclo*-[3-(4- β -N-
methylaminoethyl)phenyloxypropanoyl-L-prolyl]



- No salt
- · - · Sodium perchlorate $1.10 \times 10^{-3} \text{ M}$
- Magnesium perchlorate $9.21 \times 10^{-4} \text{ M}$
- Calcium perchlorate $1.45 \times 10^{-3} \text{ M}$

XBL 7710-4693




(—1—) Cyclo[3-(4- β -N-methyl aminoethylphenoxy)
propanoyl L-prolyl] ($c = 1.33 \times 10^{-3} \text{M}$)

(-2-) Cyclo[3-(4- β -aminoethylphenoxy)
propanoyl L-prolyl] ($c = 1.31 \times 10^{-3} \text{M}$)


(—3—) Cyclo[3-(4- β -N-methyl aminoethylphenoxy)
propanoyl L-leucyl] ($c = 1.47 \times 10^{-3} \text{M}$)

(-4-) Cyclo[3-4- β -aminoethyl phenoxy]
propanoyl L-leucyl] ($c = 1.09 \times 10^{-3} \text{M}$)

XBL 7710-4686





- (1) Cyclo[3-(4-β-N-methylaminoethylphenoxy)
propanoyl L-prolyl] (c = 1.75 × 10⁻⁴M)
- (2) Cyclo[3-(4-β-aminoethylphenoxy)
propanoyl L-prolyl] (c = 1.31 × 10⁻⁴M)
- (3) Cyclo[3-(4-β-N-methylaminoethylphenoxy)
propanoyl L-leucyl] (c = 1.55 × 10⁻⁴M)
- (4) Cyclo[3-(4-β-aminoethylphenoxy)
propanoyl L-leucyl] (c = 6.68 × 10⁻⁵M)



- No salt
- - - Sodium perchlorate $9.45 \times 10^{-3} M$
- · - Potassium hexafluorophosphate
 $8.61 \times 10^{-3} M$
- · · · · Lithium perchlorate $8.33 \times 10^{-3} M$
- Magnesium perchlorate $9.17 \times 10^{-4} M$
- · - Calcium perchlorate $1.45 \times 10^{-3} M$

XBL 7710-4692



XBL 799-11732

1.4 NMR Spectral Data for Cyclopeptide Alkaloid Model Compounds


A. Assignment of Lines


The 270 MHz ^1H NMR spectra of the cyclopeptides 6a and 6b in CDCl_3 (Figures 1.9 and 1.10 respectively) show perhaps not surprising patterns of complex overlapping multiplets. The only lines that were readily assigned were the aromatic C_{12} - C_{13} - C_{15} - C_{16} four spin multiplet in both cases, and the N-methyl group of 6a (2.95 δ , Figure 1.9).

The traditional approach to making assignments involves exhaustive homonuclear decoupling experiments. The time-consuming part of this task is determining the exact H_2 irradiation frequency in complex regions such as 2.7 - 3.1 δ of the N-H case (6b, Figure 1.10) where integration indicates the presence of four lines. When lines are this closely spaced, the careful application of the correct level of decoupler power is also required. Too high a power level would lead to the collapse of multiplets coupled to adjacent lines, making assignments very difficult. Thus both exact frequency information and a proper decoupling power level are important.


The two-dimensional J-spectra technique described in Section 1.2 was used to make the assignment procedure very simple. Two-dimensional homonuclear J-spectra were acquired with f_1 (J) frequency ranges of ± 25 Hz and f_2 frequency ranges spanning the normal spectral region. The 2-D spectrum of the N-methyl case 6a is shown along with the normal spectrum in Figure 1.11, and the N-H case 6b is similarly shown in Figure 1.12. It may be seen that the complex multiplets of the normal one-dimensional spectra (lower) are effectively rotated by 45° into the


- Figure 1.9 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -N-methylaminoethyl)phenyloxypropanoyl-L-prolyl], 6a , with tentative line assignments
- Figure 1.10 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenyloxypropanoyl-L-prolyl], 6b , with tentative line assignments
- Figure 1.11 270 MHz ^1H Two-dimensional Homonuclear J-spectrum of *Cyclo*-[3-(4- β -N-methylaminoethyl)phenyloxypropanoyl-L-prolyl], 6a , Shown with Normal Spectrum (bottom)
- Figure 1.12 270 MHz ^1H Two-dimensional Homonuclear J-spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenyloxypropanoyl-L-prolyl], 6b , Shown with Normal Spectrum (bottom)





XBL 7910-12517





two-dimensional frequency plane. It is also evident that line widths in the f_1 (J) dimension are narrower because the refocusing process eliminates the line broadening effect of field inhomogeneities.

Even when the 2-D spectra of Figures 1.11 and 1.12 are plotted on greatly enlarged scales, there is a great deal of complexity which is difficult to interpret. Nagayama, et al, have described the utility of cross-sections and projections in the presentation of 2-D data.²⁷ Two visualization schemes were used to aid in the interpretation of the 2-Ds; they are best illustrated with the enlarged 2-D of the 2.6 - 3.1 δ region of the N-H cyclopeptide 6b (Figure 1.12). Figure 1.13 shows this enlargement with a corresponding 0° projection, and Figure 1.14 shows a contour map of exactly the same area. The doublet of doublets of doublets centered at 2.747 ppm is strikingly clear in both, but visualization of the complex non-first-order pattern at 2.82 - 3.03 δ is made much easier by the contour plot. Contour plots of the other spectral regions of 6a and 6b are shown in Figures 1.15 - 1.32. With the aid of these contour plots, a number of tentative assignments could be made (*e.g.* doublets of doublets of doublets were identified as being parts of particular four-spin patterns, based on the size of the couplings).

The full assignment of all lines in both spectra was completed with a series of homonuclear decoupling experiments that made use of the 2-D 45° projections. It has been pointed out that coupled multiplets are effectively rotated 45° into the 2-D frequency plane, thus coupled lines always appear on a 45° diagonal (frequency).²⁸ A 45° projection of a 2-D plot will align coupled multiplets as single lines, yielding what has been termed a proton decoupled proton spectrum.²⁹ 45° and 0°

projection sums of the cyclopeptide 2-D spectra (Figures 1.11 and 1.12) are shown in Figures 1.33 - 1.36, plotted above the respective normal spectra. The simplification that is seen in the proton decoupled proton spectra is quite dramatic; although intensity information is misleading, the collapse of all multiplets into single lines makes the selection of decoupling frequencies simple. Homonuclear decoupling experiments using the irradiation frequencies obtained in this manner led to the rapid, direct assignment of all lines in both normal spectra (6a and 6b). The tentative assignments are shown in Figures 1.9 and 1.10; line frequencies and relative intensities are tabulated in Tables 1.1 and 1.2.

B. Comparison of Assignments

1. Aromatic Four-Spin Multiplet

Although the four-spin patterns exhibited by the C₁₂-C₁₃-C₁₅-C₁₆ aromatic protons would not be expected to change very much in going from the N-methyl cyclopeptide to the N-H case, they appear to be quite different (Figures 1.9 and 1.10). The reason for the difference is revealed by the 45° projections (Figures 1.33 and 1.35). The N-methyl case 6a has three protons that are closely spaced (7.0-7.2 δ) and a single different proton at higher field (6.85 δ), while the N-H case 6b has two protons at \sim 7.2 δ and two at 6.85 δ . This is a rather striking change to result from simply adding a methyl group some distance away.

2. Proline Seven-Spin Systems

The proline $C_5-C_{17}-C_{17}'-C_{18}-C_{18}'-C_{19}-C_{19}'$ seven-spin systems were assigned with the aid of published data on the *cyclo*-tri-L-prolyl system.³⁰ The two C_{19} protons were virtually identical in the N-H and N-methyl cases, but slight differences in chemical shift are noticeable for one of the C_{18} protons, and fairly pronounced differences are seen for the C_{17} and C_5 protons. Examination of models reveals that changes in these shifts are consistent with a reorientation of the C_4 carbonyl. A summary of differences is tabulated in Table 1.3

3. $C_8 - C_9$ Four-Spin System

The $C_8-C_8'-C_9-C_9'$ four-spin patterns should be important indicators of ring conformation, as twisting about the C_8-C_9 axis should be reflected in changes in the vicinal coupling constants. Examination of these four-spin patterns shows minimal differences in the chemical shifts of the lines, but the line patterns are different -- reflecting changes in the coupling constants.

4. $C_1 - C_2$ Four-Spin System

As in the case of the C_8-C_9 four-spin system, the $C_8-C_8'-C_9-C_9'$ four-spin pattern would be expected to reflect rotation about the C_1-C_2 bond axis. Examination of models reveals that inversion of the C_4 carbonyl should markedly change the C_1-C_2 proton dihedral angles and hence the vicinal coupling constants.

One of the C_2 protons moves upfield by 220 Hz in going from the

N-methyl case 6a to the N-H cyclopeptide 6b; otherwise the main difference is in the appearance of the multiplets. Analysis of this region is made difficult by the overlap with lines from other spin systems, and the closeness of the chemical shifts of the two C₁ and one C₂ protons, resulting in highly non-first-order patterns.

Figure 1.13 270 ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl) phenyloxypropanoyl-L-prolyl] in the 2.60 - 3.10 δ Spectral Region --

Top: Two-dimensional Homonuclear J-Spectrum,
J Range = ± 25 Hz

Lower: 0° Projection



Figure 1.14 Contour Plot of the 2.62 - 3.08 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-H Cyclopeptide, 6b

Figure 1.15 Contour Plot of the 7.34 - 6.88 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide, 6a


Figure 1.16 Contour Plot of the 6.98 - 6.52 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide, 6a

Figure 1.17 Contour Plot of the 5.00 - 4.54 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide, 6a


Figure 1.18 Contour Plot of the 4.46 - 4.00 δ Spectral Region of the Two-dimensional Homonuclear J-Spectrum of the N-Methyl Cyclopeptide, 6a





XBL 799-11391




XBL 798-11167









XBL 798-11165




- Figure 1.19 Contour Plot of the 3.71 - 3.25 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-Methyl Cyclopeptide, $\underset{\sim}{6a}$
- Figure 1.20 Contour Plot of the 3.14 - 2.68 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-Methyl Cyclopeptide, $\underset{\sim}{6a}$
- Figure 1.21 Contour Plot of the 2.89 - 2.43 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-Methyl Cyclopeptide, $\underset{\sim}{6a}$
- Figure 1.22 Contour Plot of the 2.46 - 2.00 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-methyl Cyclopeptide, $\underset{\sim}{6a}$
- Figure 1.23 Contour Plot of the 2.06 - 1.60 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-Methyl Cyclopeptide, $\underset{\sim}{6a}$
- Figure 1.24 Contour Plot of the 7.38 - 6.92 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-H Cyclopeptide, $\underset{\sim}{6b}$
- Figure 1.25 Contour Plot of the 7.06 - 6.60 δ Spectral Region of
the Two-dimensional Homonuclear J-Spectrum of the
N-H Cyclopeptide, $\underset{\sim}{6b}$




XBL 798-11162






XBL 798-11161


XBL 798-11154






CHEMICAL SHIFT, PPM


XBL 798-11158





CHEMICAL SHIFT, PPM


XBL 798-11159




CHEMICAL SHIFT, PPM


XBL 798-11157

XBL 798-11153







CHEMICAL SHIFT, PPM




CHEMICAL SHIFT, PP δ



XBL 798-11169





CHEMICAL SHIFT, PPM

XBL 798-11170

Figure 1.33 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -N-methylamino-ethyl)phenoxypropanoyl-L-prolyl], $\underset{\sim}{6a}$

Top: 0° Projection Sum of the 2-D Homonuclear J-Spectrum

Lower: Normal Spectrum

Figure 1.34 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -N-methylamino-ethyl)phenoxypropanoyl-L-prolyl], $\underset{\sim}{6a}$

Top: 45° Projection Sum of the 2-D Homonuclear J-Spectrum

Lower: Normal Spectrum

Figure 1.35 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenoxypropanoyl-L-prolyl], $\underset{\sim}{6b}$


Top: 0° Projection Sum of the 2-D Homonuclear J-Spectrum

Lower: Normal Spectrum


Figure 1.36 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenoxypropanoyl-L-prolyl], $\underset{\sim}{6b}$

Top: 45° Projection Sum of the 2-D Homonuclear J-Spectrum


Lower: Normal Spectrum




XBL 7910-12173



XBL 7910-12174



XBL 7910-12175



XBL 7910-12176

Table 1.1 NMR Line Assignments for *Cyclo*-[3-(4- β -N-methylamino-ethyl)phenoxypropanoyl-L-prolyl], 6a

Frequency (Hz)	Height	Assignment
1923.64	508.22	Complex multiplet, part of aromatic (12,13,15,16-H) 4-spin pattern
1921.47	463.05	
1917.36	620.41	
1915.38	658.96	
1913.24	700.48	
1910.91	644.90	
1903.15	649.96	
1894.90	280.46	
1832.18	453.60	Apparent doublet of doublets, part of aromatic (12,13,15,16-H) 4-spin pattern
1830.02	423.53	
1824.04	378.73	
1821.78	326.85	
1336.69	171.70	Complex multiplet containing 5-H, 2-H, and 9-H
1333.04	164.82	
1324.71	336.52	
1314.12	656.58	
1306.44	757.58	
1295.19	463.07	
1289.29	123.70	
1284.18	324.52	
1151.93	282.38	Doublet of doublet of doublets, 9'-H, part of 8,8',9,9'-H 4-spin pattern
1150.95	301.42	
1146.69	316.43	
1145.40	305.38	
1140.32	282.62	
1139.09	282.42	
1134.95	268.73	
975.73	159.47	19-H, part of prolyl 7-spin system (5,17,17',18,18',19,19'-H)
971.29	204.94	
966.08	380.76	
961.47	387.09	
956.55	286.94	
951.95	242.31	
933.50	233.40	19'-H, part of prolyl 7-spin system (5,17,17',18,18',19,19'-H)
926.91	360.32	
926.08	359.37	
918.32	337.96	
912.77	119.18	
911.90	117.63	
909.41	180.71	

Table 1.1 [continuation]

Frequency (Hz)	Height	Assignment
832.38	164.56	2-H, appearance approximates a doublet of doublets
829.73	225.77	
823.28	411.21	
816.52	309.45	
809.35	443.47	
796.49	3645.98	N-CH ₃
762.89	125.09	1-H, 1'-H, complex multiplet, part of 1,1',2,2'-H 4-spin pattern
754.51	329.52	
747.59	546.43	
742.15	549.86	
738.64	546.70	
735.12	705.05	
732.47	650.02	
727.90	338.29	8-H, appearance close to that of a quartet, part of 8,8',9,9'-H 4-spin pattern
721.07	464.09	
716.15	394.62	
714.83	417.29	
710.84	187.78	
706.21	372.33	
697.45	210.56	18-H, part of prolyl 7-spin system (5,17,17',18,18',19,19'-H)
693.90	184.86	
690.90	216.57	
688.09	275.08	
685.31	246.56	
682.01	236.72	
678.92	263.70	
676.38	231.55	
672.84	171.87	
672.17	175.33	
669.50	176.27	
667.28	147.17	
661.00	106.62	
608.35	382.25	8'-H, dominated by doublet of doublets, part of 8,8',9,9'-H 4-spin pattern
602.90	372.36	
597.90	372.36	
596.81	117.27	
595.60	121.22	
591.92	315.41	
586.57	277.28	

Table 1.1 [continuation]

Frequency (Hz)	Height	Assignment
536.97	130.41	18'-H, 17-H, broad multiplet, part of prolyl 7-spin system (5,17,17',18, 18',19,19'-H)
536.16	132.42	
532.85	166.07	
532.21	168.35	
521.81	404.77	
520.53	405.84	
513.80	501.58	
511.05	532.25	
502.22	363.77	
484.66	190.89	17'-H, complex multiplet, part of prolyl 7-spin system (5,17,17',18, 18',19,19'-H)
481.25	144.91	
474.84	315.30	
471.83	225.74	
466.84	318.36	
462.74	226.13	
457.53	185.63	
454.22	192.06	
445.18	102.12	

Table 1.2 NMR Line Assignments for *Cyclo*-[3-(4- β -aminoethyl)phenyloxypropanoyl-L-prolyl], 6b

Frequency (Hz)	Height	Assignment
1948.33	237.39	(Triplet) Part of aromatic 4-spin (12,13,15,16-H) pattern
1946.94	344.27	
1945.53	253.01	
1939.91	613.11	(Triplet) Part of aromatic 4-spin (12,13,15,16-H) pattern
1938.37	781.76	
1937.15	569.86	
1929.90	756.86	Part of aromatic 4-spin (12,13,15,16-H) pattern
1922.36	218.17	Part of aromatic 4-spin (12,13,15,16-H) pattern
1857.78	174.86	Complex multiplet dominated by triplet, part of aromatic 4-spin (12,13,15, 16-H) pattern
1855.99	203.35	
1850.82	1574.15	
1849.48	1990.71	
1848.13	1298.90	
1843.33	101.10	
1841.81	103.99	
1841.17	102.99	
1838.58	108.40	
1716.71	159.18	(Broad) 3-H (N-H)
1258.74	255.34	9-H, part of 8,8',9,9'-H 4-spin pattern, has appearance of a doublet of doublets
1248.90	387.82	
1248.27	373.31	
1247.28	381.09	
1246.62	388.46	
1237.13	387.37	
1162.89	729.11	5-H and 9'-H; 5-H is part of the prolyl 7-spin system (5,17,17',18,18',19,19'), 9'-H is part of 8,8',9,9'-H 4-spin pattern
1161.79	697.38	
1157.36	462.87	
1151.24	284.61	
1149.77	271.15	
1145.61	276.19	
1144.15	189.60	
1037.37	120.49	2-H, part of 1,1',2,2'-H 4-spin pattern
1033.86	91.38	
1030.38	179.07	
1029.71	218.61	
1026.26	166.60	

Table 1.2 [continuation]

Frequency (Hz)	Height	Assignment
1024.70	210.73	(continued)
1023.19	156.68	2-H, part of 1,1',2,2'-H 4-spin pattern
1021.55	112.39	
1018.65	193.03	
1016.94	188.98	
1015.22	89.97	
952.53	139.01	19-H, part of 5,17,17',18,18',19,19'-H prolyl 7-spin system (doublet of triplets)
950.25	156.05	
943.30	322.82	
941.16	308.50	
934.26	231.16	
932.09	208.43	
907.08	191.66	19'-H, part of 5,17,17',18,18',19,19'-H prolyl 7-spin system (doublet of triplets)
899.84	254.45	
897.27	335.42	
890.07	351.39	
887.54	158.88	
880.35	139.36	
817.63	73.25	
814.27	124.63	
811.62	257.30	1-H, 2'-H, 1'-H
810.21	207.93	
805.01	294.11	Complex multiplet that includes the 1, 1', and 2-H protons of the
799.27	300.12	
796.37	562.73	1,1',2,2'-H 4-spin pattern
794.15	343.18	
791.03	281.71	
788.88	260.84	
786.01	282.87	
784.30	357.82	
782.96	413.61	
782.19	355.95	
778.00	287.44	
776.83	322.34	
775.56	292.93	
774.12	244.95	
771.97	194.70	
770.96	266.65	
766.40	214.16	
764.97	195.23	
760.89	110.23	

Table 1.2 [continuation]

Frequency (Hz)	Height	Assignment
755.95	218.65	8-H, part of 8,8',9,9'-H 4-spin pattern (quartet of doublets)
754.50	209.76	
745.82	230.84	
744.35	202.14	
739.18	274.56	
737.72	261.76	
729.06	264.05	
727.60	196.77	
641.09	178.25	17-H, part of 5,17,17',18,18',19,19'-H prolyl 7-spin system (broad quartet)
634.52	217.49	
628.95	190.63	
622.27	224.67	
604.77	368.97	8'-H, part of 8,8',9,9'-H 4-spin pattern (broad doublet of doublets)
602.90	135.02	
599.27	370.15	
595.81	115.08	
592.68	148.84	
590.44	209.23	
588.05	357.49	
582.44	389.76	
578.39	126.60	18-H, part of 5,17,17',18,18',19,19'-H prolyl 7-spin system (complex multiplet)
574.64	133.48	
573.92	133.87	
571.83	155.04	
569.61	107.49	
568.74	95.17	
565.12	76.65	
562.84	76.14	
562.04	70.66	
540.47	87.59	
538.69	123.99	18'-H, part of 5,17,17',18,18',19,19'-H prolyl 7-spin system (complex multiplet)
536.60	90.08	
533.54	158.00	
531.66	224.05	
529.55	142.79	
526.44	149.24	
524.49	155.52	
521.75	110.92	
519.23	124.00	

Table 1.2 [continuation]

Frequency (Hz)	Height	Assignment
440.29	100.85	
432.81	160.03	17'-H, part of 5,17,17',18,18',19,19'-H prolyl 7-spin system (complex multiplet)
428.19	199.54	
425.21	119.30	
420.77	258.40	
416.10	97.19	
413.19	149.51	
408.84	84.37	

Table 1.3 Cyclopeptide Prolyl Seven-Spin Systems: Summary of
 Chemical Shift Differences Between N-Methyl (6a)
 and N-H (6b) Cases

Proton	N-Methyl Case Chemical Shift	N-H Case Chemical Shift
C ₅	1310.30 Hz	1157.63 Hz
C ₁₇	511.64 Hz	630.68 Hz
C _{17'}	464.70 Hz	420.13 Hz
C ₁₈	683.37 Hz	577.92 Hz
C _{18'}	518.01 Hz	520.30 Hz
C ₁₉	963.70 Hz	946.25 Hz
C _{19'}	921.48 Hz	893.37 Hz

1.5 Analysis of NMR Data

A. Summary of Approach

Non-first-order coupling patterns are best analyzed by comparison with calculated simulations. In the analysis of the cyclopeptide NMR data, the procedure used may be briefly outlined:

1. Approximate simulation -- Using frequencies derived from the 2-D 45° projections and coupling constants derived from the 2D contour plots, a first rough simulation was calculated. Coupling constants were then adjusted in successive simulations to produce the closest possible match to the actual spectrum. N-Methyl cyclopeptide 6a calculations also employed the dihedral angle data available from the x-ray structure, along with Karplus relations to give relative magnitudes of vicinal couplings.
2. Transition assignment and iterative fit -- After a calculation produced a simulated spectrum that appeared to be close enough to the actual spectrum for transitions to be assigned, an iterative fit was performed using a version of the LAOCN3 program of Castellano and Bothner-By.³¹ Both coupling constants and line frequencies were allowed to vary in the least-squares procedure.
3. Incorrect assignment of transitions -- If the iterative fit routine did not converge, it was assumed that incorrect assignment of some transitions had been made. A different simulation was calculated, and new transition assignments were made.

4. Convergence of iterative fit -- When the iterative fit routine

converged and produced a low RMS error of fit, good values for coupling constants and chemical shifts were presumably obtained. The calculated chemical shifts and couplings were then used to generate a new simulated spectrum, which was plotted for comparison with the actual spectrum. Convergence does not insure that the correct assignments had been made, this is where there is a possibility of ambiguity. In highly non-first-order systems, conclusions are regarded as tentative. In systems approaching first-order (such as the $C_8-C_8'-C_9-C_9'$ four-spin system), it is quite straightforward to make transition assignments, and greater confidence may be placed in the values of couplings so obtained.

B. $C_8 - C_9$ Four-Spin System

Dihedral angles for the protons of the $C_8-C_8'-C_9-C_9'$ four-spin system of the N-methyl cyclopeptide are listed in Table 1.4. The Karplus relations (2) and (3), although not necessarily valid for a $O-CH_2-CH_2-CO-$ configuration, were used to give order of magnitude estimates of the vicinal coupling constants. Combined with estimates of couplings derived from the contour plots of Figure 1.18, 1.21, and 1.22, a first approximate coupling pattern was derived for use in a simulation (Figure 1.37). Of the calculated transitions generated by this pattern, all but two were readily assigned (Table 1.5 lists these assignments). The calculated shifts and coupling constants (Figure 1.38) produced the simulated spectrum of Figure 1.39 (shown with the

normal spectrum for comparison); the two are in excellent agreement.

The $C_8-C_8'-C_9-C_9'$ four-spin system of the N-H cyclopeptide 6b was analyzed in a similar fashion. The approximate values derived from the 2-D 45° projection and contour maps (Figure 1.40) were used to calculate transition frequencies (Table 1.6). All but three transitions were assigned, leading to an RMS error of fit of 0.200 Hz in the iterative solution. The calculated best-fit shifts and coupling constants (Figure 1.41) led to the simulated spectrum of Figure 1.42, again an excellent match to the actual spectrum.

Both bond angles and coupling constants were available for the N-methyl case, therefore different constants for the Karplus relations (1) - (3) could be calculated. These constants would then reflect the unique substitution patterns of the C_8-C_9 carbons in the cyclopeptides, and they would be directly applicable to the calculation of bond angles in the N-H case. Attempts were made to fit the J -- angle data to relationships such as:

$$^3J(\theta) = A \cos^2\theta + B \cos\theta + C \quad (8)$$

and

$$^3J(\theta) = A \cos^2\theta + B \cos\theta + C \sin\theta \quad (9)$$

but the best estimator of the vicinal coupling in this case turned out to be:

$$^3J(\theta) = 9.512 \cos^2\theta + 0.681 \quad (10)$$

Application of (10) to the N-H data gave the following bond angles:

<u>Vicinal Coupling</u>	<u>Calculated J</u>	<u>Calculated Angle</u>
$H_{\alpha}-8-9-H_{\alpha}$	1.457 Hz	73°
$H_{\alpha}-8-9-H_{\beta}$	10.313 Hz	180°
$H_{\beta}-8-9-H_{\alpha}$	5.903 Hz	42°
$H_{\beta}-8-9-H_{\beta}$	-0.053 Hz	90°

C. $C_1 - C_2$ Four-Spin System

The exact same method of analysis used for the $C_8-C_8'-C_9-C_9'$ four-spin system was applied to the analysis of the $C_1-C_1'-C_2-C_2'$ four-spin system. Approximate coupling patterns (Figures 1.43 and 1.46) were derived partly with the aid of 360 MHz 1H spectra, because of the overlap with interfering lines at 270 MHz. Calculated transition frequencies were considerably more difficult to assign in these cases (Tables 1.7 and 1.8), but the calculated simulations (Figures 1.45 and 1.48) are not unreasonable approximations of the observed spectra.

Consideration of the exact couplings derived from the calculations for the N-methyl case (Figure 1.44) led to an expression for 3J in this case:

$$^3J(\theta) = 13.78 \cos^2\theta - 1.14 \quad (11)$$

Application to the calculated couplings for the N-H case (Figure 1.47) led to the following bond angles:

<u>Vicinal Coupling</u>	<u>Calculated J</u>	<u>Calculated Angle</u>
$H_{\alpha}-1-2-H_{\alpha}$	4.999	48°
$H_{\alpha}-1-2-H_{\beta}$	4.296	51°
$H_{\beta}-1-2-H_{\alpha}$	4.261	51°
$H_{\beta}-1-2-H_{\beta}$	15.341	0°

D. Proline Seven-Spin Systems

The structural rigidity of the five-membered ring proline accounts for the similarity of appearance in the seven-spin patterns of the two cyclopeptides. Using chemical shift frequencies derived from the 2-D 45° projections and published coupling constants (Figure 1.49),³⁰ the simulated spectrum of Figure 1.50 was generated. While not an exact match, it was sufficiently close to the actual spectrum that an iterative solution was deemed not worthwhile (because of the large number of transition assignments that would have been necessary).

E. Aromatic Four-Spin Systems

The same procedure outlined above was applied to the analysis of the $C_{12}-C_{13}-C_{15}-C_{16}$ aromatic four-spin systems. Approximate values (Figures 1.51 and 1.54) used for an iterative fit yielded calculated values and spectra (Figure 1.52 - 1.56, Tables 1.9, 1.10).

Despite the availability of x-ray data for the N-methyl cyclopeptide, it was not possible to unambiguously assign any of the aromatic resonances. The patterns are quite unusual though, and a significant change occurs

in chemical shifts in going from the N-methyl cyclopeptide to the N-H case.

F. Analysis of Spectral Data of *Cyclo*-[3-(4- β -aminoethyl)phenoxy-4-methylpentanoyl-L-prolyl], $\underset{\sim}{6c}$

No mention has been made yet of the 9-isopropyl substituted N-H cyclopeptide $\underset{\sim}{6c}$. Examination of the normal and the two-dimensional homonuclear J-spectra and projections reveal a great deal of similarity to the N-H case $\underset{\sim}{6c}$ (Figures 1.57 - 1.60). Coupling patterns and chemical shifts are virtually identical except for the substitution of an isopropyl group for the 9- α -H. Thus the configuration at C₉ appears to be the same as in the natural product Ceanothine-B, $\underset{\sim}{4}$. The NMR spectrum indicates that the cyclization that produced $\underset{\sim}{6c}$ resulted in the predominant formation of this isomer (essentially the only isomer by NMR). This is significant from a synthetic standpoint because no resolution of isomers was done prior to the cyclization; only the desired product was formed.

G. Metal Affinity of N-Methyl Cyclopeptide

The NMR data presented above clearly demonstrate a marked difference in conformation between the N-H cyclopeptide $\underset{\sim}{6b}$ and the N-methyl case $\underset{\sim}{6a}$. It was hypothesized that the arrangement of the C₄ and C₇ carbonyls might be a key factor in metal binding, and that because the carbonyls face opposite sides of the ring in the N-methyl case, distinctly different metal affinities would be displayed.

The circular dichroism spectrum of the N-methyl cyclopeptide in the

presence of Na^+ and Mg^{2+} is shown in Figure 1.61. Clearly there is a negligible change in the CD upon the addition of Na^+ , and there is a slight change upon the addition of Mg^{2+} . As expected then, the N-methyl case 6a exhibits only slight affinity for the divalent Mg^{2+} in comparison to the N-H case 6b (Figure 1.7). This is consistent with the proposal that the C_4 - C_7 carbonyls face the same side of the ring in the N-H cyclopeptide.

Table 1.4 N-Methyl Cyclopeptide -- Dihedral Angles from
X-ray Data

Atoms	Angle	8.5 $\cos^2\phi$ - 0.28 (0-90°)
		9.5 $\cos^2\phi$ - 0.28 (90-180°)
H_α -9-8- H_α	84.3°	-0.2
H_α -9-8- H_β	33.8°	5.589
H_β -9-8- H_α	-157.6°	7.840
H_β -9-8- H_β	84.4°	-0.199
H_α -1-2- H_α	-42.5°	4.340
H_α -1-2- H_β	76.0°	0.217
H_β -1-2- H_α	-160.7°	8.182
H_β -1-2- H_β	-42.5°	4.340

Figure 1.37 Approximate Chemical Shift and Coupling Constant
Pattern of the $C_8-C_8'-C_9-C_9'$ Four-Spin Pattern,
N-Methyl Cyclopeptide (6a), Used in Simulation
Calculations





Table 1.5 Calculated Transition Frequencies and Assignment of
Observed Lines for the Simulation of the
C₈-C₈'-C₉-C₉' Four-Spin System, N-Methyl
Cyclopeptide

Transition From To	Calculated Frequency	Calculated Intensity	Observed Frequency	Obs-Calc Error
10 15	586.270	.853	586.560	.289
5 11	586.576	.862	586.560	-.016
3 9	591.738	.866	591.820	.081
1 4	592.046	.876	591.820	-.226
13 16	602.624	1.126	602.770	.145
8 14	602.931	1.118	602.770	-.161
6 12	608.091	1.152	608.190	.098
2 7	608.400	1.144	608.190	-.210
10 13	704.336	1.115	704.350	.013
3 6	705.998	1.129	706.090	.091
5 8	714.506	1.141	714.340	-.166
1 2	716.165	1.155	716.100	-.065
15 16	720.690	.841	720.830	.139
9 12	722.351	.843	722.340	-.011
11 14	730.861	.885		
4 7	732.519	.887		
7 14	1133.140	.934	1133.060	-.080
4 11	1134.798	.932	1134.680	-.118
2 8	1138.609	.912	1138.690	.080
1 5	1140.268	.910	1140.370	.101
12 16	1144.944	1.093	1144.900	-.044
9 15	1146.604	1.079	1146.590	-.014
6 13	1150.411	1.075	1150.470	.058
3 10	1152.072	1.061	1152.090	.017
7 12	1284.014	1.097	1283.900	-.114
2 6	1284.324	1.098	1284.290	-.033
4 9	1294.182	1.056	1294.050	-.132
1 3	1294.491	1.057	1294.790	.298
14 16	1295.818	.938	1295.610	-.208
8 13	1296.125	.935	1296.250	.124
11 15	1305.989	.909	1305.830	-.159
5 10	1306.294	.906	1306.520	.225

Figure 1.38 Calculated Chemical Shift and Coupling Constant
 Pattern of the $C_8-C_8'-C_9-C_9'$ Four-Spin Pattern,
 N-Methyl Cyclopeptide (6a)

Figure 1.39 Top: Calculated Spectrum for the $C_8-C_8'-C_9-C_9'$
 Four-Spin Pattern, N-Methyl Cyclopeptide (6a)

 Lower: Normal Spectrum



XBL 7910-12568




Figure 1.40 Approximate Chemical Shift and Coupling Constant
Pattern of the $C_8-C_8'-C_9-C_9'$ Four-Spin Pattern,
N-H Cyclopeptide (6b), Used in Simulation
Calculations




Table 1.6 Calculated Transition Frequencies and Assignment of
Observed Lines for the Simulation of the
C₈-C₈'-C₉-C₉' Four-Spin System, N-H Cyclopeptide

Transition From To	Calculated Frequency	Calculated Intensity	Observed Frequency	Obs-Calc Error
5 11	582.354	.881	582.490	.135
10 15	582.358	.875	582.490	.131
1 4	588.226	.897	588.080	-.146
3 9	588.228	.891	588.080	-.148
13 16	599.104	1.103	599.270	.165
8 14	599.104	1.097	599.270	.165
6 12	604.970	1.129	604.820	-.150
2 7	604.972	1.123	604.820	-.152
10 13	727.237	1.090	727.350	.112
3 6	728.740	1.102	728.740	-.000
5 8	737.487	1.124	737.490	.002
1 2	738.987	1.136	738.900	-.087
15 16	743.983	.862	744.040	.056
9 12	745.482	.864	745.530	.047
11 14	754.236	.908	754.220	-.016
4 7	755.733	.911	755.620	-.113
7 14	1143.721	.890	1143.879	.158
4 11	1145.218	.896	1145.300	.081
2 8	1149.589	.862	1149.600	.010
1 5	1151.090	.869	1151.050	-.040
12 16	1155.078	1.137	1154.240	-.838
9 15	1156.577	1.117		
6 13	1160.944	1.123	1161.460	.515
3 10	1162.447	1.102	1162.560	.112
7 12	1237.032	1.144	1236.940	-.092
2 6	1237.034	1.150	1236.940	-.094
1 3	1247.281	1.095		
4 9	1247.283	1.090	1247.260	-.023
14 16	1248.389	.896		
8 13	1248.389	.889	1248.360	-.029
5 10	1258.638	.863	1258.760	.121
11 15	1258.642	.869	1258.760	.117




Figure 1.41 Calculated Chemical Shift and Coupling Constant
 Pattern of the $C_8-C_8'-C_9-C_9'$ Four-Spin Pattern,
 N-H Cyclopeptide (6b)
 

Figure 1.42 Top: Calculated Spectrum for the $C_8-C_8'-C_9-C_9'$
 Four-Spin Pattern, N-H Cyclopeptide (6b)
 
 Lower: Normal Spectrum



XBL 7910-12567




Figure 1.43 Approximate Chemical Shift and Coupling Constant
 Pattern of the $C_1-C_1'-C_2-C_2'$ Four-Spin Pattern,
 N-Methyl Cyclopeptide (6a), Used in Simulation
 Calculations





Table 1.7 Calculated Transition Frequencies and Assignment of
Observed Lines for the Simulation of the
C₁-C₁'-C₂-C₂' Four-Spin System, N-Methyl
Cyclopeptide


Transition From To	Calculated Frequency	Calculated Intensity	Observed Frequency	Obs-Calc Error
3 9	726.879	.849		
10 15	726.886	.864		
7 14	727.061	.883		
2 8	727.071	.907		
5 11	732.107	1.164	732.320	.212
1 4	732.117	.603	732.320	.202
4 11	733.746	.735	735.110	1.363
1 5	733.756	1.323	735.110	1.353
6 12	737.833	1.174	738.650	.816
13 16	737.840	1.072	738.650	.809
8 14	743.130	1.160	742.180	-.950
2 7	743.140	1.111	742.180	-.960
6 13	748.314	.850	747.400	-.914
12 16	748.322	.999	747.400	-.922
3 10	754.949	1.179	754.460	-.489
9 15	754.956	1.119	754.460	-.496
10 13	805.057	1.245		
3 6	811.692	.973		
15 16	816.012	.933		
5 8	816.477	1.205		
9 12	822.646	.751		
1 2	823.162	1.070	822.850	-.312
11 14	827.499	.950	829.300	1.800
4 7	834.185	.824	832.670	-1.515
7 12	1297.486	1.069		
2 6	1302.793	1.050		
4 9	1309.024	1.000	1306.520	-2.504
1 3	1314.263	1.003	1313.920	-.343
14 16	1318.747	.993	1318.130	-.617
8 13	1324.036	.974	1324.970	.933
11 15	1330.235	.949	1331.720	1.484
5 10	1335.456	.913	1336.510	1.053

Figure 1.44 Calculated Chemical Shift and Coupling Constant
 Pattern of the $C_1-C_1'-C_2-C_2'$ Four-Spin Pattern,
 N-Methyl Cyclopeptide (6a)

Figure 1.45 Top: Calculated Spectrum for the $C_1-C_1'-C_2-C_2'$
 Four-Spin Pattern, N-Methyl Cyclopeptide (6a)

 Lower: Normal Spectrum





XBL 7910-12518

Figure 1.46 Approximate Chemical Shift and Coupling Constant
Pattern of the $C_1-C_1'-C_2-C_2'$ Four-Spin Pattern,
N-H Cyclopeptide (6b), Used in Simulation
Calculations \sim






Table 1.8 Calculated Transition Frequencies and Assignment of
Observed Lines for the Simulation of the
C₁-C₁'-C₂-C₂' Four-Spin System, N-H Cyclopeptide

Transition From To	Calculated Frequency	Calculated Intensity	Observed Frequency	Obs-Calc Error
9 13	745.004	.041		
10 13	761.800	.548		
3 6	762.644	.697		
5 8	764.239	.524		
7 14	764.854	.202		
1 2	767.629	.755		
4 11	768.244	.762		
15 16	770.222	.486	772.170	1.947
9 12	771.066	.962	772.170	1.103
11 14	778.202	1.937	775.710	-2.492
12 16	779.849	1.488		
9 15	780.693	.106		
2 8	780.896	1.262		
4 7	781.591	.874		
1 5	784.286	1.550	783.630	-.656
10 12	787.862	1.237		
6 13	792.897	1.493	793.170	.272
3 10	793.741	2.204	794.770	1.028
5 11	795.457	2.053	794.770	-.687
10 15	797.488	1.380	797.470	-.018
13 16	805.911	1.041	804.240	-1.671
5 7	808.804	.048		
8 14	809.419	.842	810.690	1.270
3 9	810.537	.108	810.690	.152
1 4	811.499	.682	810.690	-.809
2 11	812.113	.080		
6 12	818.959	.228		
2 7	825.460	.359		
4 10	1001.967	.118		
7 12	1008.238	.590		
2 6	1014.739	1.052		
7 15	1017.864	.476		
4 9	1018.762	.920	1019.620	.857
1 3	1019.724	1.010	1019.620	-.104
11 12	1021.585	.461		
14 16	1023.232	.983	1021.310	-1.922
8 13	1026.740	.942	1025.940	-.800
5 10	1029.179	.842	1029.880	.700
11 15	1031.211	.494	1032.490	1.278
5 9	1045.975	.081		






XBL 7910-12519

Figure 1.49 Approximate Chemical Shift and Coupling Constant
Pattern of the $C_5-C_{17}-C_{17}'-C_{18}-C_{18}'-C_{19}-C_{19}'$
Seven-Spin Prolyl Pattern, N-Methyl Cyclopeptide
(6a), Used in Simulation Calculations

Figure 1.50 Top: Approximate Calculated Spectrum for the
 $C_5-C_{17}-C_{17}'-C_{18}-C_{18}'-C_{19}-C_{19}'$ Seven-Spin
Prolyl Pattern, N-Methyl Cyclopeptide (6a)
Using Values from Figure 1.49

Lower: Normal Spectrum






Figure 1.51 Approximate Chemical Shift and Coupling Constant
 Pattern of the C₁₂-C₁₃-C₁₅-C₁₆ Four-Spin
 Aromatic Pattern, N-Methyl Cyclopeptide (6a),
 Used in Simulation Calculations






Table 1.9 Calculated Transition Frequencies and Assignment of
Observed Lines for the C₁₂-C₁₃-C₁₅-C₁₆ Four-Spin
Aromatic System, N-Methyl Cyclopeptide

Transition From To	Calculated Frequency	Calculated Intensity	Observed Frequency	Obs-Calc Error
9 12	1821.265	.884	1821.380	.114
15 16	1821.292	.887	1821.380	.087
4 7	1823.746	.838	1823.810	.063
11 14	1823.867	.935	1823.810	-.057
4 6	1824.945	.097		
3 7	1828.468	.106		
10 13	1829.544	1.053	1829.590	.045
3 6	1829.668	.946	1829.590	-.078
5 8	1832.063	1.117	1831.990	-.073
1 2	1832.092	1.120	1831.990	-.102
9 15	1892.901	.553	1892.800	-.101
12 16	1892.928	.460	1892.800	-.128
6 13	1894.810	.502	1894.860	.049
3 10	1894.933	.608	1894.860	-.073
6 14	1899.897	.043		
4 11	1900.975	1.287	1901.020	.044
7 14	1901.097	1.254	1901.020	-.077
2 8	1902.978	1.709	1903.140	.161
1 5	1903.007	1.538	1903.140	.132
14 16	1910.411	1.364	1910.590	.178
8 13	1910.472	1.699	1910.590	.117
11 15	1912.986	1.349	1912.970	-.016
5 10	1912.991	1.430	1912.970	-.021
13 16	1915.498	1.287	1915.610	.111
8 14	1915.559	1.128	1915.610	.050
6 12	1917.380	1.347	1917.230	-.150
2 7	1917.441	.565	1917.230	-.211
7 12	1918.579	.225	1918.340	-.239
2 6	1918.640	.846	1918.340	-.300
4 9	1921.060	.683	1921.190	.129
1 3	1921.065	.430	1921.190	.124
10 15	1923.750	.984	1923.800	.049
5 11	1923.755	.987	1923.800	.044
3 9	1925.783	.754	1925.850	.066
1 4	1925.787	.909	1925.850	.062

Figure 1.52 Calculated Chemical Shift and Coupling Constant
 Pattern of the C_{12} - C_{13} - C_{15} - C_{16} Four-Spin
 Aromatic Pattern, N-Methyl Cyclopeptide (6a)

Figure 1.53 Top: Calculated Spectrum for the C_{12} - C_{13} - C_{15} - C_{16}
 Four-Spin Aromatic Pattern, N-Methyl
 Cyclopeptide (6a)
 Lower: Normal Spectrum





XBL 7910-12569

Figure 1.54 Approximate Chemical Shift and Coupling Constant
Pattern of the C₁₂-C₁₃-C₁₅-C₁₆ Four-Spin
Aromatic Pattern, N-H Cyclopeptide (6b),
Used in Simulation Calculations ~






Table 1.10 Calculated Transition Frequencies and Assignment of
Observed Lines for the C₁₂-C₁₃-C₁₅-C₁₆ Four-Spin
Aromatic System, N-H Cyclopeptide

Transition From To	Calculated Frequency	Calculated Intensity	Observed Frequency	Obs-Calc Error
12 16	1848.592	1.941	1848.290	-.302
7 12	1848.727	1.947	1848.290	-.436
11 15	1848.998	1.962	1849.480	.481
4 11	1849.133	1.970	1849.480	.346
6 13	1849.871	1.991	1849.480	-.391
2 6	1849.974	1.995	1849.480	-.494
5 10	1850.276	2.020	1850.720	.443
1 5	1850.382	2.027	1850.720	.337
15 16	1920.499	.473	1920.940	.440
9 14	1920.891	.484	1920.940	.048
11 12	1920.905	.485	1920.940	.034
4 7	1921.311	.525	1920.940	-.371
5 8	1921.810	.051		
10 13	1928.961	1.545	1929.390	.428
3 8	1929.355	1.449	1929.390	.034
5 6	1929.366	1.447	1929.390	.023
1 2	1929.774	1.453	1929.390	-.384
13 16	1936.626	1.563	1936.640	.013
3 6	1936.911	.053		
8 14	1937.889	1.493	1937.860	-.029
6 12	1937.906	1.504	1937.860	-.046
2 7	1939.153	1.436	1939.080	-.073
10 15	1945.088	.492	1944.960	-.128
3 9	1946.353	.490	1946.350	-.003
5 11	1946.367	.505	1946.350	-.017
1 4	1947.616	.507	1947.680	.063

Figure 1.55 Calculated Chemical Shift and Coupling Constant
Pattern of the C₁₂-C₁₃-C₁₅-C₁₆ Four-Spin
Aromatic Pattern, N-H Cyclopeptide (6b)

Figure 1.56 Top: Calculated Spectrum for the C₁₂-C₁₃-C₁₅-C₁₆
 Four-Spin Aromatic Pattern, N-H Cyclopeptide
 (6b)
 Lower: Normal Spectrum





CHEMICAL SHIFT, HZ

XBL 7910-12520

Figure 1.57 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenyloxy-4-methylpentanoyl-L-prolyl] (δC),
Showing Line Assignments

Figure 1.58 270 MHz ^1H Two-dimensional Homonuclear J-Spectrum of
Cyclo-[3-(4- β -aminoethyl)phenyloxy-4-methylpentanoyl-L-prolyl] (δC)

Figure 1.59 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenyloxy-4-methylpentanoyl-L-prolyl] (δC)


Top: 0° Projection Sum of the 2-D Homonuclear J-Spectrum


Lower: Normal Spectrum

Figure 1.60 270 MHz ^1H NMR Spectrum of *Cyclo*-[3-(4- β -aminoethyl)phenyloxy-4-methylpentanoyl-L-prolyl] (δC)


Top: 45° Projection Sum of the 2-D Homonuclear J-Spectrum

Lower: Normal Spectrum





XBL 7910-12172







Figure 1.61 Metal Binding Behavior of *Cyclo*-[3-(4- β -N-methylamino-ethyl)phenoxypropanoyl-L-prolyl], 6a, as shown by Circular Dichroism



- No salt
- - - Sodium perchlorate 10.92×10^{-3} M
- - - Magnesium perchlorate 12.04×10^{-3} M

1.6 Summary

A. Cyclopeptide Alkaloid Model Systems Conformations

The employment of the technique for obtaining two-dimensional homonuclear J-spectra has led to conformational data on model systems for the cyclopeptides. NMR data has shown that the solution conformations of *cyclo*-[3-(4- β -N-methylaminoethyl)phenyloxypropanoyl-L-prolyl] 6a and *cyclo*-[3-(4- β -aminoethyl)phenyloxypropanoyl-L-prolyl] 6b differ markedly, and the change in bond angles as reflected by the vicinal coupling constants are consistent with a major difference in the orientation of the C₄ carbonyl. This difference in conformation accounts for the contrasting metal affinities of the two compounds.

The differences in bond angles are summarized in Table 1.11. The approximate nature of these calculations is reflected by the less than desired internal consistency of the changes. In the absence of large changes in the geminal couplings, it would be expected that H _{α} -H _{α} angles and H _{β} -H _{β} angles would change in the same direction. When changes in the geminal couplings are observed, however, it becomes more difficult to draw conclusions.

B. Methodology

The utility of 2-D homonuclear J-spectra in conformational analysis has been demonstrated for small polypeptides that have relatively few discrete conformational possibilities. The approach that has been used may be briefly reviewed:

1. If a series of compounds is to be studied, it is beneficial to have an x-ray crystal structure for at least one member of the group. The availability of both bond angles and coupling constants allows the derivation of $J - \theta$ relations that reflect the local substituent pattern.
2. Acquire the normal and the 2-D homonuclear J-spectrum.
3. Construct 45° projections of the 2-D spectrum, and plot contour maps of the spectral regions of interest.
4. Use the frequencies derived from the 45° projection for comprehensive homonuclear decoupling experiments to make all assignments.
5. Simulate spin systems where changes in orientation are likely to be reflected by changes in vicinal coupling constants.
6. Using the accurate coupling constant data with angle information from the x-ray, develop predictors for the vicinal coupling constants. This permits the calculation of dihedral angles where the x-ray is not available.

Table 1.11 Comparison of Bond Angles -- N-Methyl Cyclopeptide
 (X-Ray Data) and N-H Cyclopeptide (Calculated)

<u>Dihedral Angle</u>	<u>N-Methyl Case <u>6a</u></u>	<u>N-H Case <u>6b</u></u>
$H_{\alpha}-8-9-H_{\alpha}$	84.3°	73°
$H_{\alpha}-8-9-H_{\beta}$	-157.6°	180°
$H_{\beta}-8-9-H_{\alpha}$	33.8°	42°
$H_{\beta}-8-9-H_{\beta}$	84.4°	90°
$H_{\alpha}-1-2-H_{\alpha}$	-42.5°	48°
$H_{\alpha}-1-2-H_{\beta}$	76.0°	51°
$H_{\beta}-1-2-H_{\alpha}$	-160.7°	51°
$H_{\beta}-1-2-H_{\beta}$	-42.5°	0°

1.7 Experimental Section

Cyclopeptide model compounds $\underset{\sim}{6a}$ - $\underset{\sim}{6c}$ were purified by chromatography on Sephadex LH-20 followed by sublimation at *ca.* 100° C at < 100 μ pressure. NMR samples were prepared with nitric acid washed pipettes and sample tubes (Wilmad 528-PP) using Aldrich "100.00%" (99.9985 atom % D) CDCl₃. Samples were prepared in a Vacuum Atmospheres argon filled glove box, and were filtered through acid washed pipettes with glass wool plugs.

NMR spectra were acquired on the 270 MHz spectrometer system described in Part II. Prior to the acquisition of a 2-D file, the 90° pulse length was measured. This was found to vary from 4.95 - 6.00 μ sec. Acquired spectral widths ranged from ± 1400 Hz to ± 2500 Hz (quadrature detection). The waiting time between the 90° and the 180° pulse was incremented by 10 msec, corresponding to a J spectral width of ± 25 Hz.

Homonuclear decoupling experiments were run as described in Section 2.3.4. In some instances a wider than normal spectral width was used because of the spectrum offset resulting from using one frequency synthesizer to drive both the observe and the decouple transmitters. The decoupler was normally run at a 4 - 20% duty cycle at a pulse delay from the DCLOCK adjusted to minimize DC offsets in the detected signal.

1.8 References

1. W.A. Thomas in "Annual Reports on NMR Spectroscopy," Vol. 6B, E.F. Mooney, Ed., Academic Press, New York, NY., 1976, pp 1-41
2. F.A. Bovey, A.I. Brewster, D.J. Patel, A.E. Tonelli, D.A. Torchia, *Accounts Chem. Res.*, 5, 193 (1972)
3. D.W. Urry and M. Ohnishi in "Spectroscopic Approaches to Biomolecular Conformation," D.W. Urry, Ed., American Medical Association Press, Chicago, IL., 1970, pp 263 ff
4. F.A. Bovey in "Chemistry and Biology of Peptides, Proceedings of the Third American Peptide Symposium," J. Meienhofer, Ed., Ann Arbor Science Publishers, Ann Arbor, MI., 1972, pp 3 ff
5. G.N. Ramachandran and V. Sasisekharan in "Advances in Protein Chemistry," Vol. 23, C.F. Afonso, Ed., Academic Press, New York, NY, 1968, pp 284 ff
6. M. Karplus, *J. Chem. Phys.*, 30, 11 (1959)
7. M. Karplus, *J. Amer. Chem. Soc.*, 85, 2870 (1963)
8. G.N. Ramachandran and R. Chandrasekaran, *Biopolymers*, 10, 935 (1971)
9. G.N. Ramachandran, R. Chandrasekaran, and K.D. Kopple, *Biopolymers*, 10, 2113 (1971)
10. J. Jeener, presented at Ampère International Summer School, Basko Polje, Yugoslavia, 1971, unpublished
11. W.P. Aue, E. Bartholdi, and R.R. Ernst, *J. Chem. Phys.*, 64, 2229 (1976)
12. A. Kumar, D. Welti, and R.R. Ernst, *J. Mag. Resonance*, 18, 69 (1975)
13. L. Müller, A. Kumar, and R.R. Ernst, *J. Chem. Phys.*, 63, 5490 (1975)
14. G. Bodenhausen, R. Freeman, R. Niedermeyer, and D. Turner, *J. Mag. Resonance*, 26, 133 (1977)
15. R. Freeman and G. Morris, *Bull. Mag. Resonance*, 1, 5 (1979)
16. K. Wüthrich, K. Nagayama, and R.R. Ernst, *Trends Biol. Sci.*, 178 (1979)
17. E.L. Hahn and D.E. Maxwell, *Phys. Rev.*, 88, 1070 (1952)

18. M. Pais, X. Monseur, X. Lusinchi, and R. Goutarel, *Bull. Soc. Chim. France*, (1964) 817
19. E.W. Warnhoff, *Fortschr. Chem. Org. Naturst.*, 28, 162 (1970)
20. R. Tschesche and E.V. Kaussman, in "The Alkaloids," Vol. XV, R.H.F. Manske, Ed., Academic Press, New York, NY, 1975, pp 165-203
21. E. Haslinger, *Tetrahedron*, 34, 685 (1978)
22. J.C. Lagarias, R.A. Houghten, and H. Rapoport, *J. Amer. Chem. Soc.*, 100, 8202 (1978)
23. Y. Ogihara, *Experientia*, 33, 1454 (1977)
24. T. Wieland in "Chemistry and Biology of Peptides," J. Meienhofer, Ed., Ann Arbor Science Publishers, Ann Arbor, MI., 1972, p. 377
25. V. Madison, M. Atreyl, C.M. Deber, and E.R. Blout, *J. Amer. Chem. Soc.*, 96, 6725 (1974)
26. F.K. Klein and H. Rapoport, *J. Amer. Chem. Soc.*, 90, 3576 (1968)
27. K. Nagayama, P. Bachmann, K. Wuthrich, and R.R. Ernst, *J. Mag. Resonance*, 31, 133 (1978)
28. 2-D plots do not always appear to deploy multiplets along 45° diagonals because of the difference in plot scaling (Hz/cm) that is usually used.
29. W.P. Aue, J. Karhan, and R.R. Ernst, *J. Chem. Phys.*, 64, 4226 (1976)
30. C.M. Deber, D.A. Torchia, and E.R. Blout, *J. Amer. Chem. Soc.*, 93, 4893 (1971)
31. S. Castellano and A.A. Bothner-By, *J. Chem. Phys.*, 41, 3863 (1964)

II. Multinuclear NMR Spectrometer for the Study of Biological Systems

"The reasonable man adapts himself to the world; the unreasonable one persists in trying to adapt the world to himself. Therefore all progress depends on the unreasonable man."

- George Bernard Shaw

2.1 Introduction

In recent years, nuclear magnetic resonance has been applied to the study of a wide range of biological problems. With the advent of commercially available pulsed Fourier transform spectrometer systems operating at ever higher magnetic field strengths, problems of great complexity have become addressable. It is the purpose of this section to describe the design principles and some of the construction details of a broadband multinuclear NMR spectrometer system that was built in this laboratory. The spectrometer had as a design objective the inclusion of features that would make it especially suited to working with biological systems.

The text that follows will describe particular aspects of the 8 - 270 MHz multinuclear NMR spectrometer at the Chemical Biodynamics Laboratory at Berkeley. The system was used for the majority of the work described in Part 1. What follows may be classified into three categories:

- Overall system design and description of major functional blocks.
- Test and characterization procedures.
- Detailed description of original design work. Parts fabricated at the Electronics Research and Development Group, Lawrence Berkeley Laboratory, are included in functional descriptions but detailed information is not provided. The same applies to commercial equipment.

Most original designs and most significant modifications to

commercial equipment have been assigned Lawrence Berkeley Laboratory (LBL) print numbers (16X series for fabrication done at LCB). All original drawings are on permanent file at the LBL Electronics Engineering Print Room.

2.2 System Design Considerations

When the construction of the 8 - 270 MHz spectrometer system was initiated, the goal was the implementation of a number of specific capabilities:

1. Multiple pulse experiments -- The system had to be capable of generating precisely timed sequences of radio frequency excitation pulses. Different pulses in the sequence might also have different relative phase shifts. This capability would be useful for proton-enhanced, spin-locking, and a variety of other experiments.
2. Broadband observe -- It was desirable to be able to observe a wide range of nuclei using the same frequency generation and observe system. The goal was to construct a system that could observe at any frequency (within the range dictated by the magnetic field strength) while decoupling at any second frequency and optionally locking at a third frequency. The only limitation on frequency (nuclei) combinations would be the availability of an appropriate probe.
3. Flexible decoupling -- for many of the proton experiments that were to be done, an efficient homonuclear decoupling capability would be essential. For high resolution and proton-enhanced ^{13}C and ^{15}N experiments, broadband noise-modulated low and high power proton decoupling would be necessary.

4. Simple structure for generating new pulse sequences -- it was desired that new pulse sequences should be generated via software microprograms rather than hardwired connection of timers. This would permit rapid reconfiguration of experiments and a degree of automation not otherwise possible.

Frequency generation on the 8 - 270 MHz spectrometer system follows the widely used scheme described by Ellett, et al.¹ Observe, decouple, and lock frequencies are generated by single sideband mixing of a 30 MHz intermediate frequency (IF) with local oscillator frequencies supplied by frequency synthesizers (Figure 2.1). Lower sidebands are selected in all three channels, and the actual single sideband conversion is realized by a term-by-term implementation of the trigonometric identity

$$\cos (\omega_1 - \omega_2) = \cos \omega_1 \cos \omega_2 + \sin \omega_1 \sin \omega_2 \quad (1)$$


using quadrature hybrids to perform Hilbert Transforms and power combiners to perform summing operations (Figure 2.2).


On the transmitter side of the spectrometer system, the use of an intermediate frequency makes the generation of pulses with different relative phases quite straightforward. Phase shifting is done at the 30 MHz IF using narrowband 30 MHz quadripole networks; these phase shifts are then frequency translated in the course of single sideband mixing. This circumvents the need for wide bandwidth phase shifters.

The use of an IF considerably simplifies the task of building a sensitive and selective receiver, as the design is just that of a

Figure 2.1 8 - 270 MHz Spectrometer System Block Diagram


Figure 2.2 Phase Shift Modulator Scheme for Single Sideband
Frequency Conversion





superheterodyne receiver. A superheterodyne receiver has good sensitivity and selectivity because the IF filter provides most of the predetection filtering. Since this filter need not be tunable, it can be relatively complex. The same local oscillator frequency used in the transmitter is used to tune the receiver (Figure 2.3). Image frequency ($\omega_{\text{obs}} \pm 2\omega_{\text{IF}}$) rejection may be accomplished either by single sideband mixing (expensive because of the need for two broadband quadrature hybrids for the signal and the local oscillator) or by suitable Rf pre-filtering prior to mixing.

Figure 2.3 Superheterodyne Receiver Structure



2.3 Major Subsystems

2.3.1 Magnet

The spectrometer system is based on a commercial Bruker Instruments 63.5 kGauss superconducting solenoid with a 5 cm room temperature bore. When the room temperature shim coils are installed, 4 cm remains for the probes. Operation at this field strength implies the NMR frequencies listed in Table 2.1.

The magnet system was used with some modification to the room temperature shim control system, the H_0 field sweep amplifiers, and the shim reference voltage supplies. These modifications are described in detail below.

2.3.2 Observe Channel Transmitter


The observe channel transmitter is outlined in Figure 2.4. The 30 MHz IF is derived from a 30 MHz Active Power Divider (10X223) which in turn receives 30 MHz from the Hewlett-Packard 5110A frequency synthesizer driver, which serves as the master frequency reference for the whole system. Phase shifts of 0° , 90° , 180° , or 270° may be imposed on the 30 MHz IF by the 10X209 Four Phase Generator, which operates under the control of the 16X892 data processor interface. Phase shifts are produced by gated selection of one of the four outputs of a Merrimac 18134 quadripole network.

The phase-shifted 30 MHz output is single sideband converted to the desired NMR frequency in the 10X206 Single Sideband Generator, using a local oscillator frequency supplied by one of the HP5105A

Table 2.1 8 - 270 MHz Spectrometer System, Principal Operating Frequencies

Isotope	Natural Abundance	NMR Frequency at 63.42 KGauss	Relative Sensitivity
^1H	99.9844%	269.9982 MHz	1.000
^{31}P	100.0 %	109.297 MHz	0.0664
^{13}C	1.108 %	67.8975 MHz	0.0159
^{15}N	0.365 %	27.360 MHz	0.0010
^2H	0.0156%	41.4513 MHz	0.0097

Figure 2.4 8 - 270 MHz Spectrometer System, Observe Channel
Transmitter Interconnect [16X9673-B3]



frequency synthesizers. Bandwidth limitations of components in the single sideband generator restrict the NMR operating frequency range to 8 - 270 MHz. Pulse gating is also done in the single sideband generator, again under data processor control via the 16X892 interface.

Pulse amplification to a suitable power level is accomplished with either of two sets of Rf power amplifiers, depending on the frequency. $^1\text{H}/^{19}\text{F}$ observe pulses are amplified by the 10X225 270 MHz solid state power amplifier employing cascaded 2N5641/42/43 common emitter stages. The output of this unit is fed to the 10X224 270 MHz Power Amplifier. This amplifier employs an Eimac 8877 (3CX1500A7) high-mu power triode. Final stage maximum output is approximately 1.2 kW PEP at 1% duty cycle.

Rf power amplification in the frequency range of 8 - 110 MHz uses an ENI 420L broadband linear power amplifier as a first stage driver (45 dB gain, 150 kHz - 250 MHz power bandwidth). The output of this amplifier is routed to the 10X229 8-110 MHz Power Amplifier for final stage gain. Optionally, 2N5643 common-emitter Rf amplifiers may be used at 68 or 109 MHz to provide higher drive levels. The 10X229 amplifier again employs an Eimac 8877 power triode, but this time the output is tunable over the specified operating range.

Transmitter signal routing is via two 16X973 Rf switches connected in a cross-bar arrangement, controlled by 50 Ω TTL level control lines. Transmitter outputs are routed to transmit/receive (T/R) switches (duplexers) that are designed for specific operating frequencies. The 16X965 T/R switch is used for 27.36 MHz (^{15}N) and 270 MHz (^1H), and the 16X966 T/R switch is used for 67.90 MHz (^{13}C) and 109.30 MHz (^{31}P). These switches isolate the receivers from the probe and transmitter during Rf pulse transmission, and they isolate the probe and receiver

from transmitter noise during the balance of the operating cycle.

Design details of both T/R switch units are described below.

270 MHz (^1H) signals are also routed through a Hewlett-Packard 778D directional coupler. A low power transmitter may be attached to one of the coupled ports for homonuclear proton decoupling, solvent suppression by presaturation, or other experiments requiring concurrent use of both a high power and a low power Rf amplifier.

2.3.3 Observe Channel Receiver

The observe channel receiver is a superheterodyne receiver employing synchronous detection (Figure 2.5). Initial Rf preamplification and band limiting is done by amplifiers contained in the T/R switch assemblies, 16X965 (270 MHz and 27.36 MHz) and 16X966 (67.90 MHz and 109.30 MHz). Rf amplifier bandwidths and gains are summarized in Table 2.2. In all cases, the image frequency is far enough removed so that image noise should not be a problem.

The appropriate preamplifier output is routed to the down-converter via a step attenuator. Frequency conversion to the 30 MHz IF is done in the 10X210 8-270 MHz to 30 MHz Linear Converter using a standard output level HP 10514A double-balanced mixer running with a saturated local oscillator signal supplied by the 10X206 Single Sideband Generator. Varying degrees of signal attenuation may be required in the case of strong signals (^1H) to minimize intermodulation distortion that will increase drastically as the signal level approaches the 1 dB compression point. Similarly, as signals approach 2 to 5 dB below the mixer's conversion compression level, Rf desensitization by large signals

Figure 2.5 8 - 270 MHz Spectrometer System Observe Channel
Receiver Interconnect [16X9673-B4]

Table 2.2 Observe Channel Receiver Characteristics

<i>NUCLEUS</i>	<i>RECEIVER CENTER FREQUENCY</i>	<i>RECEIVER RF BANDWIDTH</i>	<i>RF AMPLIFIER GAIN</i>
¹ H	269.9985 MHz	40 MHz	36 dB
² H	41.4513 MHz	3 MHz	18 dB
¹³ C	70.0000 MHz	30 MHz	30 dB
³¹ P	120.0000 MHz	40 MHz	30 dB
¹⁵ N	27.3600 MHz	3 MHz	36 dB

may cause increased small signal conversion loss.

A saturating local oscillator drive level is used to drive the mixer diodes into the linear regions of the I-V curve and reduce the percentage of time spent in the nonlinear region about $V = 0$. Such a mode of operation improves odd order harmonic suppression and two-tone third order suppression. Use of a higher level mixer that incorporates series diodes could in principle further improve the receiver dynamic range, but at the expense of requiring a higher local oscillator drive level.

After conversion to the 30 MHz IF, the signal is amplified by up to 80 dB and is pre-detection filtered. The IF amplifier output is routed to another step attenuator so that the output level to the synchronous detector can be adjusted, again to minimize odd order harmonic distortion and effects caused by excessive signal levels.

Synchronous detection is performed in the 10X211 Double Phase Sensitive Detector using a 30 MHz reference from the 10X223 30 MHz Active Power Divider. The reference 30 MHz is split into quadrature components with a quadrature hybrid, and the incoming signal is phase detected twice using these two references. The two resulting signals (in phase quadrature) are low-pass filtered and amplified by up to 40 dB by the 16X902 Remote Controlled Filter. After this anti-aliasing filter (4-pole Butterworth characteristic), the two signals are digitized and signal averaged by a Nicolet Instruments Corp. Model 1180 Data System.

2.3.4 Decouple Channel Transmitter


The decouple channel transmitter is structured in the same fashion as the observe channel transmitter. Two interconnection drawings show the configurations for homonuclear ^1H decoupling (Figure 2.6) and heteronuclear noise-modulated proton decoupling for ^{13}C , ^{31}P , or ^{15}N observation. (Figure 2.7).


Efficient homonuclear proton decoupling utilizes a pulsed Rf irradiation at the decoupling frequency that is gated at the sampling rate of the digitizer (of the data acquisition system) as set by the acquisition spectral width (software set). The timing of this pulse train is synchronized with the sampling window of the sample and held on the digitizer input so that decoupling pulses occur only during the null time between samples. Pulse width and delay relative to the sampling window are set using a pulse generator (Hewlett-Packard 8013A or equivalent) triggered by the D CLOCK signal (dwell clock) available on the 1180 data system rear panel after suitable modification to the sweep control board (board 13) of the processor. Sweep control board modifications are summarized on Print 16X983.

At present, there are two Hewlett-Packard 5105A frequency synthesizers in the spectrometer system. One is normally used to supply the observe channel local oscillator, and the other is used for the lock channel local oscillator. The need for a third synthesizer for homonuclear decoupling is circumvented by using the observe channel local oscillator to also drive the decoupler. This requires that decoupling always take place at the center of a spectrum when quadrature phase detection is employed. This is generally not a problem for proton

Figure 2.6 8 - 270 MHz Spectrometer System Decouple Channel
Transmitter Interconnect for Homonuclear
Decoupling (^1H) [16X9673-B7]

Figure 2.7 8 - 270 MHz Spectrometer System Decouple Channel
Transmitter Interconnect for Heteronuclear
Proton Decoupling [16X9673-B8]





spectra when a wide spectral width and enough data points can be used.

The observe channel local oscillator is obtained from the "pass-through" 38 - 300 MHz LO output of the observe channel 10X206 Single Sideband Generator. The 30 MHz IF is gated by the decouple channel 10X207 Four Phase Generator, and it is then passed through a 0 - 10.9 dB (0.1 dB step) attenuator for fine power regulation. After single sideband conversion in the 10X206 Single Sideband Generator, the gated signal is passed through a 16X980 high isolation Rf switch (90 dB on/off ratio) prior to amplification by a Hewlett-Packard/Boonton Radio 230A tuned Rf power amplifier.

The low power homodecoupling Rf pulse train is duplexed with the high power observe channel pulses using a Hewlett-Packard 778D Dual Directional Coupler. This coupler isolates the output stages of the two transmitters, preventing adverse interactions. The 778D coupler has produced excellent results for proton homodecoupling, yielding spectra that are surprisingly free of spurious lines.

The heteronuclear proton decoupling connection shown in Figure 2.7 uses primarily the 10X221 30 - 300 MHz Frequency Multiplier as a local oscillator source. This requires that the magnetic field, as adjusted by the coarse H_0 field sweep amplifier (part of the 16X935 63 kG Magnet Shim Control) be adjusted so that the center of the decoupling bandwidth will occur in the desired part of the proton spectrum. The observe and lock local oscillator frequencies are then adjusted accordingly. Alternatively, the lock channel may be disabled and the lock synthesizer used to supply the decoupling center frequency.

For noise-modulated proton decoupling, the 30 MHz IF is 0/180°

phase modulated using a Hewlett-Packard 10514A double-balanced mixer as a modulator. A Hewlett-Packard 3722A shift register type pseudo-random noise generator is used to provide a binary pulse string to drive the modulator. Modulation bandwidth is adjusted by varying the clock period; the radiated power spectral density may be checked with a spectrum analyzer operating at either the 30 MHz intermediate frequency (looking at the output immediately after the modulator) or at the 270 MHz decoupling frequency (looking at the signal before the 10X225 270 MHz Power Amplifier or at the probe connection using a dummy load and -40 dB pick-off). Again the 0.1 dB step attenuator inserted into the 30 MHz line after the 10X209 Four Phase Generator is used for fine power control. Decoupling power is monitored with a Bird Thru-Line Rf directional wattmeter with either a standard plug-in element or a -40 dB directional coupler element. Decoupler power levels of < 10 watts RMS are normally sufficient for high resolution heteronuclear proton decoupling experiments.

2.3.5 Lock Channel

Although superconducting magnets have extremely stable fields, it has been found that some form of internal locking (*i.e.* H_0 field stabilization based on a resonance internal to the sample) is of great utility. Any nucleus may be selected for locking; deuterium is the choice in most high resolution work because of the pervasive use of deuterated solvents in NMR sample preparation. The principal benefits of operating in an internal locked mode are:

- Minimization of field drift during long term signal averaging


that is required for experiments such as 2D data acquisition.

- A combination of monitoring the lock channel signal amplitude and viewing the lock signal lineshape and ringing pattern provides a convenient method of optimizing the H_0 field homogeneity.

The lock channel block diagram is shown in Figure 2.8. The channel may be subdivided into several functional units:

1. Lock channel transmitter -- operating under gating control of the 16X894 Lock Channel Control. The lock channel transmitter is very similar to the observe channel transmitter (2.3.2) except that the output power level is far lower. Frequency generation is again via single sideband mixing with a synthesizer supplied local oscillator.
2. Lock channel receiver -- is also similar to the observe channel receiver. After conversion to the 30 MHz intermediate frequency and predetection filtering, the lock signal is phase detected and supplied to the 16X894 control unit.
3. Lock channel control -- detects the lock signal and generates an error signal that is applied to the H_0 field sweep amplifier that is incorporated in the room temperature shim system controller (16X935).
4. Magnet shim control -- contains the H_0 field sweep amplifier to which error signals are applied. This unit also handles the display functions associated with field homogeneity adjustments.

Figure 2.8 8 - 270 MHz Spectrometer System ^2H Lock System
Interconnect [16X9673-B6]



The lock channel operates in the pulsed time-shared mode, with a gating frequency of 1.25 KHz. This implies that the lock nucleus is resonated at 1.25 KHz below the actual Larmor frequency, using the lower modulation sideband. The 1.25 KHz is derived in the 16X894 Lock Channel Control from a 3.0 MHz reference (from the HP 5110A frequency synthesizer driver) by direct digital division. Quadrature 5 KHz reference frequencies are also derived for use later in detection.

The 1.25 KHz pulse train (with adjustable delay and duty cycle) gates the 16X919 Lock Channel Single Sideband Generator. This unit is functionally equivalent to the 10X206 Single Sideband Generator used in the observe and decouple channels except that a phase-amplitude trim network is set to optimize conversion (harmonic suppression) at an output operating frequency of 41.451 MHz (the normal deuterium frequency).

Transmit/receive switching is handled by the 16X964 T/R Switch assembly. This unit differs from the T/R switches used for observation because of the difference in power levels. A 3 dB resistive pad is used for transmitter isolation; most of the noise blanking is performed by the final high isolation Rf switch in the 16X919 SSB generator. Diode switches were used originally in the T/R switch, but they were found to result in extreme non-linear behavior of the signal DC offset when the transmitter power level was varied.

The lock receiver is identical to the observe channel receiver; a 16X978 41.45 MHz amplifier is cascaded before the linear converter to provide 0 - 80 dB (adjustable) additional Rf gain. The 10X211 phase detector output is a low-pass signal whose amplitude is proportional to the lock signal strength.

The low-pass lock signal is bandpass filtered in the 16X894 lock control unit, and then it is phase detected at 5 KHz. Detection is done at four times the transmitter gating frequency to increase the efficiency of transmitter blanking. Because of a large amount of transmitter leakage, detection at 1.25 KHz would lead to large DC offsets in the detected output.⁴ Detection at four times the gating frequency has the added advantage that both absorption and dispersion mode signals may be conveniently generated; the dispersion signal is used to generate the actual control signal and the absorption signal is used for display purposes.

Proportional plus integral control action is used. The control signal that is fed to the H_0 field sweep amplifier, $m(t)$, is generated from the detected dispersion signal, $e(t)$, in a form such as

$$m(t) = K_p e(t) + \frac{K_p}{T_i} \int_0^t e(t) dt \quad (2)$$

K_p , the proportional sensitivity, is the controller gain; it may be adjusted at a number of points in the control loop. T_i is the integral time, it may also be internally adjusted. Proportional plus integral control action is used because pure porportional control leads to steady-state errors which may be eliminated by addition of the integral term.

2.3.6 Probes

Probes are key elements of the 8 - 270 MHz spectrometer system. Their function of course is to couple radiofrequency excitations to the

sample and detect the resulting behavior of the nuclear magnetizations.

Most high resolution probes provide at least the following capabilities:

- Observation at either a single frequency or over a set bandwidth.
- Proton decoupling for the observation of ^{13}C , ^{31}P , ^{15}N , etc.
- Provision for heteronuclear internal locking, usually on deuterium.
- Monitoring of sample temperature, with provision for heating or cooling.

The observation and locking functions are usually fulfilled by one coil, and decoupling is done with a second orthogonal coil. The solenoid configuration of a superconducting magnet imposes another constraint: if probes are to remain in place when samples are interchanged, the samples must be removed axially. Since the Rf coils must generate fields orthogonal to the H_0 field direction, this dictates the use of Helmholtz style coils rather than solenoids.

Factors influencing probe design have been discussed extensively in the literature.⁵⁻⁹ The typical arrangement consists of a parallel resonant circuit with elements R , L , and C . The sample is contained in the coil, L , and the capacitance C is used to set the resonance frequency:

$$\omega_0 = 1 / \sqrt{LC} \quad (3)$$

Losses are modeled by the parallel resistance R , and the figure of merit is defined by

$$Q = R / \omega_0 L = R \omega_0 C = R \sqrt{C/L} \quad (4)$$

To maximize power transfer from the transmitter to the resonant circuit (and conversely from the sample to the receiver), a passive (loss-less) two-port network may be used to effect an impedance transformation (Figure 2.9). Synthesis of these two port networks and general probe design principles have been described elsewhere.⁹⁻¹³

Five probes presently exist for the LCB 8 - 270 MHz spectrometer system. Associated sample tube sizes and operating frequencies are listed in Table 2.3. Circuit configurations are summarized in Figures 2.10 - 2.18. The 5 mm ^1H , 10 mm ^{13}C , 10 mm ^{31}P , and 15 mm ^{15}N probes represent reconstructions of commercial Bruker Instruments, Inc. probes. Reconstruction generally involved replacement of one or both coils, replacement of all tunable capacitors, and the repositioning of several key capacitors to make remote tuning possible. The 10 mm $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ probe was fabricated using Bruker probe housings; the balance of the probe consisted of locally fabricated parts. Remote tuning capabilities were included for all frequencies. Figures 2.19 and 2.20 show two views of this probe.


As is apparent from the drawings, tapped parallel tuned circuits (or a variation) were used in all cases. It was found that the parallel tuning capacitor could be used to adjust the resonant frequency with minimal changes to the probe input impedance, therefore when remote tuning capabilities were included, this was the capacitor that was made adjustable.

The relationship of coil geometry to performance has been discussed in detail by Hoult and Richards.¹⁴ Helmholtz coils of a geometry prescribed by Lyddane and Ruark,¹⁵ and others,^{16,17} were used for low frequencies (< 100 MHz). Above 100 MHz, coils of a design similar


Table 2.3 8 - 270 MHz Spectrometer System: Sample Probe Specifications

<i>PROBE</i>	<i>SAMPLE TUBE SIZE</i>	<i>OPERATING FREQUENCIES</i>		<i>REMOTE TUNING CAPABILITY</i>
^1H	5 mm	^1H	269.9985 MHz	
		^2H	41.4513 MHz	
^{13}C	10 mm	^{13}C	67.8975 MHz	
		^1H	269.9985 MHz	×
		^2H	41.4513 MHz	
^{31}P	10 mm	^{31}P	109.2970 MHz	
		^1H	269.9985 MHz	×
		^2H	41.4513 MHz	
^{15}N	15 mm	^{15}N	27.3600 MHz	×
		^1H	269.9985 MHz	×
		^2H	41.4513 MHz	
$^{13}\text{C}/^1\text{H}/^3\text{H}$	10 mm	^{13}C	67.8975 MHz	×
		^1H	269.9985 MHz	Wide range tuning
		^2H	41.4513 MHz	×


Figure 2.9 Matching of a Parallel Resonant Circuit with a
Passive Lossless Two-port




- Figure 2.10 5 mm $^1\text{H}/^2\text{H}$ Probe: 270 MHz Port and 41.45 MHz Port
Circuit Details
- Figure 2.11 10 mm ^{13}C Probe: 270 MHz Port Circuit Detail
- Figure 2.12 10 mm ^{13}C Probe: 67.90 MHz Port and 41.45 MHz Port
Circuit Details
- Figure 2.13 10 mm ^{31}P Probe: 270 MHz Port Circuit Detail
- Figure 2.14 10 mm ^{31}P Probe: 109.30 MHz Port and 41.45 MHz Port
Circuit Details
- Figure 2.15 15 mm ^{15}N Probe: 270 MHz Port Circuit Detail
- Figure 2.16 15 mm ^{15}N Probe: 27.36 MHz Port and 41.45 MHz Port
Circuit Details
- Figure 2.17 10 mm $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ Probe: 270 - 286 MHz Port Circuit
Detail
- Figure 2.18 10 mm $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ Probe: 67.90 MHz Port and 41.45 MHz
Port Circuit Details




All variable capacitors: Johanson 5761 paralleled with appropriate value ATC 100B series porcelain chip capacitors




All variable capacitors: Johanson 5761 paralleled with appropriate value ATC 100B series porcelain chip capacitors




All variable capacitors: Johanson 5761 paralleled with appropriate value ATC 100B series porcelain chip capacitors




All variable capacitors: Johanson 5761 paralleled with appropriate value ATC 100B series porcelain chip capacitor




All variable capacitors: Johanson 5761 paralleled with appropriate value ATC 100B series porcelain chip capacitors




All variable capacitors: Johanson 5761 or equivalent, paralleled with appropriate value ATC100B series porcelain chip capacitors



All variable capacitors: Johanson 5761 or equivalent, paralleled with appropriate value ATC 100B series porcelain chip capacitors




All variable capacitors: Johanson 5341 0.8 - 10 pF variable, paralleled with appropriate value ATC 100B series porcelain chip capacitor



All variable capacitors: Johanson 5341 0.8 - 10 pF variable, paralleled with appropriate value ATC 100B series porcelain chip capacitors

Figure 2.19 $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ 10 mm Probe, Insert Detail and 67.89 MHz/
41.45 MHz Matching Network Detail (Cover removed)

Figure 2.20 $^1\text{H}/^3\text{H}/^{13}\text{C}/^2\text{H}$ 10 mm Probe, Insert Detail and 270 MHz
Matching Network Detail (Cover removed)







Figure 2.21 Typical High Frequency Coil Cutting Pattern,
Dimensions Shown for Placement on the Inside
of a 13 mm O.D. Insert





to one of Dadok were used (visible in Figures 2.19 and 2.20). These were cut from 0.001" OFHC copper foil (Electronic Space Products, Los Angeles, CA) in a pattern such as is shown in Figure 2.21. The principal advantages of this configuration are:

1. Once the coil is soldered together, it is self-supporting and no adhesives are required to fasten it to the insert glass.
2. Reduced sensitivity to bulk susceptibility changes in the sample were noticed; the reason is not known.
3. Vastly improved heat dissipation when the coils were used for broadband noise-modulated decoupling, because of increased surface area and a consequent increase in heat transfer efficiency to the cooling gas stream.
4. Lower distributed inductance made the attainment of higher operating frequencies easier; also higher Q's were attained.

The most significant problem with this design coil was that the large surface area resulted in significant capacitive coupling to the other coil (if one was present). Probe design, construction, and tuning were done with the aid of the three network analysis methods described in Section 2.4.1.

2.3.7 Data Processing and Control

A Nicolet Instruments Corp. Model 1180 Data Processor is one of the key elements of the LCB 8 - 270 MHz spectrometer system. Operating under a Nicolet Technology Corp. supplied software package (NTCFT-1180),

this unit handles all data acquisition and control functions. This section will outline the digital interconnections between peripheral devices and the processor. Connections may be subdivided into four categories:

1. High speed mass I/O
2. Low speed I/O and display drivers
3. Direct I/O bus connections
4. Nicolet 293A Programmable Pulser I/O connections


All digital interconnections are shown in Figure 2.22. The first two categories are covered in detail in the Nicolet 1180 Service Manual,¹⁸ and only brief mention of connections is made. The third category includes devices connected via the 16X904 I/O Bus Intercept; among these are the 16X950 Temperature Control, the 16X956 CalComp 565 Interface, and the 16X902 Remote Controlled Low-Pass Filter. Devices in the last category include the 16X892 System Control Logic, and devices driven by the 293A digital-to-analog converters (DACs).

.1 High Speed Mass I/O

Two types of disk memory systems may be used with the 1180 -- a dual platter Diablo Model 44 drive using IBM 5440 type media, or a single platter (removable) Diablo Model 31 drive with IBM 2315 type media. The LCB system uses the latter; both drives share a common interface board with only minor timing variations.

When using Diablo 31 drives, bits 14 and 13 of the sector register specify the drive unit number. Up to four drives may be connected in daisy-chain fashion, using ribbon cables with integral ground planes.

Figure 2.22 8 - 270 MHz Spectrometer System Digital Interconnect
 [16X9673-B1]



Each drive has jumper pins on its motherboard that select the unit number that it will recognize, therefore drives may be chained in any order.

Only one of the drives was originally delivered with the 1180; this was assigned Unit 0. Drives 1 and 2 were obtained used; after reconditioning, they were installed with locally modified power supplies (described by print 16X9673-B2).

.2 Low Speed I/O and Display Drivers

The 1180 Slow I/O board is used to interface low speed serial data format peripherals such as terminals. Serial I/O is handled by universal asynchronous receiver/transmitter chips (UARTs); two are provided (Channels A and B). The Channel A UART was connected to a Digital Equipment Corp. LA36 Decwriter II terminal, while the channel B UART was replaced with a high speed IM6402CPD (Harris 6402) UART which permitted operation at 38.4 KBaud. Channel B was connected to either a VAX 11/780 or to a Nicolet NIC-80 processor. Discussion of these connections is in Appendix 1.

The other type of slow I/O goes via the Display Board. Display board DACs are used to drive an X-Y-Z display (Hewlett-Packard 1222A) for the display of memory contents (spectra, etc.).

.3 Direct I/O Bus Connections

Locally designed 1180 peripherals that require 20 bit parallel I/O are connected directly to the 1180 I/O Bus by way of the 16X904 I/O Bus Intercept. This implies that the device must be assigned a unique address (that is not used internally by the 1180), and access

to the device is via newly defined I/O instructions that incorporate this address. The three devices presently attached to the bus in this fashion were assigned the same device addresses that Nicolet Technology Corp. assigns to similar peripherals, thus minimal software changes were necessary. Details of each device will be provided in Section 2.5.

.4 Nicolet 293A I/O Connections

The bulk of the 293A I/O connections are for the 16X892 System Control Logic. The 16X892 print set and associated descriptions in Section 2.5 cover this subject thoroughly.

The other type of device connected to the 293A utilizes the DAC board in the 293A. These devices and their functional relationships are:

1. 16X950 Temperature Control -- uses DAC2 as part of the control loop; DAC2 drives the programmable power supply that powers the variable temperature gas heater.
2. 16X957 293A DX Extender -- provides for multiple connections to the 293A DX connector.
3. Hewlett-Packard 7004B recorder -- is used for generating plots when the CalComp 565 Incremental Plotter is not connected. DAC0 (DX-2) drives the X-axis, and DAC1 (DX-3) drives the Y-axis.

2.4 Test and Characterization Procedures

2.4.1 Complex Impedance, Reflection Coefficients, and VSWR

Sample probes for the 8 - 270 MHz spectrometer system were characterized by measuring the complex impedance of the various ports as a function of frequency. Three overlapping techniques were used; applicability was determined by the frequency range of interest, and by whether a swept frequency display was preferred or if it was necessary to get complete phase information at a single frequency.

When a transmission line is terminated with a load that is not identical to the characteristic line impedance, a standing wave pattern is set up along that line. This pattern is composed of an incident voltage wave E_i and a reflected wave E_r . The ratio of these two is the reflection coefficient, $\bar{\rho}$, which is a vector quantity that includes the relative phase between the two (Figure 2.23)

$$\bar{\rho} = \rho \angle \theta_{\rho} = \frac{E_r}{E_i} \quad (5)$$

The angle θ_{ρ} varies as the distance along the line from the load. The reflection coefficient is determined by the impedance of the load (Z) and the characteristic impedance of the line (Z_0):

$$\frac{Z}{Z_0} = \frac{1 + \bar{\rho}}{1 - \bar{\rho}} \quad (6)$$

The voltage standing wave ratio (VSWR) and return loss don't carry

phase information:

$$\text{VSWR} = \frac{1 + |\bar{\rho}|}{1 - |\bar{\rho}|} \quad (7)$$

$$\text{Return Loss} = -20 \log |\bar{\rho}| \quad (8)$$

Reference points are an open circuit ($\rho = 1 \angle 0^\circ$) or a short circuit ($\rho = 1 \angle 180^\circ$).

2.4.2 Measurement Techniques

A. Complex Impedance, 500 KHz to 109 MHz

Impedance measurements in this frequency range were made with a Hewlett-Packard 4815A Vector Impedance Meter. This device measures the Rf current flowing through a device for a given low excitation voltage, and converts this information to a complex impedance reading (magnitude and phase angle).

B. Complex Impedance, 100 MHz to 1 GHz

Impedance measurements in this frequency range were made using a Hewlett-Packard 8405A Vector Voltmeter in conjunction with a Hewlett-Packard 778D Dual Directional Coupler in the test arrangement of Figure 2.24. The directional coupler is used to provide accurate samples of the incident and reflected wave amplitudes, and the phase angle between the two. Prior to making measurements, a short ($\rho = 1 \angle 180^\circ$) is placed at the end of the coupler, and a phase reference of 180° is set on the vector voltmeter using the phase offset


and by stretching the sliding line.

The measured reflection coefficient is readily converted to complex impedance using a Smith chart. Alternatively, a short BASIC program was written to perform this conversion (Appendix 2).


C. VSWR, 1 MHz to 1 GHz

Voltage standing wave ratio measurements were made in a swept frequency mode using a Wiltron Model 610B swept frequency generator and an Anzac RB-3-50 VSWR Bridge. The test setup is illustrated in Figure 2.25. Since the VSWR does not contain phase information, this method was most useful when the network was already fairly well characterized; it was extremely useful for final probe tuning adjustments where the impedance match was already quite close.

- Figure 2.23 Angle θ and Magnitude E_r/E_i of Incident and Reflected Waves as Determined by Load Impedance
- Figure 2.24 Complex Impedance Measurement Using HP 8405A Vector Voltmeter and HP 778D Dual Directional Coupler
- Figure 2.25 Swept Frequency VSWR Measurement Using Wiltron 610B Frequency Sweeper and VSWR Bridge



$$\frac{E_r}{E_i} = \overline{\rho}$$




COMPLEX IMPEDANCE MEASUREMENT TEST SET-UP USING DIRECTIONAL COUPLER

XBL 7910-12483

WILTRON 610B SWEPT FREQUENCY GENERATOR

X-Y DISPLAY SCOPE



IMPEDANCE MISMATCH
MEASUREMENTS USING
VSWR BRIDGE

XBL 7910-12481

2.4.3 Evaluation of System Noise Performance

The usual form of the signal-to-noise ratio after a 90° pulse is given by Abragam:¹⁹

$$\frac{v_s}{v_n} = K\eta M_0 \left(\mu_0 Q \omega_0 V_c / 4FkT_c \Delta f \right)^{1/2} \quad (9)$$

where K is a numerical factor (≈ 1) related to coil geometry; η is the "filling factor," *i.e.* the fraction of the coil volume occupied by sample; V_c is the coil volume; M_0 is the nuclear magnetization; Q is the quality factor of the coil; ω_0 is the Larmor frequency; μ_0 is the permeability of free space; F is the receiver noise figure at bandwidth Δf (Hz); k is Boltzmann's constant; and T_c is the probe temperature. It has been pointed out by several authors that this is not a fundamental equation, but is useful for order of magnitude calculations.^{14,19,20}

Of the four unknowns in equation (9), K , η , Q , and F , only the latter two are easily measurable. This section is devoted to discussion of the last factor -- receiver noise performance, and other system noise contributions.

One of the traditional measures of spectrometer system performance has been the one-shot signal-to-noise ratios for certain standard samples. The difficulty with this approach is that it does not discriminate between the effects of magnetic field inhomogeneity and receiver noise performance. These effects were separated during the construction of the 8 - 270 MHz spectrometer system by making separate noise figure measurements on receiver system components, both individually and in cascade.

2.4.4 Noise Sources

Noise at an NMR receiver input may be attributed to two primary sources:

- Thermal noise in the probe, and radiative pick-up
- Transmitter noise

The probe is an RLC one-port network, therefore the mean square noise voltage at its output is given by Nyquist's formula:

$$\langle v_n^2(t) \rangle = 2 kT_c \int_{-\infty}^{\infty} R(f) df \quad (10)$$

where $R(f)$ is the real part of the complex impedance seen looking into the port. Thus for a narrow-band tuned probe, the principal contribution is thermal noise.

The nature of the output stages of the Eimac 8877 power triode amplifiers used for the transmitters suggests that transmitter noise is a combination of hot thermal noise and shot noise. The normal isolation procedure is to minimize current flow due to this noise voltage by using series crossed diodes. Since the current flow is related to the voltage drop V across each diode by

$$I = I_s \left(\exp \left(\frac{eV}{kT} \right) - 1 \right) \quad (11)$$

where I_s is the reverse saturation current, enough diodes are used in series so that the drop across each is relatively small and I is accordingly small. Further discussion is included with the description

of the transmit/receive switches (16X964, 16X965, 16X966).

2.4.5 Noise Figures

The most useful measure of a device's noise performance is the noise figure, F . For the ℓ th device in the cascade of Figure 2.26, the noise figure F_ℓ is

$$\left(\frac{S}{N} \right)_\ell = \frac{1}{F_\ell} \left(\frac{S}{N} \right)_{\ell-1} \quad (12)$$

or in terms of decibels:

$$F_{\text{dB}} = 10 \log_{10} F_{\text{ratio}} \quad (13)$$


For an ideal device that introduces no noise, $F_\ell = 1$ (0 dB). It can be shown that if the available internal noise power of the device is P_{int} ,

$$F_\ell = 1 + \frac{P_{\text{int}}}{G_a k T_s B} \quad (14)$$

where G_a is the gain of the device, T_s the internal noise source temperature, and B the input noise source bandwidth. Thus increasing the gain tends to diminish the importance of the internal noise contribution, driving the F_ℓ towards unity.

The noise figure of a cascade of devices is given by Friis' formula:

Figure 2.26 Cascade of Subsystems Making Up a System, with Signal-to-Noise Ratios Defined at Each Point



$$F = F_1 + \frac{F_2 - 1}{G_{a_1}} + \frac{F_3 - 1}{G_{a_1} G_{a_2}} + \dots \quad (15)$$

where F_i is the noise figure of stage i and G_{a_i} is the corresponding gain of that stage. Passive devices such as mixers have insertion loss, *i.e.* a fractional gain < 1 .

Inspection of the superheterodyne receiver structure of Figure 2.3 illustrates some of the important concerns with regard to the design of the observe channel receiver (Figure 2.5) for optimum noise performance. The Rf preamplifier noise figure (F_1 in this instance) will be the dominant term in the system noise figure, and the higher its gain G_{a_1} , the less important the noise figure of the mixer (which is typically 9 dB for the Hewlett-Packard 10514A units used in the 10X210 8 - 270 MHz to 30 MHz Linear Converters). The IF amplifier noise figure F_3 is still a crucial factor because G_{a_2} represents a fractional gain (HP 10514A single sideband conversion loss is typically 9 dB).

2.4.6 Noise Figure Measurement Procedure

Noise figure measurements were made with either an Ailtech Noise Figure Meter or a Hewlett-Packard 160 Noise Figure Meter. Both are equipped with broadband calibrated noise sources that are attached to the input of the receiver under test.

The 30 MHz IF output of the 10X210 8 - 270 MHz to 30 MHz Linear Converter is returned to the 30 MHz input on the noise figure meter,

and a direct reading of cascaded system noise figure may be made.

Receiver absolute performance was measured with the transmitter physically disconnected. Transmitter noise isolation was then measured by reconnecting the powered on/gated off transmitter and observing the degradation in system noise figure (if any).

Table 2.4 Observe Channel System Noise Performance

<i>NUCLEUS</i>	<i>LOCAL OSCILLATOR FREQUENCY</i>	<i>IF AMPLIFIER GAIN</i>	<i>NOISE FIGURE</i>
¹ H	299.9985000 MHz	60 dB	2.3 dB
³¹ P	139.2970000 MHz	60 dB	2.2 dB
¹³ C	97.8975000 MHz	80 dB	2.3 dB
² H	71.4513500 MHz	60 dB	5.0 dB
¹⁵ N	57.3600000 MHz	80 dB	1.5 dB

2.5 System Components

This section provides design summaries and operating characteristics of major components of the LCB 8 - 270 MHz spectrometer system. Only significant original designs are included.

2.5.1 16X892 System Control Logic

The LCB 8 - 270 MHz spectrometer system has as a key element a Nicolet Instruments Corp. Model 1180 Data System. This system includes a Model 293A Programmable Pulser, a unit that includes seven software settable timers with 32 nsec resolution, four digital-to-analog converters (DACs) with 12 bit resolution, programmable levels, and sense lines. When the 1180 system is used with a 293A and one or more disk memory subsystems, a large body of Nicolet supplied software that is applicable to NMR may be used. In particular, Nicolet Technology Corp. supplies a software package called NTCFT-1180, which will provide all of the control signals normally required to operate an NMR spectrometer system, in addition to handling data acquisition and reduction. Since the NTCFT software package is designed for use with NTC supplied spectrometer systems, there existed a large number of incompatibilities that had to be resolved before this package could be used on the LCB 8 - 270 MHz system.

The function of the 16X892 system control logic is to make the LCB 8 - 270 MHz system look like an NTC system. This approach was taken in preference to making software modifications for several reasons:

- Major software packages such as NTCFT tend to evolve over a period of years. The release of new software versions would necessitate local modification of each version as it is released.
- Software modification is a task that is perhaps an order of magnitude more difficult than designing compatible hardware.

Changing selected areas of large programs tends to be difficult, since it is hard to predict the global effects of these changes.

- After a large number of NTC systems are shipped, NTC will be obliged to insure compatibility of new software with existing hardware, therefore future hardware changes are likely to be minimal because of the cost of retrofitting.

Design Summary

The NTCFT -1180 software package running on the 1180/293A makes use of the 293A inputs and outputs that are listed in Table 2.5. The hardware implementation is shown in Figure 2.27.

The 8 - 270 MHz spectrometer system uses 10X209 Four-Phase Generators to produce phase shifts in the observing frequency. This unit has four gate control lines for 0° , 90° , 180° , and 270° ; these lines drive Rf switches on the four outputs of a Merrimac quadripole network. Since NTCFT and the 293A signal changes in the phase rather than producing gate control signals that are high for the duration of the desired phase, a simple four-state machine was used as an interface.

The state diagram for this machine is shown in Figure 2.28. State definitions are as follows:

<u>State (Q1,Q0)</u>	<u>Condition</u>
0 0	0° phase shift on
0 1	90° phase shift on
1 0	180° phase shift on
1 1	270° phase shift on

Table 2.5 Summary of 1180/293A I/O Functions

293A Timer Outputs

Output	Function	293A Connection(GND)
SP0	Observe Rf on	ZB-5 (ZB-18)
SP1	Programmable via NX command	ZB-2 (ZB-15)
SP2	Programmable via NX command	ZB-6 (ZB-19)
SP3	Programmable via NX command	ZB-3 (ZB-16)
SP4	Programmable via NX command	ZB-4 (ZB-17)
SP5	1180 ADC trigger	ZB-1 (zb-14)

293A Programmable Levels

Level	Function	293A Connection
LEV0	Decoded by 16X892 logic	YC-17
LEV1	Decoded by 16X892 logic	YC-16
LEV2	Decoded by 16X892 logic	YC-14
LEV3	Decoded by 16X892 logic	YC-15
LEV4	Decoupler off	YC-4
LEV5	Aux. logic trigger to 16X892	YC-3
LEV6	Step plotter left (unused)	YC-2
LEV7	Release plotter x (unused)	YC-21
LEV8	Delay flag	YC-7
LEV9	Decoupler CW	YC-6
LEV10	unused	YC-5
LEV11	Step plotter right	YC-8
LEV12	Tune mode on (unused)	YA-2
LEV13	Decoupler level 2 (unused)	YA-3
LEV14	Used by NTCCOR	YA-1
LEV15	Release plotter y (unused)	YA-4
PENLIFT	Lower pen	DX-7

Table 2.5 (Continued)

293A DACs		
DAC	Function	293A Connection(GND)
DAC0	Recorder X-axis	DX-2 (DX-15)
DAC1	Recorder Y-axis	DX-3 (DX-16)
DAC2	16X950 Temperature Control	DX-4 (DX-17)
DAC3	16X950 Temperature Control monitor	DX-5 (DX-18)

293A Sense Lines


Sense	Function	293A Connection
SENSE0	Plotter off left limit	YB-15
SENSE1	Plotter off right limit	YB-5
SENSE2	Lock on signal from 16X892	YB-20
SENSE3	External add/subtract flag	YB-3
SENSE4	Spinner rate (unused)	
SENSE5	Spinner rate (unused)	

Decoding of LEV0 - LEV3 by 16X892 logic

Levels (3210)	Function
0 0 0 0	Turn 90° phase shift off
0 0 0 1	Turn 90° phase shift on
0 0 1 0	Turn 180° phase shift off
0 0 1 1	Turn 180° phase shift on
0 1 0 0	No-op
1 0 1 0	Select FT mode
1 1 0 0	Low power FT mode


Figure 2.27 270 MHz Spectrometer System/Nicolet Instrument 293A
Interface

Figure 2.28 270 MHz Spectrometer System/Nicolet Instrument 293A
Interface, State Diagram and Transition Table



TRANSITION TABLE

PRESENT STATE		NEXT STATE		INPUTS								J-K IMPLEMENTATION							
Q_1	Q_0	Q_1	Q_0	J_1	J_0	K_1	K_0	L_1	L_0	M_1	M_0	J_1	J_0	K_1	K_0	L_1	L_0	M_1	M_0
0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	0	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	1	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	1	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
0	1	1	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	0	1	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	1	0	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	1	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0
1	1	1	1	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0



STATE DEFINITIONS:

00 = 0° PHASE
 01 = 90° PHASE
 10 = 180° PHASE
 11 = 270° PHASE

INPUT DEFINITIONS:

PRG LVL (PRG) → 3 2 1 0

0	0	0	0	90° P SHIFT OFF
0	0	0	1	90° P SHIFT ON
0	0	1	0	180° P SHIFT OFF
0	0	1	1	180° P SHIFT ON
1	1	1	0	DISABLE MANUAL ADJ CONTROL
1	1	1	1	ENABLE MANUAL ADJ CONTROL

170 Mhz Spectrometer System/Manual		UNIVERSITY OF CALIFORNIA	
Instrument 230 Infrared		UNIVERSITY OF CALIFORNIA	
State Diagram and Transition Table		UNIVERSITY OF CALIFORNIA	
Date: 10/1/83		UNIVERSITY OF CALIFORNIA	
Author: 16X8923-H1		UNIVERSITY OF CALIFORNIA	

The transition table (Table 2.6) shows programmable levels LEV3, LEV2, LEV1, and LEV0 as inputs, and LEV5 provides a trigger to clock the state transitions. The transition table shows an implementation with J-K flip-flops; Karnaugh maps and minimum implementations are shown in Figures 2.29 through 2.32.

Hardware implementation

A. Next State Generator (Board 2)

The next state generator board contains the two J-K flip-flops that hold the current phase state (7473s, Q0 and Q1), plus a third J-K flip-flop whose output is high when the manual pulse width control is enabled. Gating for manual pulse width control has not been installed. Combinatorial logic is implemented with 7420 NANDs, and a buffered LEV5 is used to clock transitions.

B. State Decoder (Board 3)

The current phase state is decoded by a 2-to-4 line decoder (74LS139); the decoder outputs are fed to the line driver board. All external control logic lines are buffered. The state decoder board also contains a comparator (LM311) that senses the presence or absence of LOCK from the 16X894 Lock Channel Control. This is used for a software interlock in NTCFT-1180.

C. Line Driver Board (Board 4)

The line driver board uses 12 standard 50 ohm complementary emitter-follower line drivers. All control lines on the system are

Table 2.6 16X892 Transition Table

Present State		Next State		Inputs				J-K Implementation			
Q1	Q0	Q1	Q0	LEV3	LEV2	LEV1	LEV0	J1	K1	J0	K0
0	0	0	0	0	0	0	0	0	X	0	X
0	0	0	1	0	0	0	1	0	X	1	X
0	0	0	0	0	0	1	0	0	X	0	X
0	0	1	0	0	0	1	1	1	X	0	X
0	0	0	0	1	1	X	X	0	X	0	X
0	1	0	0	0	0	0	0	0	X	X	1
0	1	0	1	0	0	0	1	0	X	X	0
0	1	0	1	0	0	1	0	0	X	X	0
0	1	1	1	0	0	1	1	1	X	X	0
0	1	0	1	1	1	X	X	0	X	X	0
1	0	1	0	0	0	0	0	X	0	0	X
1	0	1	1	0	0	0	1	X	0	1	X
1	0	0	0	0	0	1	0	X	1	0	X
1	0	1	0	0	0	1	1	X	0	0	X
1	0	1	0	1	1	X	X	X	0	0	X
1	1	1	0	0	0	0	0	X	0	X	1
1	1	1	1	0	0	0	1	X	0	X	0
1	1	0	1	0	0	1	0	X	1	X	0
1	1	1	1	0	0	1	1	X	0	X	0
1	1	1	1	1	1	X	X	X	0	X	0

Figure 2.29 J1 Karnaugh Map

Q1 Q0

L1 L0

L3 L2 = 0 0

	00	01	11	10
00	0	0	X	X
01	0	0	X	X
11	1	1	X	X
10	0	0	X	X

Q1 Q0

L1 L0

L3 L2 = 1 0
0 1
1 1

	00	01	11	10
00	0	0	X	X
01	0	0	X	X
11	0	0	X	X
10	0	0	X	X

Minimum implementation: $J1 = L0 L1 \overline{L2} \overline{L3}$

Figure 2.30 K1 Karnaugh Map

Q1 Q0

L1 L0

L3 L2 = 0 0

	00	01	11	10
00	X	X	0	0
01	X	X	0	0
11	X	X	0	0
10	X	X	1	1

Q1 Q0

L1 L0

L3 L2 = 1 0
0 1
1 1

	00	01	11	10
00	X	X	0	0
01	X	X	0	0
11	X	X	0	0
10	X	X	0	0

Minimum implementation: $K1 = \overline{L0} L1 \overline{L2} \overline{L3}$

Figure 2.31 J0 Karnaugh Map

L1 L0		Q1 Q0			
		00	01	11	10
L3 L2 = 0 0	00	0	X	X	0
	01	1	X	X	1
	11	0	X	X	0
	10	0	X	X	0

L1 L0		Q1 Q0			
		00	01	11	10
L3 L2 = 1 0 0 1 1 1	00	0	X	X	0
	01	0	X	X	0
	11	0	X	X	0
	10	0	X	X	0

Minimum implementation: $J0 = L0 \overline{L1} \overline{L2} \overline{L3}$

Figure 2.32 K0 Karnaugh Map

Q1 Q0

L1 L0

L3 L2 = 0 0

	00	01	11	10
00	X	1	1	X
01	X	0	0	X
11	X	0	0	X
10	X	0	0	X

Q1 Q0

L1 L0

L3 L2 = 1 0
0 1
1 1

	00	01	11	10
00	X	0	0	X
01	X	0	0	X
11	X	0	0	X
10	X	0	0	X

Minimum implementation: $K0 = \overline{L0} \overline{L1} \overline{L2} \overline{L3}$

50 ohm coaxial cable, and all gate inputs are 50 ohm terminations. The use of low impedance transmission lines insures that long cable runs may be used without excessive noise pickup or reflection problems. The transistors used in the drivers are MHQ2222 or MHQ2907 quad NPN and PNP arrays.

2.5.2 16X894 Lock Channel Control Unit

The 16X894 Lock Channel Control provides the control signals and the lock-in detection necessary for the implementation of a pulsed time-shared heteronuclear lock system. The unit is used in conjunction with a Hewlett-Packard 5110A Frequency Synthesizer Driver (reference frequency source), and it provides control signals to the 16X919 Lock Channel SSB Generator, the 10X210 8-270 MHz to 30 MHz Linear Converter, and the 16X935 63 KG Magnet Shim Control.

The unit gates the observing transmitter at 1.25 KHz, and lock-in detection is performed at 5 KHz. This scheme permits dual detection of a pair of signals that are in phase quadrature, and allows the simultaneous generation of both an absorption mode signal (for display purposes) and a dispersion mode signal (for field correction).


Front panel controls permit adjustment of transmitter duty cycle, audio phase, field sweep rate, field sweep width, and lock threshold. The lock threshold detector is used to communicate the presence or absence of lock to the Nicolet 1180 data system via the 16X892 interface and the Nicolet 293A.


Design Summary

A. Frequency Generation and Control Signals (Figure 2.33)

The 1.25 KHz transmitter gating rate is derived from a buffered 3.0 MHz from the HP 5110A synthesizer driver. The 5110A is designed to drive one 50 ohm line per output, thus the 3.0 MHz input is buffered with an emitter-follower to increase the input impedance. This


Figure 2.33 270 MHz Spectrometer System, Lock Channel 5 KHz Four
Phase Generator (Board 1), Gate Generator
(Board 2)





ALL RESISTORS 1% METAL FILM
 CAPACITORS: SILVER MICA (300pF)
 OR CERO (0.1)
 ALL TRANSISTORS: CERAMET

270 MHz LOCAL OSCILLATOR	5 MHz SIGNAL INPUT	712	147K	300pF	147K	10K	100	1.0K	1.0M	2N4340	147K	300pF	147K	10K	100	1.0K	1.0M	2N4340	712	270 MHz LOCAL OSCILLATOR
LARSEN & TUBER UNIVERSITY OF CALIFORNIA 16X24X3-2																				



UNLESS INDICATED:
ALL FEEDINGS - 1/4 CUP
CARBOHYDRATE - CAGE CRACKING

[illegible]

buffered signal is converted into a TTL pulse train by the LM306 comparator, and the pulse train is divided by 150 to yield 20 KHz. The 20 KHz is routed to four successive T flip-flops to produce four square waves with relative phases of 0° , 90° , 180° , and 270° . The 0° signal is fed to two successive monostables (74121); the first is used to generate an adjustable time delay for audio phase shifting, and the second is used to vary the pulse width of the transmit gate for duty cycle control. The resultant pulse train is buffered by a complementary emitter-follower, and the signal is fed to the rear panel XMIT GATE for the 16X919 Lock Channel SSB Generator. The complementary output from the last monostable is also buffered and fed to the 10X210 Linear Converter for use as a receiver gate.

B. 5 KHz Signal Active Bandpass Filter (Figure 2.34)

The signal from the 10X211 Double Phase Sensitive Detector is fed via a rear panel BNC to board 3. On this board are two identical voltage controlled voltage source (VCVS) active filters with a center frequency of 5 KHz. The last stage is an AD741 voltage follower that provides buffering.

C. Dispersion Signal Phase Detector (Figure 2.35)

Both the dispersion and the absorption signal phase detectors are based on MC1496/1596 series monolithic balanced modulators. The MC1496 consists of a differential amplifier driving a dual differential amplifier; operation of the device is based on the ability to deliver the product of two input voltages when the magnitudes of these voltages are within the limits of linear operation.

Figure 2.34 270 MHz Spectrometer System, Lock Channel 5 KHz
Bandpass Filter (Board 3)

Figure 2.35 270 MHz Spectrometer System, Lock Channel Dispersion
Signal Phase Detector and Lock Control Signal
Processing (Board 4)

In this application, the dual differential amplifier is driven into saturation by the reference 5 KHz signal; the 5 KHz lock signal is fed via a trimpot to the differential amplifier, and the input amplitude is adjusted so that the device operates in the center of the linear region. The differential output is converted to a single-ended signal by a 747 differential amplifier, and the signal is then low-pass filtered (1-pole) with a 747-based active filter.

The control algorithm used for field control is the position form of the proportional plus integral control. The proportional term is taken directly from the low-pass filter output, and the integral term is generated with an AD503 FET-input operational amplifier connected as an integrator. The two terms are summed with an AD741 connected as a summing amplifier; trimpots are provided for the adjustment of the relative weightings of the proportional and integral terms. In addition, an offset adjustment pot is brought to the front panel so that DC offsets may be conveniently bucked.


D. Absorption Signal Phase Detector (Figure 2.36)

The absorption signal phase detector board is identical to the dispersion signal board, except there is no integrator and no summing amplifier. In addition, the low-pass filter bandwidth is larger (pole at a higher frequency) to make ringing easier to observe. Again, an offset adjustment is brought to the front panel to allow for the bucking of DC components.

E. Sweep Control Logic (Figure 2.37)

The sweep control logic board provides for the front panel


- Figure 2.36 270 MHz Spectrometer System, Lock Channel Absorption
Signal Phase Detector (Board 5)
- Figure 2.37 270 MHz Spectrometer System, Lock Channel Sweep Control
Logic (Board 7)
- Figure 2.38 270 MHz Spectrometer System, Lock Channel Analog Buffer
(Board 6) and Sweep Generator (Board 9)




UNLESS INDICATED:
ALL RESISTORS 5% 1/4 W CARBON
ALL CAPACITORS CMOS CERAMIC

[illegible]

XBL 7910-12492





selection of various lock channel modes of operation. Two 74LS253 multiplexers are used to select relative phases of the reference 5 KHz for the absorption and dispersion signal phase detectors. The inputs are the 0° , 90° , 180° , and 270° relative phase TTL pulse trains from board 1, and the outputs are fed to boards 4 (dispersion) and 5 (absorption).

The LM311 comparator is used to compare the absorption channel signal amplitude (when the unit is locked) to a variable (front panel adjustable) threshold. When the unit is unlocked and sweeping through a resonance, if the lock request line is grounded (by turning the front panel lock switch on), when the lock level sense comparator detects signal amplitude in the absorption channel that is above threshold, the field control relay is latched. After a locked state is attained, if the lock level falls below threshold, the field control relay is dropped. The field control relay is debounced with the NE555 timer and the 2N2222 inverter. In addition to driving the relay (using a 75451 peripheral driver), the logic also drives two 7440 high current NAND gates that are used to signal presence or absence of lock by driving front panel LEDs or external lines (to the Nicolet 1180 data system).

F. Analog Buffer and Sweep Generator (Figure 2.38)

Board 6 contains three LM310 voltage followers acting as buffers, and an LM741 operational amplifier connected as a low-pass filter and buffer. The LM310 buffers are connected to the absorption and dispersion signal phase detector outputs, and are present to insure adequate current drive capability. The low-pass filter is connected to the absorption signal; the output is fed to the threshold comparator on

board 7.

An LM 747 dual operational amplifier is used on board 9 as an adjustable amplitude sawtooth generator. This sawtooth drives the field sweep amplifier in the 16X935 63 KG Magnet Shim Control. A separate fixed amplitude ramp is also provided to drive the horizontal sweep of a display oscilloscope (thus the display sweep width is invariant with field sweep width).


2.5.3 16X902 Remote Controlled Low-Pass Filter



The 16X902 Remote Controlled Low-Pass Filter is used to provide 0 - 40 dB (switch selectable) audio gain between the 10X211 Double Phase Sensitive Detector and the 1180 data system analog-to-digital converter (ADC) inputs. The two filter channels are exact replicates (direct and quadrature channels), and have four-pole Butterworth characteristics. Turnover frequency may be selected either via remote programming via the 1180 I/O bus and the 16X904 I/O Bus Intercept, or by front panel switches.

In remote operation, the 1180 addresses the filter as device 322. An I/O instruction of generic form ACOUT (620000) with IOP2 set (623224) causes the lowest 10 bits of the accumulator to be loaded into a D register in the filter; this number sets the turnover frequency in 100s of Hertz.

The active filter networks are voltage-controlled voltage source realizations of low-pass networks. Turnover frequency is adjusted by switching the resistances with analog switches. The complete filter unit is described in Figures 2.39 - 2.43.

- Figure 2.39 Remote Controlled Low-Pass Filter, Analog Signal Processing [16X9023-S1]
- Figure 2.40 Remote Controlled Low-Pass Filter, Remote/Local Function Selection [16X9023-S2A]
- Figure 2.41 Remote Controlled Low-Pass Filter, Device Select and Remote Controller [16X9023-S3A]
- Figure 2.42 Remote Controlled Low-Pass Filter, Differential Amplifier Buffer [16X9023-S4]
- Figure 2.43 Remote Controlled Low-Pass Filter, LED Remote Frequency and Mode Indicator [16X9023-S5]



[illegible]

2.5.4 16X903 Low Frequency Rf Preamplifier

The 16X903 Low Frequency Rf Preamplifier is a low noise, moderate gain, FET cascode amplifier that is used at 27.36 MHz and at 41.45 MHz on the 8 - 270 MHz spectrometer system. The circuit is shown in Figure 2.44, with construction details illustrated by the photograph of Figure 2.45. The low level signal input is protected by crossed 1N4447 diodes, and a series resonant 270 MHz trap shunts any stray proton decoupling frequencies to ground.


Gain is provided by a matched pair Siliconix E420 dual J-FET. The identical drain saturation current I_{DSS} for the two FETs insures that neither will have a forward biased gate if direct interstage coupling is used. The input "el" matching network is set to transform the 50 ohm source impedance to the optimum impedance for the lowest noise figure. At 27.36 MHz, the optimum source impedance is about 1850 ohms; at 41.45 MHz the optimum is about 1700 ohms.

Typical measured performance at 27.36 MHz showed a gain of 18 dB, bandwidth of 3.5 MHz, and measured noise figure of 1.5 dB.

Figure 2.44 Low Frequency Rf Preamplifier [16X-9033]

Figure 2.45 Construction Details, Low Frequency Rf Preamplifier


ENCLOSURE 2 JPD BUREAU 333



* VALUES SUBJECT TO MEASUREMENT
ALL RESISTORS IN MEG OHM UNLESS NOTED
0.01, 0.1 CAPACITORS -- KAPOT CERAMIC TYPE

LOW FREQUENCY RF PREAMPLIFIER	
DATE	9/2/77
SCALE	1:1
NO.	16X9033
UNIVERSITY OF CALIFORNIA	
ENGINEER	W. L. SHIN
DATE	
REVISION	
CHANGES	

XBB 797-9250



2.5.5 16X904 Nicolet 1180 I/O Bus Intercept

A Nicolet Instruments Corp. Model 1180 Data System is used for data acquisition and control on the LCB 8 - 270 MHz spectrometer system. This processor has three major data busses: memory (M), data (D), and world (W). Only the world lines are available for connections outside of the main processor; and they function as part of a general purpose I/O bus.

The 1180 world lines are bidirectional three-state lines. Data may be transmitted or received by any device attached to the lines. Each I/O device that is attached to the world lines is identified by a unique device address. The processor may address any device; it does this by presenting the device address on address lines DV3 - DV10. The processor may request actions by that device by simultaneously presenting certain control signals: IOP0, IOP1, IOP2, TALK, LISTEN, or ACK.

Peripheral devices attached to the I/O bus may signal a request for service to the processor using the interrupt lines INT0- through INT6-. These lines are wired-or lines, that is all devices that output to the lines have open-collector outputs and the normal line state is held up by pull-up resistors. When a device wishes to signal the processor, it pulls down one of these lines (depending on the interrupt level to which that device is assigned). When the processor senses that an interrupt line is pulled down, it branches to a software interrupt service routine whose address is stored in the interrupt vector, and it deposits the return address at the first location of the service routine for an indirect return jump upon completion.

The three input/output pulses IOP0, IOP1, and IOP2, are usually used to clock the transfer of data to or from the 1180 accumulator to a peripheral device. All three pulses are 400 nsec long, and their occurrence is governed by bits 0 - 2 of the I/O instruction. Bits 3 - 10 of the I/O instruction set DV3 - DV10, the device address selection.

When the I/O bus leaves the 1180 chassis, connections are made with a 175 conductor ribbon cable. 86 signal lines are used; the remaining lines are interleaved grounds. Normally the 1180 is only connected to a 293A Programmable Pulser and a Nicolet I/O Tee. A bus termination is provided within the 293A. Each three-state line is terminated with a 220 ohm pull-up to + 5 VDC and 330 ohms to ground. Wired-or lines are terminated with a 4.7 Kohm pull-up to + 5 VDC.

Locally constructed devices that require connection to the 1180 I/O bus are connected via the 16X904 I/O Bus Intercept. This unit has two 86 conductor edge connectors for the 1180 I/O bus daisy-chain, and four 44 conductor edge connectors for use with Control Logic Corp. (CLC) 40 conductor interleaved ground cables. Only the signal lines that are commonly used are tied to the CLC connectors; this simplifies the construction of peripheral devices and reduces cabling costs. All world lines (W0 - W19), device lines (DV3 - DV10), IOPs, ACK, TALK, LISTEN, and INIT- lines are connected. Of the interrupt lines, only INT6- is connected. In general, locally constructed peripherals will only use interrupt level 6, as higher levels are reserved for the disk memory, the 293A, and internal processor functions such as clocks. Devices attached to the CLC bus are also daisy-chained, so each device must have two 44 conductor connectors on the rear panel. The last device in the chain must have a bus termination board in the extra slot.

Because the terminator board employs powered resistive terminations, the I/O peripheral must be capable of supplying at least 0.5 amps at + 5 VDC to the rear connectors. Since any device may be the last in the chain, all devices should be capable of supplying this current.

Signal definitions and bus conductor assignments are summarized in Table 2.7. Interface details for particular devices are contained in the documentation for the respective device.

Table 2.7 Nicolet 1180 I/O Bus Intercept Wire-Wrap Cross Reference

<i>SIGNAL DEFINITION</i>	<i>1180 I/O BUS ASSIGNMENT</i>	<i>1-43A/B BUS ASSIGNMENT</i>	<i>CLC I/O BUS ASSIGNMENT</i>
GND	1	1B	
GND	2	1A	
+5 VDC	3	2B	
+5 VDC	4	2A	
	5	3B	
DV10	6	3A	15
	7	4B	
DV9	8	4A	14
	9	5B	
DV6	10	5A	11
	11	6B	
DV8	12	6A	13
SKIP-	13	7B	W
W16	14	7A	L
DV3	15	8B	8
W17	16	8A	M
DV4	17	9B	9
W18	18	9A	N
DV5	19	10B	10
W19	20	10A	P
ACK	21	11B	U
W12	22	11A	F
IOP2	23	12B	17
W13	24	12A	H
IOP1	25	13B	S
W14	26	13A	J
IOP0	27	14B	R
W15	28	14A	K
DV7	29	15B	12
	30	15A	
	31	16B	
	32	16A	
	33	17B	
	34	17A	
INIT-	35	18B	A
	36	18A	
TALK	37	19B	T
	38	19A	
LISTEN	39	20B	16
W8	40	20A	B

Table 2.7 (continuation)

<i>SIGNAL DEFINITION</i>	<i>1180 I/O BUS ASSIGNMENT</i>	<i>1-43A/B BUS ASSIGNMENT</i>	<i>CLC I/O BUS ASSIGNMENT</i>
UPPER8-	41	21B	
W9	42	21A	C
SETNEG-	43	22B	
W10	44	22A	D
INT0	45	23B	
W11	46	23A	E
INT1-	47	24B	
	48	24A	
INT2-	49	25B	
	50	25A	
INT3-	51	26B	
	52	26A	
INT4-	53	27B	
W4	54	27A	6
INT5-	55	28B	
W5	56	28A	7
INT6-	57	29B	V
W6	58	29A	18
SPARE1	59	30B	
W7	60	30A	19
SPARE2	61	31B	
W0	62	31A	2
SPARE3	63	32B	
W1	64	32A	3
SPARE4	65	33B	
W2	66	33A	4
SPARE5	67	34B	
W3	68	34A	5
SPARE6	69	35B	
	70	35A	
SPARE7	71	36B	
	72	36A	
SPARE8	73	37B	
	74	37A	
SPARE9	75	38B	
	76	38A	
SPARE10	77	39B	
	78	39A	
SPARE11	79	40B	
	80	40A	

Table 2.7 (continuation)

<i>SIGNAL DEFINITION</i>	<i>1180 I/O BUS ASSIGNMENT</i>	<i>1-43A/B BUS ASSIGNMENT</i>	<i>CLC I/O BUS ASSIGNMENT</i>
<hr/>			
SPARE12	81	41B	
	82	41A	
	83	42B	
	84	42A	
GND	85	43B	20,22
GND	86	43A	X,Z

2.5.6 16X919 Lock Channel Single Sideband Generator

The 16X919 Lock Channel Single Sideband Generator performs frequency conversion of a synthesized local oscillator to the desired lock observation frequency. The unit also provides for lock channel transmitter gating and output level controls.

The 8 - 270 MHz spectrometer system uses a 30 MHz intermediate frequency, hence with lower sideband mixing the local oscillator frequency is the desired operating frequency + 30 MHz. Frequency conversion is again via a term-by-term realization of the trigonometric identity (1). Quadrature components of the IF and local oscillator frequencies are generated using quadrature hybrids. The products are then realized using double-balanced mixers, and the sum is formed with an in-phase power combiner.

Design Summary (Figure 2.46)

The IF and local oscillator frequencies are fed in via rear panel TNC connectors. Since a 30 MHz Active Power Divider (10X223) is available on the system, no pass-through 30 MHz output is provided. The local oscillator is fed to a reactive power divider (Anzac TU-50); one output is then returned to the rear panel and the other output is fed to a 3 dB pad followed by a Watkins-Johnson A7 hybrid wideband amplifier. The 3 dB pad adjusts the output level of the amplifier as it feeds the Anzac JH-14 40-400 MHz broadband quadrature hybrid. The attenuation value of this pad was found to give the maximal output level within the constraint of desired spectral purity at the output. The quadrature components of the local oscillator available at the quadrature hybrid

outputs are then fed through matched cables to dual phase/amplitude trim networks. It was found that the quadrature hybrid outputs did not have perfectly matched characteristics, therefore to increase the spectral purity of the output when operating at 41.451 MHz (the deuterium observe frequency) it was necessary to add this network. The trim network outputs are fed via another pair of matched cables to the R ports of matched Hewlett-Packard 10514A double balanced mixers.

The input 30 MHz IF is fed to a Merrimac QH-2-30 30 MHz quadrature hybrid. The outputs of this unit were found to have excellent phase and amplitude balance, therefore no trim network was necessary. The outputs are fed via matched cables to the L ports of the two mixers.

The X port outputs were routed via another pair of matched cables to an Anzac TU-50 in-phase power combiner. The power combiner output was fed to a Wavetek 5080 (TNC) 0 - 79 dB step attenuator and another Watkins-Johnson A7 amplifier. The amplifier output was gated by a Watkins-Johnson/Relcom S1 Rf switch, which was controlled from the rear panel switch drive input by a driver circuit (variation of 16X980). The S1 switch was placed last in the signal chain because it has an excellent on/off ratio (approaching 90 dB at low frequencies). Thus when the unit is gated off, there is minimal transmitted noise at the rear panel output. This was significant because the unit generated sufficient output power for internal deuterium locking on the 8 - 270 MHz spectrometer system, and the absence of a higher power transmitter followed by additional non-linear switching devices (*e.g.* series crossed diodes) was found to minimize phase shifts and DC level changes as the lock power level was varied.

Construction Details


Internal construction details of the 16X919 Lock Channel SSB Generator are similar to a number of other devices on the system. Figures 2.47 - 2.49 illustrate the construction methods used.


Figure 2.46 Lock Channel Single Sideband Generator [16X9193]

Figure 2.47 16X919 Lock Channel SSB Generator, Internal Construction
Detail (1)

Figure 2.48 16X919 Lock Channel SSB Generator, Internal Construction
Detail (2)


Figure 2.49 16X919 Lock Channel SSB Generator, Internal Construction
Detail (3)






XBB 797-9252

XBB 797-9249





XBB 797-9251

2.5.7 16X935 63.42 KGauss Magnet Room Temperature Shim Control


The 16X935 63.42 KGauss Magnet Room Temperature Shim Control replaces a unit originally supplied with the magnet system. Principal changes are in the four Z shim controls; the original +/- range switch plus potentiometer arrangement was changed to a continuous (500 center) 10 - turn potentiometer adjustment. In the original configuration, it was not possible to reach zero shim current.

Fine controls are provided for all four Z shims, and coarse and fine H_0 field offset controls are brought to the front panel. The Z/ H_0 summing amplifier board has incorporated minor changes. A high stability LM108A operational amplifier is used for the field sweep summing amplifier, and an Analog Devices AD741 has been substituted for the original LM741 in the Z summing amplifier.

The unit also incorporates switching and display metering of the lock channel signal amplitude, for convenience in shimming adjustments. An X-Y oscilloscope display, when connected to the 'scope vertical and horizontal outputs, may be switched to display a dispersion, absorption, or undetected 5 KHz signal (along the ordinate) coupled with either a sawtooth or the field correction signal amplitude (along the abscissa). Circuit details are summarized in Figures 2.50 and 2.51.


Figure 2.50 63 KG Magnet System -- Room Temperature Shim Control,
Front Panel Wiring [16X9353-S1A]

Figure 2.51 63 KG Magnet System -- Room Temperature Shim Control,
Z/H₀ summing Amplifier Board [16X9353-S2]



ALL POTENTIOMETERS : 10-TURN LINEAR, 5%
2 WATT

A SHIN		7/1/67 REDRAWN W/ ADDITIONS		63 KG MAGNET SYSTEM -- ROOM	
				TEMPERATURE SHUNT CONTROL	
				FRONT PANEL WIRING	
				SHIN	
				DATE 7/1/67	
				LAWRENCE BERKELEY LABORATORY	
				UNIVERSITY OF CALIFORNIA	
				16X9353-SIA	
				SHIN	
				DATE 7/1/67	
				LAWRENCE BERKELEY LABORATORY	
				UNIVERSITY OF CALIFORNIA	
				16X9353-SIA	
				SHIN	
				DATE 7/1/67	
				LAWRENCE BERKELEY LABORATORY	
				UNIVERSITY OF CALIFORNIA	
				16X9353-SIA	



2.5.8 16X950 Nicolet 1180 Driven Probe Temperature Control

The 16X950 Probe Temperature Control unit is part of a direct digital control system that is designed to maintain software selected temperatures in the probes used in the 63 KG Bruker Instruments solenoid that is the heart of the LCB 8 - 270 MHz Spectrometer System. The complete control system is outlined in Figure 2.52 [16X9503-B1]; it consists of the 16X950 unit attached to the Nicolet Instruments 1180 Data System, a Nicolet 293A Programmable Pulser, a programmable power supply (16X949 or equivalent), and the required heater and gas transfer line assemblies.

Actual probe temperature is measured by a copper/constantan thermocouple that is an integral part of each probe. The temperature control unit generates periodic interrupts on 1180 Interrupt Level 6, and an interrupt service routine that is written as an overlay to the Nicolet Technology Corp. NTCFT-1180 software package reads this temperature, compares it to the desired setpoint, and adjusts the output of the programmable power supply/heater accordingly. Because the thermocouple is physically located in the gas stream just ahead of the sample tube, temperature measurements actually only reflect the probe temperature, which may not be the same as the actual sample temperature. In instances when the decoupling is not in use and the sample has been equilibrated for 15 minutes or more, these temperatures are generally quite close. When proton decoupling is in use with aqueous samples, dielectric heating will increase the sample temperature above the indicated reading by a significant amount, depending on decoupler power level. If the actual sample temperature must be known quite accurately,

an external calibration should be run. The measurement error may be minimized to some extent by using a high variable temperature gas flow rate. When proton decoupling is used with non-aqueous samples, the dielectric heating problem is less serious. It is not expected that homonuclear decoupling should cause any measurement errors, because of the much lower Rf power levels involved.

Design Summary

A. Interface Hardware

The 16X950 unit utilizes a Fluke 2170A Digital Thermometer equipped with the Option -02 Digital Output Unit. The thermocouple connections on the 2170A are connected to a rear panel thermocouple plug with a copper/constantan pair that is wound through several ferrite beads. This aids in the suppression of decoupling frequencies that are picked up in the thermocouple. Rectified common-mode input can otherwise lead to significant measurement errors.

The 2170A makes approximately 2.5 readings per second. At the end of each measurement period, the output register in the digital output unit is updated. While updating is in progress, the BUSY line is held high to indicate invalid data.

Interface circuitry is contained on two Control Logic Corp. (CLC) breadboards -- the Device Select and Interrupt Generator Board (Figure 2.53) and the 20-Bit D-Register Board (Figure 2.54). The 2170A is connected via a 34 conductor ribbon cable to the D-register board, and both boards plug into a CLC cardcage that is connected to rear panel

44 conductor edge connectors that are compatible with the 16X904 extended 1180 I/O Bus (refer to 16X904 description, Section 2.5.6).

The falling edge of the BUSY pulse is used to clock a 74LS196 divide-by-ten counter that drives a 74LS221. The 74LS221 monostable issues a 150 nsec pulse approximately every four seconds. This pulse is used to set the interrupt flip-flop (74LS74) using the preset line. When the Q output is high, the device is requesting a Level 6 interrupt. The 74LS38 NAND on the preset line checks to insure that an interrupt service is not currently in progress. The 74LS38 NAND on the Q output has an open collector output, which is compatible with the wired-or nature of the INT6- line of the 1180.

The temperature controller is assigned device address 311. When the INT6- line is pulled down, the processor jumps to the interrupt 6 service routine. Since all slow I/O devices are also on Level 6, the processor must first determine which device caused the interrupt. This is done by polling all devices on Level 6 in a skip-chain. Each possible device code (among those known to be attached to the system as set in the configuration number) is presented on device address lines DV3 - DV10 along with an ACK pulse. If device code 311 is presented, it will be decoded by the 8836s and the 74LS30. The NAND of this signal with ACK (74LS132) are compared to the output state of the interrupt flip-flop with the 74LS38 NAND. If the interrupt flip-flop was set, the SKIP- (wired-or) line is pulled down, signalling the processor that it was the temperature controller that caused the interrupt.

The processor reads out the temperature on to the World lines

W0 - W19 by issuing a device address 311 with a TALK and an IOP0 pulse. This causes the generation of a LOAD BUS signal at the output of the DS8836. This is inverted and routed to the OUTPUT CONTROL inputs of the five 74LS173s that make up the 20 bit D-register. Upon enabling these three-state outputs, the current temperature (in BCD format) and sign and check bits are placed on the I/O bus and loaded into the 1180 accumulator (read temperature control instruction = RTCIN = 603111).

The 20 bit D-register is loaded from the 2170A-02 by the rising edge of the BUSY signal. Thus at the end of each updating cycle in the 2170A, the outputs are dumped into the 20 bit D-register.

The temperature control unit may be switched between ACTIVE and STANDBY states by the processor. A device select 311 coupled with an IOP1 pulse (1180 instruction TCON = 643112) turns on the ACTIVE/STANDBY flip-flop, and a device select 311 with an IOP2 (1180 instruction TCOFF = 643114) turns the unit off. This control line is also brought to the rear panel so that the programmable power supply may also be controlled (using a solid state relay). The ACTIVE/STANDBY flip-flop must indicate that the controller is active for the 74LS38 on the INT6-line to be capable of generating interrupts.

The processor INIT- line is used to clear all registers on power up or in the event of a front panel (1180) clear command. The INIT-line is used in the controller to generate a pulse that clears both the ACTIVE/STANDBY flip-flop and the interrupt flip-flop. This insures that abnormal conditions will force the controller into a STANDBY state.

The processor calculates the desired setpoint for the programmable power supply, and outputs a voltage on DAC2 of the 293A. This signal

is available on the DX rear panel connector. The DX lines are routed to a daisy-chain pair of Canon DP-25 rear panel connectors on the 16X950 chassis, and the DAC2 signal is routed to the Analog Output Buffer Board (Figure 2.55, 16X9503-S3). A differential amplifier picks up this signal, and after passing through a front panel 10-turn pot (LOOP GAIN control) and an LM310 voltage follower, it is routed to the programmable power supply via a rear panel BNC connector.

B. Interrupt Service Software

Interrupt service software modifications to NTCFT-1180 Version 90622 are contained in the system ASCII file VT.ASC on the 1180 system. This file is assembled using either FASM, the Fast 1180 Assembler, or DASM, the old 1180 Disk Assembler.


The overlay compensates for the incompatibilities between the 16X950 unit and hardware supplied by Nicolet Technology Corp. Differences are summarized in the internal program documentation. The assembled listing is shown in Figure 2.56, appended to the end of this section.


The direct digital control algorithm used is the velocity form of the proportional-integral-differential (PID) 3-mode control algorithm. Term weightings are determined empirically.


The interrupt service routine also outputs the measured temperature on DAC3 of the 293A. It was found to be useful to monitor the temperature with an analog Y-T recorder when setting control loop gains.

Provisions have been made for adding an alarm board to sense the absence of gas flow or other malfunctions. Error conditions are transmitted to the 1180 by holding bit 19 (W19) high.

- Figure 2.52 8 - 270 MHz Spectrometer System Probe Temperature
Control System Interconnect [16X9503-B1]
- Figure 2.53 Temperature Control Device Select and Interrupt
Generator [16X9503-S1]
- Figure 2.54 Temperature Control 20 Bit D-Register Board
[16X9503-S2]
- Figure 2.55 Temperature Control Analog Buffer and Status Indicator
Board [16X9503-S3]










Figure 2.56 Assembled Listing of Temperature Control Interrupt
Service Overlay to NTCFT-1180

RUN FASM
FASM, E58-90215
@VT.ASC/VT.BIN,-TT:F

FAST 1100 ASSEMBLER, E58-90215
08/24/79 16:25 @VT.ASC/VT.BIN,-TT:F
110 SOURCE STATEMENTS
27 USER-DEFINED SYMBOLS
MEMORY NEEDED: 4571 TO 4661

/OVERLAY TO NTCFT V#90603 FOR 16X950 TE

```

/OVERLAY TO NTCFT V#90603 FOR 16X950 TEMPERATURE CONTROL
.UNDENT
/
/
/
/      WILLY C. SHIH
/E90704 *****      4 JULY 1979      *****
/
/      This version uses the velocity form of the FID control
/      algorithm, and reads temperatures to 0.2 deg C.
/
/
/MISCELLANEOUS COMMAND DEFINITIONS
620621      LDAC2=620621      /LOAD DAC2
620631      LDAC3=620631      /LOAD DAC3
603111      RTCIN=603111      /READ TEMPERATURE CONTROL REGISTER ONTO BUS
643112      TCON=643112      /TURN ON TEMPERATURE CONTROL
643114      TCOFF=643114      /TURN OFF TEMPERATURE CONTROL
/
/LOCATION EQUIVALENCES
4047      TSETI=4047
4050      DTIM1=4050
4051      EQUILT=4051
4052      TVOLTS=4052
4053      TPASS=4053
4054      TSTORE=4054
4055      TM0=4055
4056      TM1=4056
4057      TM2=4057
4060      TPT=4060
4061      TCT=4061
4060      DTI=TPT
2666      SLOIOX=2666
564      GOFLAG=564
605      QCOUNT=605
206      TTYOUT=206
/
4571      *4571
4571 603111 TCINT6,RTCIN      /READ TC REGISTER ONTO BUS
4572 735322      TRLSH 22      /ROTATE ALL BUT SIGN AND CHECK TO MQ
4573 670400      SKIP AC0      /SIGN BIT=1 => NEG TEMP
4574 1503060     ZERMZ DTI      /TEMP>0 => STORE 0 @ DTI AND SKIP
4575 1473060     MONM DTI      /TEMP<0 => STORE -1 @ DTI
4576 262002     ANDA 12      /CHECK BIT 19 FOR ERROR CONDITION
4577 44607      JPZ TC1      /JUMP IF NO ERRORS
4600 1116564     MEMZ GOFLAG    /ACQUISITION IN PROGRESS?
4601 1436605     MPOM QCOUNT   /YES, STOP IT
4602 216207      MEMA 1207      /RING BELL
4603 2000206     JMS TTYOUT
4604 4603      JMP #-1
4605 722622      TLLSH 22      /RESTORE MQ TO AC VIA ROTATE
4606 2666      JMP SLOIOX      /IGNORE READING AND QUIT
4607 216003 TC1, MEMA 13      /LOAD DIGIT COUNT FOR TEMPERATURE

```

/OVERLAY TO NTCFT V#90603 FOR 16X950 TE

```

/THIS VERSION TRUNCATES THE TENTHS OF DEGREES
/ FOR THE SAKE OF SIMPLICITY IN GETTING
/ OVERLAY WORKING AND GAIN ADJUSTED
4610 1413061 ACCM TCT /STORE IN TCT
4611 1603050 ZERMA DTIM1 /CLEAR TEMP WORD
4612 722602 TLLSH 2 /GET FIRST 1/2 DIGIT
4613 701003 TC2, LASH 3 /BCD TO BINARY CONVERSION OCCURS
4614 1221050 AFMA DTIM1 / BY SHIFTING LEFT 3 (*8) AND
4615 1421050 AFMM DTIM1 / ADDING DIGIT TO ITSELF TWICE
4616 722604 TLLSH 4 /GET NEXT DIGIT
4617 262017 ANDA (17 /MASK FOUR LSB's
4620 1621050 AFMA DTIM1 /ADD TO RUNNING SUM
4621 1535061 MMOMZ TCT /DECREMENT DIGIT COUNT,SKIP IF DONE
4622 4613 JMP TC2 /CONTINUE IF NOT DONE
4623 1117060 MEMZ DTI /TEMPERATURE < 0?
4624 1433050 MNGM DTIM1 /YES, COMPLEMENT DTIM1
4625 1217050 MEMA DTIM1 /OUTPUT TEMPERATURE ON DAC3 +
4626 220030 AFMA (30 / 30 OFFSET
4627 620631 LDAC3 / FOR MONITORING OR LOGGING
4630 1217054 TC2A, MEMA TSTORE /PUSH CURRENT READING ONTO STACK
4631 1413060 ACCM TPT /TSTORE POINTS TO TOP OF STACK
4632 216003 MEMA (3 /STACK IS ONLY LAST 3 READINGS
4633 1413061 ACCM TCT / JUST MOVE 3 LOCATIONS
4634 1217050 MEMA DTIM1 /GET TEMPERATURE READING
4635 3657060 SWP ~ TPT / AND PUSH ONTO STACK
4636 1535061 MMOMZ TCT /DECREMENT DIGIT COUNT,SKIP IF DONE
4637 4635 JMP #-2
4640 1735053 TC2B, MMOMAZ TPASS /TIME TO EVALUATE?
4641 62666 JPF SLOIOX /NO,EXIT INTERRUPT SERVICE

```

/DIRECT DIGITAL CONTROL ALGORITHM:

```

/
/ DV = 4*(E + 10*E' - E'')
/
/ WHERE: DV = CHANGE IN DAC OUTPUT
/ E = CURRENT TEMPERATURE ERROR
/ E' = 1ST DERIVATIVE OF ERROR
/ E'' = 2ND DERIVATIVE OF ERROR
/
/ E(0) = T(5) - T(0)
/ E'(0) = E(0) - E(-1)
/ = T(-1) - T(0)
/ E''(0) = E'(0) - E'(-1)
/ = 2*T(-1) - T(0) - T(-2)
/

```

```

4642 1217055 MEMA TMO /T0
4643 1223056 AMMA TM1 /T0 - T1
4644 701002 LASH 2 /4*T0 - 4*T1
4645 1221055 AFMA TMO /5*T0 - 4*T1
4646 701001 LASH 1 /10*T0 - 8*T1
4647 1223057 AMMA TM2 /10*T0 - 8*T1 - T2
4650 1225047 MMAA TSETI /TS - 10*T0 + 8*T1 + T2
4651 4661 JMP TC3
4661 *4661

```

/OVERLAY TO NTCFT V#90603 FOR 16X950 TE

4661 701002 TC3, LASH 2 /*4 (GAIN CONTROL)

-3-

/OVERLAY TO NTCFT V#90603 FOR 16X950 TE

DTI	4060	DTIM1	4050	EQUILT	4051	GOFLAG	564
LDAC2	620621	LDAC3	620631	QCOUNT	605	RTCIN	603111
SLOIOX	2666	TCINT6	4571	TCOFF	643114	TCON	643112
TCT	4061	TC1	4607	TC2	4613	TC2A	4630
TC2B	4640	TC3	4661	TMO	4055	TM1	4056
TM2	4057	TPASS	4053	TPT	4060	TSETI	4047
TSTORE	4054	TTYOUT	206	TVOLTS	4052		

* 2

2.5.9 16X955 270 MHz Preamplifier

The 16X955 270 MHz Preamplifier is designed for a 270 MHz center frequency with a 1 dB bandwidth of approximately 50 MHz. The unit employs four grounded-gate FET stages, capacitively coupled, yielding about 35 dB of Rf gain. Measured noise figure (cascaded with the 10X210 8-270 MHz to 30 MHz Linear Converter) is 2.2 dB.

Circuit details are summarized in Figure 2.57. The unit employs Siliconix U310 with a common-gate forward transconductance of 20,000 μmho , with a figure of merit (g_m/C) of greater than 2.35×10^9 .

Figure 2.57 270 MHz Preamplifier [16X9553]

2.5.10 16X956 Nicolet 1180/CalComp 565 Incremental Plotter Interface

The CalComp 565 incremental plotter (California Computer Products) is an industry standard drum-type digital plotter. Control inputs are required at a rear panel Cannon SK-19 connector. These inputs control carriage step left and right, drum step up and down, and pen lift and lower. Inputs are ten volt pulses with a minimum duration of four microseconds.

The Nicolet Instrument Corp. 1180 Data System is normally intended to operate with a Nicolet Zeta Research Model 160 Incremental Plotter system. Connection is made via the supplementary "I/O Tee". This I/O tee is attached to the 1180 I/O bus, which is a three-state bus providing 20 data lines, 8 device address lines, 3 IOP lines, and assorted other control signals. The 16X956 Nicolet 1180/CalComp 565 Interface is connected to the I/O bus via the 16X904 I/O Bus Intercept (refer to Section 2.5.6 16X904 Nicolet 1180 I/O Bus Intercept).

A. Design Summary

The interface unit is connected to the 16X904 I/O Bus Intercept via two 40 conductor ribbon cables (3M Scotchflex 3365). Two rear panel I/O connectors are provided so that the unit may be conveniently daisy-chained. If the unit is the last connection on the bus, a bus terminator board must be inserted in the unused connector.

The rear panel connectors are connected to board 1. The plotter is assigned device address 374. When a plot command is issued, the lowest six bits of the 1180 accumulator contain the desired plot action (Table 2.8).

Table 2.8 16X956 Nicolet 1180/CalComp 565 Plot Command Summary

<u>Accumulator Bit</u>	<u>Desired Action</u>
AC0 = 1	Step carriage right
AC1 = 1	Step carriage left
AC2 = 1	Step drum up
AC3 = 1	Step drum down
AC4 = 1	Lift pen
AC5 = 1	Lower pen


Interface circuit details are summarized in Figures 2.58 and 2.59 [16X9563-S1 and 16X9563-S2]. Upon issuing the plot I/O command (1180 instruction DZETA = 623741), the accumulator is output to the 20 world lines (W0 - W19). The 74LS14 and 74LS30 decode the 374 device select, and with the issuance of the IOP0 pulse, W0 - W5 are clocked on the leading edge into the 74LS174 6 bit D-register. The outputs of this register then reflect which of the 75451 peripheral drivers should issue a pulse.


The device select/IOP0 pulse also triggers two timers (NE556). One timer generates a ten microsecond pulse for the peripheral drivers (in general, only one of the six will have been selected). The other timer generates a delay to allow sufficient time for the plotter to respond. Before issuing a plot command, programs either calculate a software delay or issue a "skip if ready" command. The "skip if ready" command (1180 instruction 647740) generates a device select 374 and an ACK pulse. If the NE556 timer indicates that sufficient time has elapsed since the last plot command, the ANDed device select and ACK will cause the 75452 peripheral driver to pull down the -SKIP line, which is wired-or. When the processor senses that the -SKIP line is pulled down, it assumes that the plotter is ready to receive the next step command.

The board 1 plotter outputs are routed to complementary emitter-followers on board 2. This insures that the plotter step pulses will be driven from a low impedance source with enough current capability for the CalComp. The board outputs are routed to a front panel MS connector, where a shielded cable makes connections to the plotter.

Figure 2.58 Nicolet 1180/CalComp 565 Interface, Plot Control Logic [16X9563-S1]

Figure 2.59 Nicolet 1180/CalComp 565 Interface, Line Driver Detail
[16X9563-S2]





B. Requisite Software Changes

The NTCFT-1180 software package supplied by Nicolet Technology Corp. is designed for use with a Nicolet Zeta Research Model 160 plotter with a step size of 0.005 inch. The CalComp 565 at LCB may use a step size of 0.005 or 0.01 inch, depending on the drive gear selection. Since the LCB CalComp is also used with the general laboratory data acquisition system with a step size of 0.01 inch, a large body of software would have to be modified if the step size were to be changed. It was also deemed to be impractical to change the gears whenever the plotter was moved to a different machine, therefore the only alternative was to modify the NTCFT-1180 software. An assembled listing of these changes is shown in Figure 2.60.

Figure 2.60 Assembled Listing of CalComp 565 Modifications to
NTCFT-1180

RUN FASM
FASM, E58-90215
@CALCMP.ASC/CALCMP.BIN,-TT:F

FAST 1180 ASSEMBLER, E58-90215
08/23/79 10:01 @CALCMP.ASC/CALCMP.BIN,-TT:F
70 SOURCE STATEMENTS
3 USER-DEFINED SYMBOLS
MEMORY NEEDED: 12640 TO 16107

CALCOMP PLOTTER PATCHES - 90615

.TITLE "CALCOMP PLOTTER PATCHES - 90615"
 .UNDENT

/PATCHES TO NTCFT-1180 - 90622 FOR DRIVING 100 STEP/INCH
 /CALCOMP PLOTTER RATHER THAN ZETA 160 (200 STEPS/INCH)

/PATCHES TO DIR MODULE

14227		*14227		
14227	216004		MEMA (4	/TO SLOW DOWN PLOTTER
14310		*14310		
14310	217730		MEMA (1730	/STEPS/25 CM
14631		*14631		
14631	114713		114713	/STEPS/1000 CM

/PATCHES TO OUT MODULE

12640		*12640		
12640	1522	PAGSIZ, 1522		/STEPS/8.5 IN (=PAGE)
15507		*15507		
15507	241157		DBMA (1157	/STEPS/PAGE (VERTICAL
15515		*15515		
15515	216100		MEMA (100	/-> BOTTOM OF PLOT

/PATCHES TO AXL MODULE

13010		*13010		
13010		1 CHARSZ,1		/DEFAULT CHARACTER SIZE
13063		*13063		
13063	216001		MEMA (1	/CHAR SIZE FOR AXIS
13073		*13073		
13073	220014		APMA (14	/AXIS Y-OFFSET
13303		*13303		
13303	3777770		-10	/SMALL TIC SIZE
13307		*13307		
13307	10		10	/ " " "
13324		*13324		
13324	3777757		-21	/LARGE TIC SIZE
13330		*13330		
13330	21		21	/ " " "
13347		*13347		
13347	232050		MNGA (50	/AXIS NUMBER Y-OFFSET
13402		*13402		
13402	222060		AMMA (60	/= 3 CHAR WIDTHS
13404		*13404		
13404	222007		AMMA (7	/HALF SPACE
13432		*13432		
13432	220064		APMA (64	/CHAR Y-OFFSET
13435		*13435		
13435	220040		APMA (40	/TEST FOR END OF AXIS
13704		*13704		
13704	7		7	/HALF CHAR SPACE
15047		*15047		

CALCOMP PLOTTER PATCHES - 90615

15047	217234		MEMA (1234	/=17 CM
15051		*15051		
15051	220450		APMA (450	/-> 24.5 CM
15250		*15250		
15250	220043		APMA (43	/FOR LEFT EDGE OF PLOT
15260		*15260		
15260	3777634		-144	/PARAMETER POSITIONING
15261	3777714		-64	
16004		*16004		
16004	1522	PGSIZ,	1522	/=8.5 IN
16046		*16046		
16046	222006		AMMA (6	/= HALF CHAR (X)
16057		*16057		
16057	222007		AMMA (7	/= HALF CHAR (Y)
16107		*16107		
16107	3777764		-14	/= 5 MM

-2-

CALCOMP PLOTTER PATCHES - 90615

CHARSZ	13010	PAGSIZ	12640	PGSIZ	16004
--------	-------	--------	-------	-------	-------

*

2.5.11 16X964 41.45 MHz Receiver - T/R Switch

The 16X964 41.45 MHz Receiver T/R Switch operates at a center frequency of 41.451 MHz for ^2H observation. The unit serves two principal functions:

- Front-end Rf amplification. The preamplifier in this unit supplies approximately 18 dB of gain. It employs a low noise junction FET amplifier which sets the noise performance of the deuterium receiver channel.
- Transmit/receive switching. When the transmitter produces a pulse, the unit protects the preamplifier from overloading. After the pulse, the unit switches to the receive mode.


Design Summary (Figure 2.61)

The unit employs a 16X903 dual J-FET cascode amplifier (refer to Section 2.5.4). Protective crossed diodes are incorporated in the input, and a series LC shunt tuned to 270 MHz is employed to suppress decoupling frequencies. A semi-rigid UT-141 0.141 quarter-wave line is used as part of the switching arrangement. A Circutel FLT/20-140-13/150-2A/2A 140 MHz low-pass filter is used for additional suppression of decoupling frequencies.

Instead of the conventional diode switch, only a 3 dB pad is used to isolate the lock channel transmitter (16X919) from the receiver. This arrangement was used because the high isolation Rf switch used in the 16X919 final stage has an excellent on/off ratio, and the addition of diode switches only resulted in non-linearities when the Rf power

level of the transmitter was varied. The unit employs internal supply voltage regulation, and is normally powered by the 16X974 pre-amplifier power supply (+24 VDC).

Figure 2.61 41.45 MHz Receiver - T/R Switch [16X9643]



2.5.12 16X965 27.36 MHz / 270 MHz Dual Frequency Receiver - T/R Switch

The 16X965 dual frequency T/R switch assembly operates at 27.36 MHz for ^{15}N and at 270 MHz for ^1H observe. The unit serves three principal functions:

- Front-end Rf amplification -- The two preamplifiers in this unit are the first gain stages encountered by ^{15}N and ^1H signals. They provide approximately 36 - 40 dB of gain, and the outputs are fed to the 10X210 8 - 270 MHz to 30 MHz Linear Converter for conversion to the 30 MHz intermediate frequency. Since they are the first gain stages, noise performance is crucial.
- T/R switching -- When the transmitter produces a high power pulse, the unit protects the sensitive front-end receivers. After the pulse, the unit switches to the receive mode.
- Transmitter noise suppression -- Even when the transmitter is gated off, it has a high noise temperature. The 16X965 unit effectively isolates this noise source from the receivers, leading to minimal degradation in system noise performance.

Design Summary

A. 27.36 MHz Section


The 27.36 MHz section uses two cascaded 16X903 low noise dual J-FET cascode amplifiers (refer to Section 2.5.4). The first amplifier has protective crossed diodes on the input, and it also employs a series LC shunt tuned to 270 MHz to suppress proton decoupling frequen-

cies (see photograph of Figure 2.45). The amplifier cascade, when coupled to the 10X210 Linear Converter, shows a system noise figure of 1.5 dB when measured at the 30 MHz IF output. The T/R switching employs a standard quarter-wave line arrangement. Transmitter isolation was found to be sufficient when four pairs of crossed diodes were placed in series with the transmitter output, and a measured line was used to connect the transmitter.

B. 270 MHz Section

The 270 MHz section employs a 16X955 four stage grounded gate J-FET amplifier to provide 36 - 40 dB of Rf gain (see Section 2.5.9, 16X955 270 MHz Preamplifier). This unit has protective crossed diodes on the input. When coupled to the 10X210 Linear Converter, the system noise figure is approximately 2.2 dB when measured at the 30 MHz IF amplifier output. The T/R switching is done with a standard quarter-wave line arrangement, using UT-141 0.141 semi-rigid coax for the line segment. The junction capacitance of the crossed diodes in the transmitter switch is tuned with a short coil to maximize the impedance at 270 MHz and cancel as much reactive component as possible.

Figure 2.62 27.36 MHz / 270 MHz Dual Frequency Receiver - T/R Switch
Assembly [16X9653]



2.5.13 16X966 67.89 MHz / 109.30 MHz Dual Frequency Receiver T/R Switch

The 16X966 dual frequency T/R switch assembly operates at 67.89 MHz for ^{13}C and at 109.30 MHz for ^{31}P observe. The unit serves three principal functions:

- Front-end Rf amplification -- The two preamplifiers in this unit are the first gain stages encountered by ^{13}C and ^{31}P signals. They provide approximately 30 - 36 dB of gain, and the outputs are fed to the 10X210 8 - 270 MHz to 30 MHz Linear Converter for conversion to the 30 MHz intermediate frequency. Since they are the first gain stages, noise performance is crucial.
- T/R switching -- When the transmitter produces a high power pulse, the unit protects the sensitive front-end receivers. After the pulse, the unit switches to the receive mode.
- Transmitter noise suppression -- Even when the transmitter is gated off, it has a high noise temperature. The 16X966 unit effectively isolates this noise source from the receivers, leading to minimal degradation in system noise performance.

Design Summary

A. 67.89 MHz Section


The 67.89 MHz section uses an RHG ICFT7030 thin-film hybrid preamplifier, protected by passive quarter-wave line switches. For this operating frequency, it was found that additional quarter-wave line sections were necessary for good noise performance (see Figure 2.63).

In the transmit mode, the two quarter wave sections that are terminated with diode pairs will have the diodes turned on, therefore they will look like open circuits $\lambda/4$ away from the diodes. In the receive mode, the diodes will be off; therefore the quarter-wave section closest to the transmitter will look like a short. This sets a boundary condition and forces the transmitter leg of the lower tee to look like an open circuit.

B. 109.30 MHz Section

The 109.30 MHz section uses an RHG ICFT12040 thin-film hybrid preamplifier, with a conventional quarter-wave/diode switching arrangement. It was found that a large number of isolated diode pairs was necessary for the transmitter switch.

Figure 2.63 67.89 MHz / 109.30 MHz Receiver - T/R Switch [16X9663]



ALL CABLE PULSES KINC K1A 59-57 OR K1A 59-86
UNLESS SPECIFIED CABLE - RG 55 B/U

[illegible]

2.5.14 Diablo 31 Disk Memory Subsystem -- Power Supply Interconnect

Three Diablo 31 series 2315 type cartridge disk memory subsystems are presently installed on the Nicolet 1180 Data System. Drive unit 0 is powered by a Diablo Model 029 8 amp \pm 15 VDC supply; drives 1 and 2 are powered by the supply arrangement shown in Figure 2.64. Because the Dynage 700 KHD series supplies are not as high in current capability as the Diablo 029, capacitor bank buffers are provided. These are attached by very short cables to the drive rear panels.

2.5.15 Rf Power Amplifiers -- Power Supply Interconnect

The 10X224 and 10X229 high power pulse amplifiers are powered by modified Universal Voltronics BRE-3-1000-LRL-1 3 KV 1 A high voltage power supplies. These supplies are buffered by 10X227 Capacitor Banks to minimize pulse droop during long pulses. Upon turn on, the current surge required to charge the capacitor banks would normally trip the current overload shutdown sensor, therefore the supplies were modified to include a time-delayed turn on shunt. Interconnections are shown in Figure 2.65

2.5.16 16X973 Rf Switch

The 16X973 Rf Switch is a relay actuated DPDT Rf switch employing Teledyne 712TN5 miniature Rf relays. The circuit layout is shown in Figure 2.66. The relay is mounted on 0.125 G-10 fiberglass-epoxy circuit board (double-sided), and connections to front panel switch inputs are via 50 ohm microstrip. The double-sided board layout is

shown in Figures 2.67 and 2.68

2.5.17 16X974 Preamplifier Power Supply

The 16X974 Preamplifier Power Supply has four switched outputs that may be used to power the 16X964 - 16X965 - 16X966 receiver - T/R switch assemblies. The unit contains two Powermate UNI-88 DC supplies set to provide + 24 VDC and + 15 VDC. Two 15,000 μ fd capacitors are used to provide high quality output filtering. Circuit details are shown in Figure 2.69.

2.5.18 16X978 41.45 MHz Amplifier

The 16X978 41.45 MHz Rf Amplifier is used to provide 0 - 80 dB adjustable Rf gain over a 5 MHz bandwidth centered at 41 MHz. The unit is used as part of the ²H internal locking system, providing adjustable Rf gain prior to frequency conversion to the 30 MHz IF.

Circuit details are shown in Figure 2.70. The unit incorporates and RHG EVT40AJ66 amplifier. Gain is regulated by controlling voltage supplied to the AGC input. A trailing Wavetek 5070 0 - 70 dB step attenuator is used to adjust the output level.

2.5.19 16X980 Rf Switch

The 16X980 Rf Switch is a general purpose, fast, high isolation Rf switch unit. The switching element is a Watkins-Johnson/Relcom S1 switch; a drive circuit is provided that accepts TTL level control signals. Circuit details are shown in Figure 2.71.

Figure 2.64 Nicolet 1180 Data System -- Disk Memory Subsystems
 DC Supplies and Signal Interconnections
 [16X9673-B2]

Figure 2.65 8 - 270 MHz Spectrometer System -- Rf Power Amplifier
 Connections [16X9673-B5]

Figure 2.66 DPDT Rf Switch [16X9733-S1]


Figure 2.67 16X973 Rf Switch PC Layout (Signal Side)


Figure 2.68 16X973 Rf Switch PC Layout (Ground Plane Side)


Figure 2.69 Preamplifier Power Supply [16X9743]

Figure 2.70 41.45 MHz Amplifier [16X9783]

Figure 2.71 Rf Switch [16X9803]










UNLESS SPECIFIED:
ALL RESISTORS - 1/4 5% CARBON,
1% RECURRENT - KINGS KA-79-28


[illegible]








INDICATOR AVAILABLE: DE N/A:

[illegible]



2.5.20 Data Processor Power Fail Interlock

The 16X981 Data Processor Power Fail Interlock employs a General Electric CR205 high power magnetic contactor to switch system AC power to the 10X225 270 MHz Rf Amplifier, the ENI 420L Rf Power Amplifier, and the programmable power supply used to heat the variable temperature gas. The unit(s) are controlled by attachment to the switched AC line of the 1180 Data System. The units protect the system probes from burn out resulting from sudden CW operation of the transmitters, or from full power heating by the programmable power supply. Such conditions may arise if the 1180 Data System is shut down, leaving system control logic in indefinite states.

Figure 2.72 Data Processor Power Fail Interlock
 [16X9813]

2.6 References

1. J.D. Ellet, Jr., M.G. Gibby, U. Haeberlen, L.M. Huber, M. Mehring, A. Pines, and J.S. Waugh, *Adv. Mag. Res.*, 5, 117 (1971)
2. D.J. Sakrison, "Communication Theory: Transmission of Waveforms and Digital Information," John Wiley & Sons, Inc., New York, NY, 1968
3. R.E. Ziemer and W.H. Tranter, "Principles of Communications," Houghton Mifflin Co., Boston, MA, 1976
4. T.C. Farrar and E.D. Becker, "Pulse and Fourier Transform NMR," Academic Press, New York, NY, 1971
5. M.E. Stoll, A.J. Vega, and R.W. Vaughan, *Rev. Sci. Instrum.*, 48, 800 (1977)
6. I.J. Lowe and C.E. Tarr, *J. Phys. E, Ser. 2*, 1, 320 (1968)
7. I.J. Lowe and M. Eisenberg, *Rev. Sci. Instrum.*, 45, 1159 (1974)
8. K.E. Kisman and R.L. Armstrong, *Rev. Sci. Instrum.*, 45, 1159 (1974)
9. P.D. Murphy and B.C. Gerstein, "Analysis and Computerized Design of NMR Probe Circuits," U.S.D.O.E. Report IS-4436, Ames, Iowa, 1978
10. E. Schwartz, *IEEE Trans. Microwave Theory and Techniques*, MTT-16, 158 (1968)
11. R.M. Fano, *J. Franklin Inst.*, 249, 57 (1950)
12. R.M. Fano, *J. Franklin Inst.*, 249, 139 (1950)
13. D.C. Youla, *IEEE Trans. Circuit Theory*, CT-11, 30 (1964)
14. D.I. Hoult and R.E. Richards, *J. Mag. Resonance*, 24, 71 (1976)
15. R.H. Lyddane and A.E. Ruark, *Rev. Sci. Instrum.*, 10, 253 (1939)
16. W. Franzen, *Rev. Sci. Instrum.*, 33, 933 (1962)
17. W.J.C. Grant and M.W.P. Strandberg, *Rev. Sci. Instrum.*, 36, 343 (1965)
18. Service Manual for NIC-1180 Data System, Nicolet Instrument Corp., 1976
19. A. Abragam, "The Principles of Nuclear Magnetism," London, Oxford University Press, 1961, p. 83
20. H.D.W. Hill and R.E. Richards, *J. Phys. E, Ser. 2*, 1, 977 (1968)

APPENDIX 1 Nicolet 1180 / VAX 11/780 Data Transmission Protocol

A1.1 Introduction

The 8 - 270 MHz spectrometer system may be connected via the 1180 data system to two other data processors: the Nicolet NIC-80 attached to the LCB 220 MHz ^1H pulsed Fourier Transform system, or to a Digital Equipment Corporation VAX 11/780. Both connections are implemented with standard RS232 serial links, the former at data rates up to 38.4 Kbaud, and the latter at data rates up to 9600 baud. The purposes of the 1180/VAX connection are:

- Archival storage of data on magnetic tape. The 1180 data system uses IBM 2315 type 2200 bpi 3.0 Mbyte disks for mass storage, which are not ideal for archival storage purposes.
- Data files stored on the VAX are readily manipulated using FORTRAN language routines. Of particular value is the use of FORTRAN display processing routines for 2D NMR data.
- Data files may be readily reformatted on the VAX, making them transportable to other systems for use with foreign analysis and display processing routines (such as the Surface Display Library at the Berkeley CDC 6400 facility).

The current data transmission protocol does not conform to the standards of the CAL Data Acquisition Facility (CALDAF) being implemented on the LCB VAX. The particular approach was chosen because of the speed with which it could be implemented.

A1.2 Data Transmission Protocol

.1 1180 Memory Organization and Interrupt Structure

The 1180 system presently includes $24K \times 20$ bits of non-ECC memory. Programs are generally written to occupy locations 00000_8 - 17777_8 , and data is usually acquired into or processed in the balance (20000_8 - 57777_8). Reserved locations are summarized in Table A1.1.

Memory locations below 400_8 are reserved for DEXTER/2 disk monitor functions. The NTCFT-1180 program package supplied by Nicolet Technology Corp. is the normal resident program. This package consists of a main program (NTCFT) occupying locations 400_8 - 12637_8 , an overlay file (NTCFT.OVL) containing 16 blocks of command processing routines that are selected after a command table look-up, and various external routines for spectrum simulation, etc. Overlays are loaded into locations 12640_8 - 16477_8 ; locations 16500_8 - 17777_8 are used as a scratch area.

VAX data transmission routines are written as part of the USR overlay. They were composed as .ASC (ASCII text) files using the DEDIT disk-based text editor, then assembled using the FASM Fast 1180 Assembler. Stored as a .BIN file (USROVL.BIN), the routines are incorporated into NTCFT via the P↑A command (load overlay). The FASM assembled listing is appended to the end of this section.

.2 1180 Interrupt Structure

The 1180 has a seven-level vectored priority interrupt structure. Interrupt lines are wired-or bus lines (refer to print 16X904) that are asserted by grounding at any point along the bus. The interrupt vector

addresses and devices associated with each priority level are summarized in Table A1.2.

When an interrupt occurs, the processor completes the current instruction and then does an indirect subroutine jump through the interrupt vector address corresponding to the priority of the interrupt. The address of a software interrupt service routine must be pointed to by this address; upon performing the JMS, the return address is deposited at the start of the service routine for an indirect jump return upon completion.

Interrupt levels are enabled or disabled using the interrupt level mask register. Devices and levels are enabled only when program dictates require that they be used. When they are no longer needed, they are disabled immediately. The VAX data transmission routines use the Level 5 interrupt for the RS232 B Channel UART (Universal Asynchronous Receiver-Transmitter). Upon entering the routine, Level 5 and the RS232B are enabled in the mask registers. On exiting, a check is done to see if data acquisition is in progress, because the 293A Programmable Pulser is also on Level 5. If data acquisition is in progress, only the RS232B bit is reset in the Slow I/O mask. If data acquisition is not in progress, both the RS232B and the Level 5 bits are reset.

NTCFT-1180 permits simultaneous execution of up to two background processes (plotting and data acquisition) and one foreground process (such as performing a FFT or phasing a spectrum). Data transmission is a foreground process because the serial RS232 format requires partitioning of data words into byte strings; interrupt driven processing would have been very difficult to implement within the

framework of NTCFT-1180.

.3 Hardware Connections, VAX Configuration

Hardware connection of the 1180 to the VAX employs the RS232 standard (2 wire + ground). The 1180 connection is to the RS232 B Channel UART (an IM6402 CPD UART was substituted for the slower COM 2502 that was originally supplied, the former has a much higher baud capability); the VAX connection is to one of the ports on the terminal multiplexer. It does not matter which VAX port is connected, but the port number must be known.

Each terminal multiplexer port on the VAX has its own associated buffer space. Communications lines connected to the ports may operate in several modes:

1. HOSTSYNC -- the VAX synchronizes the flow of input into the terminal multiplexer buffer by issuing CTRL-Q to start the data stream or CTRL-S to stop. Default mode is NOHOSTSYNC.
2. TTSYNC -- Device attached to assigned terminal multiplexer port synchronizes output from VAX by issuing CTRL-Q to start data flow and CTRL-S to stop. Default mode is TTSYNC.

A VAX process may accept data from several sources:

1. Interactively, from the terminal that is controlling the process.
2. From a file stored on the system.

3. From an allocated external device attached to an assigned port.

VAX / 1180 data communications programs were written to accept data in the last mode, that is through an allocated port on the terminal multiplexer. This obviates the need to reconnect terminals and computers whenever such a process is to be run. Prior to running any of these programs, it is necessary to ALLOCATE the port, ASSIGN logical files to the port, and to set the port in HOSTSYNC mode. An example sequence that would configure port TTB6 correctly would be:

ALLOCATE TTB6 FOR001	/ALLOCATE TTB6, ASSIGN LOGICAL
	/ NAME FOR001 TO IT
ASSIGN TTB6: FOR002	/ALSO ASSIGN THE LOGICAL NAME
	/ FOR002 TO IT
SET TERM TTB6/HOSTSYNC	/SET HOSTSYNC CONTROL

.4 Data Format

The VAX interprets the input byte stream in terms of ASCII characters and control signals. Because the bit patterns of the 1180 data are random, it is entirely possible that certain data bytes would be misinterpreted as control characters, thereby causing undesired VAX actions. To avoid this possibility, the 1180 segments each 20 bit word into four 5 bit strings, which are then transmitted in an unpacked format with a leading 001. Thus each byte has the format 001XXXXX, where X represents a data bit. The VAX FORTRAN programs FILECOPY and ARRAY then repack each string of four bytes into a 20 bit word, and they then write these to files as 32 bit integers. The ASCII equivalents of 0010 0000 through 0011 0000 do not include any control characters.

Alternatively, the VAX terminal multiplexer port could be set to PASSALL mode, which would permit the input of byte strings without interpretation. This mode is considered dangerous, however, because its use requires I/O privileges.

.5 Timing

When data transmission is running under HOSTSYNC control, the VAX interrupts the 1180 whenever the terminal multiplexer buffer is full and it is desired to stop transmission. Similarly when the buffer can accept more bytes, the 1180 is interrupted. Data transmission on the 1180 is a foreground process; the flag CTRLCH is inspected prior to any four byte transmission. This flag is set under interrupt control by the VAX control signals. Timing is well illustrated by the flow chart of Figure A1.1.

An 1180 assembly language listing for the data transmission routine, and a VAX FORTRAN IV PLUS listing of the routines FILECOPY and ARRAY are appended to the end of this section.

.6 1180/NIC-80 Data Transmission

The 1180/NIC-80 data link was used to support data storage for the NIC-80, which had no mass storage peripheral. The link was set to operate at 38.4 Kbaud to a NIC-80 modified per drawing 16X943. The 1180 data transmission routine listing is appended to the end of this section.

Table A1.1 1180 Data System -- Resident Monitor Reserved
Memory Locations

<i>OCTAL LOCATION</i>	<i>CONTENTS</i>
0	JMP to 112, Keyboard Monitor Restart
1-7	Levels 0 - 6 Interrupt Vectors
10	Reserved for DBUG11, NICBUG
11	Software copy of interrupt level enable mask
12	Software copy of slow I/O device enable mask
13	Pointer to overlay exit routine
14	MONR subroutine
17	PUSH subroutine
22	POP subroutine
25	OVRLD overlay loader
30	RDISK disk read service routine
32	WDISK write disk subroutine
34-377	Miscellaneous monitor functions
360007	Paper tape binary loader for high speed reader
360010	Paper tape binary loader for low speed reader
360114	Disk monitor cold start from sector 0, disk drive unit 0
360112	Disk monitor backup bootstrap from sector 6250, disk drive unit 0

Table A1.2 Interrupt Vector Addresses and Device Priority Level
Assignment -- Nicolet 1180 Data System

<i>LEVEL</i>	<i>ADDRESS</i>	<i>ASSIGNED DEVICES</i>
0	1	ADC - in 256 word mode
1	2	ADC - end of sweep
2	3	Clock 1 - general purpose
3	4	Disk
4	5	Clock 2 - Display
5	6	RS232 - B Channel UART, 293A Programmable Pulser
6	7	RS232 - A Channel UART, High speed paper tape reader/punch, 16X950 Variable Temperature Controller

Figure A1.1 Assembled Listing of USR Overlay to NTCFT-1180
 for Data Transmission to VAX 11/780

RUN FASM

FASM, E58-90215
@USROVL.ASC/USROVL.BIN,-TT:F

FAST 1180 ASSEMBLER, E58-90215

08/24/79 16:16 @USROVL.ASC/USROVL.BIN,-TT:F

200 SOURCE STATEMENTS
38 USER-DEFINED SYMBOLS
MEMORY NEEDED: 4751 TO 16343

/USER OVERLAY TO NTCFT V#90615

/USER OVERLAY TO NTCFT V#90615

/UNDENT

/A90823 WILLY C. SHIH

/23 AUGUST 1979

/U1 =>DATA TRANSMISSION TO VAX11/780 VIA TERMINAL
/MULTIPLEXER. TERMINAL MULTIPLEXER PORT SHOULD
/BE SET AS FOLLOWS:

/ALLOCATE TTB6 FOR001
/ (ASSIGNS PORT TTB6 TO FORTRAN LOGICAL UNIT 001)
/ALLOCATE TTB6 FOR002
/ (ASSIGNS PORT TTB6 TO FORTRAN LOGICAL UNIT 002)
/SET TERM TTB6/HOSTSYNC
/ (TERMINAL MUX CONTROLS TRANSMISSION)
/SET TERM TTB6/SPEED=(9600,9600)
/ (9600 BAUD TRANSMISSION)

/U1 WILL TRANSMIT THE DISPLAYED AREA, USE THE FS COMMAND
/TO SET TO A PARTICULAR NUMBER OF POINTS. U1
/MAY BE USED IN A LINK FOR SENDING SETS OF FILES

/U2 =>SAME AS U1, EXCEPT WITH ACCOMPANYING PARAMETER TABLE
/ (LOCATIONS 640-1377, 430 OCTAL WORDS)

466	DSFSIZ=466	/SIZE OF DISPLAYED REGION
467	DSFSTA=467	/START OF DISPLAYED REGION
11	ZILMSK=11	/COPY OF INT LEVEL MASK
12	SPDMSK=12	/COPY OF SLOW I/O MASK
10061	GADCHK=10061	/CHECK FOR AIP IN DISPLAY REGION
564	GOFAG=564	/ACQUISITION IN PROGRESS FLAG
2301	ZDISPLA=2301	/DISPLAY MEMORY SUBROUTINE
2533	ZUPKBUF=2533	/MESSAGE OUTPUT ROUTINE
2627	ZWFCHAR=2627	/CHARACTER INPUT
2636	ZTTYBUF=2636	/ZWFCHAR+ZACCBUF
545	DIFLAG=545	
10152	BRETRN=10152	/RETURN TO NTCFT INTS SERVICE

/ADDITION TO INT5- SERVICE ROUTINE

10121	*10121		
10121	44751	JPZ RSBINT	/JMP ON ZERO (i.e. 293A DID NOT /CAUSE INTERRUPT
4751	*4751		
4751	644720	RSBINT, RSINF	/RS232B READY TO READ
4752	4751	JMP #-1	/NO
4753	600721	RSIN	/READ ONE BYTE
4754	126021	SNQ (21	/VAX11 START CHAR = CTRL-Q (021)

/USER OVERLAY TO NTCFT V#90615

```

4755 1403761      ZERM CTRLCH      /DEPOSIT 0 FLAG IF START RECEIVED
4756 126023      SNQ (23          /VAX11 STOP CHAR = CTRL-S (023)
4757 1471761      DNEM CTRLCH     /DEPOSIT 1 FLAG IF STOP RECEIVED
4760 10152        JMP BRETRN      /RETURN
4761      1 CTRLCH, 1
13000      *13000
/Routine FOR UNPACKED DATA TRANSMISSION TO VAX 11 TERMINAL MUX
13000      0 VAXDMP, 0          /BYTE FORMAT = 001XXXXX
13001 2010061      JMS GADCHK      /CHECK IF ACQ IN PROGRESS IN
                        /CURRENT DISPLAY AREA
13002 2002533      JMS ZUPKBUF     /WRITE MESSAGE
13003      13176      MSG1         /MESSAGE POINTER
13004 2002301      JMS ZDISPLA     /DISPLAY MEMORY ONCE
13005 1403172      ZERM CHKSUM     /ZERO CHECKSUM
13006 1216467      MEMA DSPSTA     /GET START OF DISPLAY REGION
13007 1413174      ACCM START      / AND PUT IT IN START
13010 1216466      MEMA DSPSIZ     /GET SIZE OF DISPLAY REGION
13011 1413175      ACCM SIZE       / AND PUT IT IN SIZE
13012 2013114      JMS INTENA     /ENABLE LVL5 INTERRUPTS
13013 2013065      JMS DUMP        /DUMP DISPLAYED REGION(CTRL-E)
13014 1217172      MEMA CHKSUM     /TRANSMIT CHECKSUM
13015 2013127      JMS VAXOUT      /CALL RS232B XMIT WORD ROUTINE
13016 1434545      MMOM DIFLAG     /REENABLE "O,"X
13017 1217164      MEMA MASK3      /DISABLE RS232B IN SLOW I/O
13020 1462012      ANDM SPIMSK     / MASK COPY AND MASK
13021 500761      LDMASK          /LOAD SLOW I/O DEVICE MASK
13022      12      SPDMASK
13023 2117171      MEMZ @ GFLAG     /ACQ IN PROGRESS
13024 1013000      JMP @VAXDMP     /RETURN WITHOUT DISABLING LVL5
                        /SINCE 293A IS ALSO ON LVL5
13025 1217165      MEMA MASK4      /OTHERWISE RESET INTS-
13026 1462011      ANDM ZILMSK     /STORE IN MASK COPY AND
13027 503771      LMASK          /LOAD INT LEVEL MASK
13030      11      ZILMSK
13031 1013000      JMP @VAXDMP     /NORMAL RETURN
/
/Routine FOR UNPACKED DATA TRANSMISSION TO VAX11 TERMINAL MUX
/ WITH ACCOMPANYING PARAMETER TABLE
13032      0 LEADER, 0          /PARAM TABLE @ 640-1377 (430 WORDS)
13033 2010061      JMS GADCHK      /CHECK IF ACQ IN PROGRESS
13034 2002533      JMS ZUPKBUF     /WRITE MESSAGE
13035      13213      MSG2         /MESSAGE POINTER
13036 1403172      ZERM CHKSUM     /ZERO CHECKSUM
13037 216640      MEMA (640        /PARAMETER TABLE STARTS AT 640
13040 1413174      ACCM START      /PUT START ADDRESS IN START
13041 216540      MEMA (540        /PARAM TABLE IS 540 WORDS (OCTAL)
13042 1413175      ACCM SIZE       /PUT IN SIZE
13043 2013114      JMS INTENA     /ENABLE LVL5 INTERRUPTS
13044 2013065      JMS DUMP        /OUTPUT PARAMETER TABLE
13045 2013075      JMS STATUS      /NOW GET SIZE OF DATA REGION
13046 2013065      JMS DUMP        /DUMP BUTTON REGISTER SELECTED AREA
13047 1217172      MEMA CHKSUM     /LOAD CHECKSUM
13050 2013127      JMS VAXOUT      /AND WRITE IT TO VAX ALSO
13051 1434545      MMOM DIFLAG     /REENABLE "O,"X

```


/USER OVERLAY TO NTCFT V#90615

13052	1217144	MEMA MASK3	/DISABLE RS232B IN SLOW I/O MASK
13053	1462012	ANDM SPDMSK	/AND IT WITH MASK COPY
13054	500761	LDMSK	/LOAD SLOW I/O MASK AND COPY
13055	12	SPDMSK	/ MASK COPY POINTER
13056	2117171	MEMZ @ GFLAG	/ACQ IN PROGRESS?
13057	1013032	JMP @ LEADER	/IF YES, RETURN WITHOUT DISABLING INT5-
13060	1217165	MEMA MASK4	/OTHERWISE RESET INT5-
13061	1462011	ANDM ZILMSK	/AND IT WITH MASK COPY
13062	503771	LMASK	/LOAD INT LEVEL MASK
13063	11	ZILMSK	/ MASK COPY
13064	1013032	JMP @ LEADER	/NORMAL RETURN
13065	0 DUMP,	0	/OUTPUT A MEMORY REGION ROUTINE
13066	2217174 LOOP2,	MEMA @ START	/LOAD FIRST POINTED TO WORD (20 BITS)
13067	1421172	APMM CHKSUM	/ADD IT TO RUNNING CHECKSUM
13070	2013127	JMS VAXOUT	/CALL RS232B XMIT WORD ROUTINE
13071	1437174	MPOM START	/INCREMENT WORD POINTER
13072	1535175	MMOMZ SIZE	/DECREMENT SIZE POINTER
13073	13066	JMP LOOP2	/JMP IF NOT DONE
13074	1013065	JMP @ DUMP	/NORMAL RETURN
13075	0 STATUS,	0	/READ SIZE AND POSITION OF CURRENT
			/DISPLAY REGION FROM BUTTON REGISTER
13076	1217170	MEMA STADDR	/LOAD BASE ADDRESS
13077	1413174	ACCM START	/AND PUT IT IN START
13100	603640	RDBUTN	/READ BUTTONS
13101	1263166	ANDIA MASK5	/SELECT START ADDRESS BITS
13102	250000	AMOA	/SEE IF ADDRESS IS ZERO
13103	53106	JFZ SETADR	/START ADDRESS IS ZERO
13104	246000	AFOA	/ADD 1 BACK
13105	701007	LASH 7	/CONVERT TO OFFSET
13106	1421174 SETADR,	APMM START	/ADD TO BASE ADDRESS
13107	603640	RDBUTN	/READ BUTTONS AGAIN
13110	1263167	ANDIA MASK6	/SELECT SIZE BITS
13111	716001	RISH 1	/CONVERT TO ABSOLUTE SIZE
13112	1413175	ACCM SIZE	/ AND STORE IT IN SIZE
13113	1013075	JMP @ STATUS	/NORMAL RETURN
13114	0 INTENA,	0	/ENABLE INT5- FOR RS232B ROUTINE
13115	1436545	MPOM DIFLAG	/DISABLE "O,"X
13116	216004	MEMA (4	/LEVEL 5 DEVICE ENABLE
13117	1460012	AOMM SPDMSK	/LOAD COPY OF SLOW I/O MASK
13120	500761	LDMSK	/LOAD SLOW I/O MASK
13121	12	SPDMSK	/ MASK COPY
13122	216040	MEMA (40	/LVL5 ENABLE
13123	1460011	AOMM ZILMSK	/LOAD INT ENABLE MASK COPY
13124	503771	LMASK	/LOAD INT ENABLE MASK
13125	11	ZILMSK	/ MASK COPY
13126	1013114	JMP @ INTENA	/NORMAL RETURN
13127	0 VAXOUT,	0	/RS232B XMIT 20 BIT WORD ROUTINE
13130	1413173	ACCM TEMP	/STORE WORD IN TEMP
13131	2013152	JMS RSXMIT	/OUTPUT LOWEST 5 BITS
13132	1217173	MEMA TEMP	/RESTORE FULL WORD
13133	716005	RISH 5	/SHIFT 5
13134	2013152	JMS RSXMIT	/XMIT NEXT 5 BITS
13135	1217173	MEMA TEMP	/RESTORE FULL WORD
13136	716012	RISH 12	/SHIFT 10

/USER OVERLAY TO NTCFT V#90615

```

13137 2013152      JMS RSXMIT      /XMIT NEXT 5 BITS
13140 1217173      MEMA TEMP      /RESTORE FULL WORD
13141 716017       RISH 17        /SHIFT 15
13142 2013152      JMS RSXMIT      /XMIT NEXT 5 BITS
13143 216015       MEMA (15)      /LOAD A VAX <CR>
13144 2117163      MEMZ @ XMTFLG  /SKIP IF VAX BUFFER NOT FULL
13145 13144        JMP #-1
13146 644730       RSOUTF         /RS232B OUT READY?
13147 13146        JMP #-1        /NO
13150 620731       RSOUT          /OUTPUT THE <CR>
13151 1013127      JMP @ VAXOUT
13152 0 RSXMIT, 0   /ROUTINE TO OUTPUT 1 BYTE
13153 262037       ANDA (37)       /AND AC WITH 00011111
13154 260040       ADMA (40)       /SET BITS 001XXXXX
13155 2117163      MEMZ @ XMTFLG  /SKIP IF VAX11 BUFFER NOT FULL
13156 13155        JMP #-1        /NO GO
13157 644730       RSOUTF         /RS232B OUT READY?
13160 13157        JMP #-1        /NO
13161 620731       RSOUT          /OUTPUT 1 BYTE
13162 1013152      JMP @ RSXMIT    /DONE
13163 4761 XMTFLG, CTRLCH
13164 3777773 MASK3, 3777773      /MASK FOR DISABLING SLOW I/O
13165 3777737 MASK4, 3777737      /MASK FOR DISABLING INTS-
13166 777 MASK5, 777              /BUTTON REGISTER MASK
13167 777000 MASK6, 777000        /BUTTON REGISTER MASK
13170 20000 STADDR, 20000         /BASE ADDRESS FOR STATUS ROUTINE
13171 564 GFLAG, GFLAG           /ACQ IN PROGRESS FLAG LOC EQUIV
13172 0 CHKSUM, 0                /CHECKSUM
13173 0 TEMP, 0                  /TEMP STORAGE
13174 0 START, 0                 /LOCATION START POINTER
13175 0 SIZE, 0                  /BLOCK SIZE
13176 222362 MSG1, .TEXT "RS232B DATA XMIT TO VAX (NO LEADER)C"
13177 636202
13200 400401
13201 240140
13202 301511
13203 244024
13204 174026
13205 13040
13206 501617
13207 401405
13210 10405
13211 225133
13212 370000
13213 222362 MSG2, .TEXT "RS232B DATA TRANSMIT TO VAX (LEADER)C"
13214 636202
13215 400401
13216 240140
13217 242201
13220 162315
13221 112440
13222 241740
13223 260130
13224 405014

```

/USER OVERLAY TO NTCFT V#90615

```

13225  50104
13226  52251
13227  333700
16341      *16341
16341  13000 POINT1, VAXDMP
16343      *16343
16343  13032 POINT2, LEADER

```

-5-

/USER OVERLAY TO NTCFT V#90615

BRETRN	10152	CHKSUM	13172	CTRLCH	4761	DIFLAG	545
DSFSIZ	466	DSFSTA	467	DUMP	13065	GADCHK	10061
GFLAG	13171	GOFLAG	564	INTENA	13114	LEADER	13032
LOOP2	13066	MASK3	13164	MASK4	13165	MASK5	13166
MASK6	13167	MSG1	13176	MSG2	13213	POINT1	16341
POINT2	16343	RSBINT	4751	RSXMIT	13152	SETADR	13106
SIZE	13175	SFDMSK	12	STADDR	13170	START	13174
STATUS	13075	TEMP	13173	VAXDMF	13000	VAXOUT	13127
XMTFLG	13163	ZDISPL	2301	ZILMSK	11	ZTTYBU	2636
ZUPKBU	2533	ZWFCHA	2627				

*

Figure A1.2 VAX-11 FORTRAN IV-PLUS Listing of the Data Transmission
 Program FILECOPY.FOR

5-Sep-1979 16:52:07

VAX-11 FORTRAN IV-PLUS V1.2-13

Page

1

FILECOPY.FOR.12

```

C RS232B ENCODED DATA RECEIVE AND PACK ROUTINE
C
C THIS ROUTINE RECEIVES ENCODED DATA FROM THE NICOLET 1180
C DATA SYSTEM VIA A TERMINAL MULTIPLEXER INPUT. 1180 DATA
C WILL BE IN THE FORM OF VECTORS OR ARRAYS, WITH A MAXIMUM
C SIZE OF 16K X 128 20 BIT WORDS. THE 20 BIT WORDS ARE
C TRANSMITTED 5 BITS/BYTE IN THE FORMAT 001XXXXX WHERE XXXXX
C ARE THE SIGNIFICANT BITS. EACH WORD IS TRANSMITTED AS
C A STRING OF 4 BYTES FOLLOWED BY A <CR> CHARACTER. THE 001
C HIGH ORDER BITS INSURE THAT NO CONTROL CHARACTERS WILL
C APPEAR AS PART OF THE DATA STREAM.
C
C VARIABLE DEFINITIONS:
C
C      A(*)          DATA POINT ARRAY, MAX SIZE 16K + 352 WORD
C                     PARAMETER TABLE
C      CHECKSUM      RUNNING CHECKSUM CALCULATED DURING DATA RECEIVE
C      TITLE         A60 HEADER BLOCK WRITTEN ON OUTPUT FILE
C      MASK1         00011111 MASK USED TO SELECT 5 LSB OF EACH BYTE
C      MASK2         11111111111111111111 MASK USED TO SELECT
C                     20 LSBs OF EACH 32 BIT WORD ON VAX
C      KK            SET TO OCTAL 21 FOR CTRL-Q TO INITIATE
C                     DATA TRANSMISSION
C      IP            NUMBER OF DATA POINTS
C      IJ(*)         VECTOR USED FOR DATA READ-IN, 4 WORDS LONG,
C                     ONE BYTE PER LOCATION
C      N             RECORD COUNT ON DATA READ FROM RS232
C      ISUM          RUNNING SUM USED IN PACKING DATA
C
C-----
0001      DIMENSION A(16736),IJ(4),TITLE(60)
0002      INTEGER A, CHECKSUM, TITLE, ANS
0003      LOGICAL CONTIN,COUNT,TYPE,LOG
C      OUTPUT WRITTEN ON FILE 'FOR009.DAT;*' WITH CARRIAGE CONTROL SUPPRESSED
C
0004      OPEN (UNIT=09, NAME='FOR009',TYPE='NEW',
0005             1CARRIAGECONTROL='LIST')
0005      CONTIN = .TRUE.
0006      COUNT = .FALSE.
0007      TYPE = .FALSE.
0008      LOG = .FALSE.
0009      MAXLIM = 100
C      SUPPRESS ERROR MESSAGES FROM INTEGER OVERFLOW CAUSED BY CHECKSUM
C      CALCULATION
0010      CALL ERRSET(70,CONTIN,COUNT,TYPE,LOG,MAXLIM)
0011      MASK1 = '37
0012      MASK2 = '3777777
0013      KK = '21
0014      WRITE (6,1100)
0015      1100      FORMAT ('ENTER NUMBER OF FILES, POINTS PER FILE')
0016      READ (5,1110) IF,IP
0017      1110      FORMAT (2I6)
C      OUTPUT FILE MAY BE WRITTEN WITH OR WITHOUT A 352 WORD NTCFT-1180
C      COMPATIBLE PARAMETER TABLE. ALL FILES IN A FILEGROUP MUST
C      HAVE THE SAME OPTION EXERCISED. ON RELOADING THE 1180 WITH
C      SUCH A PARAMETER TABLE, LOAD ADDRESS IS 17240. NTCFT THEN
C      EXAMINES THE LEADING TWO KEYWORDS TO INSURE THAT THE FILE IS

```

FILECOPY\$MAIN 5-Sep-1979 16:52:07 VAX-11 FORTRAN IV-PLUS V1.2-13 Page
2
FILECOPY.FOR.12

```

      C      AN NTCFT WRITTEN FILE.
0018      WRITE (6,1112)
0019      1112  FORMAT ('STORE WITH PARAMETER TABLES (Y OR N)?')
0020      50    READ (5,1114) ANS
0021      1114  FORMAT(1A1)
0022      IF (ANS.EQ.1HY) THEN
0023      IF = IP + 352
0024      GO TO 60
0025      ELSE IF (ANS.EQ.1HN) THEN
0026      GO TO 60
0027      ELSE
0028      WRITE (6,1116)
0029      1116  FORMAT ('PLEASE ANSWER Y OR N')
0030      GO TO 50
0031      END IF
      C      WRITE A TWO LINE HEADER BLOCK (DEFAULT STANDARD FORMAT)
0032      60    WRITE (6,1120)
0033      1120  FORMAT ('ENTER TWO LETTER ID, 5 DIGIT EXP #, AND 2 DIGIT
      1 BLOCK # --')
0034      READ (5,1130) IX1,IX2,IX3
0035      WRITE (9,1130) IX1,IX2,IX3
0036      1130  FORMAT (A2,I5,I2)
0037      WRITE (6,1140)
0038      1140  FORMAT (' EXPERIMENT IDENTIFICATION (A60) --')
0039      READ (5,1150) TITLE
0040      WRITE (9,1150) TITLE
0041      1150  FORMAT (60A1)
      C      NOW LOOP THROUGH ONCE FOR EACH FILE.  START BY TRANSMITTING A
      C      CTRL-Q TO INSURE THAT 1180 IS DUMPING DATA.  CTRL-Q SETS INT5-
      C      ON THE 1180 AND CAUSES THE INTERRUPT SERVICE TO SET A DATA
      C      TRANSMIT FLAG.
0042      DO 230, J=1,IF
0043      WRITE (2) KK
0044      CHECKSUM = 0
      C      FOR EACH FILE, LOOP THROUGH IF TIMES TO RECONSTRUCT IF WORDS FROM
      C      IF*4 BYTES.  A 20 BIT WORD IS RECEIVED IN THE FOLLOWING BYTE
      C      SEQUENCE:
      C
      C      AC4-0,AC9-5,AC14-10,AC19-15 <CR>
      C
      C      RECONSTRUCTION (PACKING) OCCURS BY APPROPRIATE LEFT-SHIFTING
      C      AND LOGICAL OR-ING OF THE FOUR BYTES STORED IN IJ(*)
0045      DO 200, K=1,IF
0046      100    READ (1,1000) N,IJ
0047      1000  FORMAT (Q,100A1)
0048      ISUM = 0
0049      DO 120 L=4,1,-1
0050      IJ(L) = IAND(IJ(L),MASK1)
0051      LSHIFT = (L-1)*5
0052      120    ISUM = IOR(ISUM,ISHFT(IJ(L),LSHIFT))
0053      CHECKSUM = CHECKSUM + ISUM
0054      200    A(K) = ISUM
      C      NOW RECEIVE AND RECONSTRUCT THE CHECKSUM CALCULATED BY THE 1180,
      C      AND COMPARE TO THE LOCALLY CALCULATED ONE.
0055      READ(1,1000) (IJ(I),I=1,4)
0056      ISUM = 0

```

FILECOPY\$MAIN 5-Sep-1979 16:52:07 VAX-11 FORTRAN IV-PLUS V1.2-13 Page
3

FILECOPY.FOR.12

```

0057      DO 220, L=4,1,-1
0058      IJ(L) = IAND(IJ(L),MASK1)
0059      LSHIFT = (L-1)*5
0060      ISUM = IOR(ISUM,ISHFT(IJ(L),LSHIFT))
0061      WRITE (9,2000) (A(KL),KL=1,IP)
0062      CHECKSUM = IAND(CHECKSUM,MASK2)
0063      ISUM = IAND(ISUM,MASK2)
0064      230  WRITE(6,1010) J,CHECKSUM,ISUM
0065      1010 FORMAT('OFILE =',I3,'OALCULATED CHECKSUM =',09/
0066      2000 1'O RECEIVED CHECKSUM =',09)
0067      FORMAT (BI10)
      END

```

FILECOPY\$MAIN 5-Sep-1979 16:52:07 VAX-11 FORTRAN IV-PLUS V1.2-13 Page

4

FILECOPY.FOR.12

PROGRAM SECTIONS

Name	Bytes	Attributes
0 \$CODE	835	PIC CON REL LCL SHR EXE RD NOWRT LONG
1 \$PDATA	313	PIC CON REL LCL SHR NOEXE RD NOWRT LONG
2 \$LOCAL	67396	PIC CON REL LCL NOSHR NOEXE RD WRT LONG

ENTRY POINTS

Address	Type	Name
0-00000000		FILECOPY\$MAIN

VARIABLES

Address	Type	Name	Address	Type	Name	Address	Type	Name
2-00010684	I*4	ANS	2-00010680	I*4	CHECKSUM	2-00010688		
L*4 CONTIN			L*4 COUNT					
2-000106D4	I*4	I	2-000106A8	I*4	IF	2-000106AC		
I*4 IP			I*4 ISUM					
2-000106E0	I*4	IX1	2-000106B4	I*4	IX2	2-000106B8		
I*4 IX3			I*4 J					
2-000106C0	I*4	K	2-000106A4	I*4	KK	2-000106D8		
I*4 KL			I*4 L					
2-00010694	L*4	LOG	2-000106D0	I*4	LSHIFT	2-0001069C		
I*4 MASK1			I*4 MASK2					
2-00010698	I*4	MAXLIM	2-000106C4	I*4	N	2-00010690		
L*4 TYPE								

ARRAYS

Address	Type	Name	Bytes	Dimensions
2-00000000	I*4	A	66944	(16736)
2-00010580	I*4	IJ	16	(4)
2-00010590	I*4	TITLE	240	(60)

LABELS

Address	Label	Address	Label	Address	Label	Address	Label
0-00000093	50	0-000000EB	60	**	100	**	
120	**	200	**	220	**		
**	230	1-000000EE	1000'	1-000000F4	1010'	1-0000000B	
1100'	1-00000035	1110'	1-0000003A	1112'			
1-00000063	1114'	1-00000066	1116'	1-0000007E	1120'	1-000000BC	
1130'	1-000000C3	1140'	1-000000E9	1150'			
1-00000134	2000'						

FUNCTIONS AND SUBROUTINES REFERENCED

ERRSET FOR\$OPEN MTH\$JSHIFT

Figure A1.3 VAX-11 FORTRAN-IV PLUS Listing of the Data Transmission
 Program ARRAY.FOR

3-SEP-1979 13:55:55

VAX-11 FORTRAN IV-PLUS V1.2-13

Page

ARRAY.FOR.23

```

C RS232B ENCODED DATA RECEIVE AND PACK ROUTINE
C
C THIS ROUTINE RECEIVES ENCODED DATA FROM THE NICOLET 1180
C DATA SYSTEM VIA A TERMINAL MULTIPLEXER INPUT. 1180 DATA
C WILL BE IN THE FORM OF VECTORS OR ARRAYS, WITH A MAXIMUM
C SIZE OF 8K X 128 20 BIT WORDS. THE 20 BIT WORDS ARE
C TRANSMITTED 5 BITS/BYTE IN THE FORMAT 001XXXXX WHERE XXXXX
C ARE THE SIGNIFICANT BITS. EACH WORD IS TRANSMITTED AS
C A STRING OF 4 BYTES FOLLOWED BY A <CR> CHARACTER. THE 110
C HIGH ORDER BITS INSURE THAT NO CONTROL CHARACTERS WILL
C APPEAR AS PART OF THE DATA STREAM.
C
C VARIABLE DEFINITIONS:
C
C      A(*,*)      INTEGER ARRAY FOR UP TO 128X128 20 BIT WORDS
C      IJ(*)        4 WORD VECTOR FOR READING 4 UNPACKED BYTES
C      CHECKSUM     RUNNING CHECKSUM CALCULATED ON RECEIVED DATA
C      TITLE(*)     A40 HEADER BLOCK FOR LABELING FILES
C      MASK1        00011111 MASK FOR PICKING OUT SIGNIFICANT BITS
C      MASK2        11111111111111111111 MASK TO PICK OUT 20 BITS
C      KK           SET TO OCTAL 21 (CTRL-Q), USED TO INITIATE
C                  DATA TRANSMISSION
C      IF           NUMBER OF FILES TO BE READ, USUALLY 128
C      IP           NUMBER OF POINTS PER FILE, USUALLY 128
C      N            SET TO CHARACTER COUNT ON EXECUTION OF READ
C                  STATEMENT 100
C      LSHIFT       LEFT SHIFT COUNT
C
C -----
0001      DIMENSION A(128,128),IJ(4)
0002      INTEGER A, CHECKSUM, TITLE(60)
0003      LOGICAL CONTIN,COUNT,TYPE,LOG
0004      DATA A/16384*0/
C  OUTPUT FILE 'FOR010.DAT' WILL CONTAIN THE DATA ARRAY
0005      OPEN (UNIT=10, NAME='FOR010',TYPE='NEW',
             1CARRIAGECONTROL='LIST')
0006      CONTIN = .TRUE.
0007      COUNT = .FALSE.
0008      TYPE = .FALSE.
0009      LOG = .FALSE.
0010      MAXLIM = 100
C  CALL ERRSET TO OVERRIDE INTEGER OVERFLOWS DURING CALCULATION
C  OF CHECKSUMS
0011      CALL ERRSET(70,CONTIN,COUNT,TYPE,LOG,MAXLIM)
0012      MASK1 = '37
0013      MASK2 = '3777777
0014      KK = '21
0015      WRITE (6,1100)
0016      1100  FORMAT ('ENTER NUMBER OF FILES, POINTS PER FILE')
0017      READ (5,1110) IF,IP
0018      1110  FORMAT (2I4)
C  READ AND WRITE OUT A HEADER BLOCK THAT IS COMPATIBLE WITH THE
C  VAX 11 FORTRAN IV PLUS CONTOUR PLOTTING ROUTINE
C  IKAR120C3NICOLET
0019      20  WRITE (6,1120)

```

ARRAY\$MAIN 3-Sep-1979 13:55:55 VAX-11 FORTRAN IV-PLUS V1.2-13 Page
2

ARRAY.FOR.23

```

0020 1120 FORMAT ('CENTER TWO LETTER ID, 5 DIGIT EXP #, AND 2 DIGIT
        1 BLOCK # --')
0021      READ (5,1130) IX1,IX2,IX3
0022      WRITE (10,1130) IX1,IX2,IX3
0023 1130 FORMAT (A2,I5,I2)
0024      WRITE (6,1140)
0025 1140 FORMAT (' EXPERIMENT IDENTIFICATION (A60) --')
0026      READ (5,1150) TITLE
0027      WRITE (10,1150) TITLE
0028 1150 FORMAT (60A1)
        C INITIATE 1180 DATA TRANSMISSION BY SENDING A CTRL-Q OUT
        C TERMINAL MULTIPLEXER PORT THAT IS ALLOCATED TO 1180
        C RS232B CHANNEL
0029      WRITE (2) KK
0030      DO 230, J=1,IF
0031      CHECKSUM = 0
0032      DO 200, K=1,IF
        C A 20 BIT WORD IS RECEIVED IN THE FOLLOWING SEQUENCE:
        C
        C AC4-0,AC9-5,AC14-10,AC19-14 <CR>
        C
        C RECONSTRUCTION OCCURS BY APPROPRIATE MASKING, LEFT SHIFTING,
        C AND LOGICAL OR-ING OF THE FOUR BYTES STORED IN IJ(*)
0033 100 READ (1,1000) N,IJ
0034 1000 FORMAT (8,100A1)
0035      ISUM = 0
0036      DO 120 L=4,1,-1
0037      IJ(L) = IAND(IJ(L),MASK1)
0038      LSHIFT = (L-1)*5
0039 120 ISUM = IOR(ISUM,ISHFT(IJ(L),LSHIFT))
0040      CHECKSUM = CHECKSUM + ISUM
0041 200 A(J,K) = ISUM
        C NOW RECEIVE THE 1180 CALCULATED CHECKSUM AND COMPARE IT TO
        C THE LOCALLY CALCULATED ONE
0042      READ(1,1000) (IJ(I),I=1,4)
0043      ISUM = 0
0044      DO 220, L=4,1,-1
0045      IJ(L) = IAND(IJ(L),MASK1)
0046      LSHIFT = (L-1)*5
0047 220 ISUM = IOR(ISUM,ISHFT(IJ(L),LSHIFT))
0048      CHECKSUM = IAND(CHECKSUM,MASK2)
0049 230 WRITE(6,1010) J,CHECKSUM,ISUM
0050 1010 FORMAT('OFILE =',I3,'OCALCULATED CHECKSUM =',O9/
        1'ORECEIVED CHECKSUM =',O9)
0051      WRITE (10,2000) ((A(M,N),N=1,128),M=1,128)
0052 2000 FORMAT (8I10)
        C EKAR120CINICOLET LOOKS FOR TAGWORDS AS DELIMITERS OF SERIALY
        C STORED DATA BLOCKS. AS MANY 128 X 128 BLOCKS AS
        C DESIRED MAY BE STORED SEQUENTIALLY, OR THEY MAY BE
        C LATER COMBINED USING THE VMS COMMAND <APPEND>.
        C
0053      WRITE (6,2100)
0054 2100 FORMAT ('OAPPEND ANOTHER ARRAY (Y OR N)?')
0055 300 READ (5,2110) IANS
0056 2110 FORMAT (A1)
0057      IF (IANS.EQ.1NY) THEN

```

```
ARRAY$MAIN      3-Sep-1979 13:55:55   VAX-11 FORTRAN IV-PLUS V1.2-13   Page
  3
                                ARRAY.FOR.23

0058             GO TO 20
0059             ELSE IF (IANS.EQ.1HN) THEN
0060             GO TO 400
0061             ELSE
0062             WRITE (6,2120)
0063 2120         FORMAT ('ANSWER Y OR N')
0064             GO TO 300
0065             END IF
0066 400          END
```

ARRAY\$MAIN 3-Sep-1979 13:55:55 VAX-11 FORTRAN IV-PLUS V1.2-13 Page
4
ARRAY.FOR.23

PROGRAM SECTIONS

Name	Bytes	Attributes
0 \$CODE	852	PIC CON REL LCL SHR EXE RD NOWRT LONG
1 \$PDATA	299	PIC CON REL LCL SHR NOEXE RD NOWRT LONG
2 \$LOCAL	65988	PIC CON REL LCL NOSHR NOEXE RD WRT LONG

ENTRY POINTS

Address	Type	Name
0-00000000		ARRAY\$MAIN

VARIABLES

Address	Type	Name	Address	Type	Name	Address	Type	Name
2-00010100	I*4	CHECKSUM	2-00010104	L*4	CONTIN	2-00010108		
L*4 COUNT		2-00010150	I*4 I					
2-00010158	I*4	IAN5	2-00010124	I*4	IF	2-00010128		
I*4 IP		2-00010144	I*4	ISUM				
2-0001012C	I*4	IX1	2-00010130	I*4	IX2	2-00010134		
I*4 IX3		2-00010138	I*4	J				
2-0001013C	I*4	K	2-00010120	I*4	KK	2-00010148		
I*4 L		2-00010110	L*4	LOG				
2-0001014C	I*4	LSHIFT	2-00010154	I*4	M	2-00010118		
I*4 MASK1		2-0001011C	I*4	MASK2				
2-00010114	I*4	MAXLIM	2-00010140	I*4	N	2-0001010C		
L*4 TYPE								

ARRAYS

Address	Type	Name	Bytes	Dimensions
2-00000000	I*4	A	65536	(128,128)
2-00010000	I*4	IJ	16	(4)
2-00010010	I*4	TITLE	240	(60)

LABELS

Address	Label	Address	Label	Address	Label	Address	Label
0-0000007D	20	**	100	**	120	**	
200	**	220	**	230			
0-000002FF	300	0-00000350	400	1-000000AA	1000'	1-00000080	
1010'	1-0000000B	1100'	1-00000035	1110'			
1-0000003A	1120'	1-00000078	1130'	1-0000007F	1140'	1-000000A5	
1150'	1-000000F0	2000'	1-000000F5	2100'			
1-00000117	2110'	1-0000011A	2120'				

FUNCTIONS AND SUBROUTINES REFERENCED

ERRSET FOR\$OPEN MTH\$JISHFT

Figure A1.4 Assembled Listing of USR Overlay to NTCFT-1180
 for Data Transmission to NIC-80

RUN FASM

FASM, E58-90215
@HR220.ASC/HR220.BIN,-TT:F

FAST 1180 ASSEMBLER, E58-90215

08/23/79 10:03 @HR220.ASC/HR220.BIN,-TT:F

150 SOURCE STATEMENTS
28 USER-DEFINED SYMBOLS
MEMORY NEEDED: 13000 TO 16343

/USR OVERLAY TO NTCFT V#90615

/USR OVERLAY TO NTCFT V#90615
 .UNDENT

```

/      LAST UPDATE: 27 JUNE 1979
/
/      For transmission of data and Programs to the NIC-80 system
/      attached to the HR-220 spectrometer system
/
/      U1 => read Programs/data from NIC-80
/      U2 => send Programs/data to NIC-80
/
/      293A SSW      UP => Program, occupying 0000-7577
/                   down => Data, occupying 100000-117777 max
/                   (button selected)

```

```

600610      TSPW=600610
13000      *13000
13000      0 STATUS, 0      /READ BUTTONS FOR MEMORY START
13001 1217172      MEMA STADDR /LOAD BASE ADDRESS
13002 1413173      ACCM START  /AND STORE IN START
13003 603640      RDBUTN      /READ BUTTONS
13004 1263170      ANDA MASK1  /SELECT START ADDRESS BITS
13005 250000      AMOA        /SEE IF START ADDRESS IS ZERO
13006 53011      JFZ SETADR    /START ADDRESS IS 0
13007 246000      APOA        /ADD 1 BACK
13010 701007      LASH 7      /CONVERT TO OFFSET
13011 1221173 SETADR, AFMA START /ADD TO BASE ADDRESS
13012 1413173      ACCM START  /STORE FINAL START ADDRESS
13013 603640      RDBUTN      /READ BUTTONS AGAIN
13014 1263171      ANDA MASK2  /SELECT SIZE BITS
13015 716001      RISH 1      /CONVERT TO ABSOLUTE SIZE
13016 1413174      ACCM SIZE   /STORE SIZE IN SIZE
13017 600610      TSPW        /LOOK AT SENSE SWITCHES
13020 1363027      ANDAZ MASK0 /PICK OUT SW 3 (BIT 15)
13021 2013023      JMS SETSIZ  /SET SIZE = 7600
13022 1013000      JMP @ STATUS
13023      0 SETSIZ, 0      /SET SIZE TO 7600
13024 1217030      MEMA FSIZE
13025 1413174      ACCM SIZE
13026 1013023      JMP @ SETSIZ
13027 100000 MASK0, 100000      /MASK TO PICK OUT SSW 3
13030 7600 FSIZE, 7600
13031      0 RSREAD, 0      /START OF RS232B READ
13032 1403175      ZERM CHKSUM /ZERO CHECKSUM
13033 2013000      JMS STATUS  /FIND ADDRESSES TO USE
13034 2013114 LOOP1, JMS XON   /TRANSMIT X-ON
13035 644720      RSINF       /RS232 READY TO LOOK FOR X-ON?
13036 13035      JMP #-1      /NOT READY
13037 600721      RSIN        /READ RS232
13040 322021      AMMAZ (021 /IS IT X-ON?
13041 13034      JMP LOOP1    /NO, KEEP LOOKING
13042 2013101 LOAD, JMS READ   /READ 8 BITS FROM RS232
13043 2413173      ACCM @ START /STORE LOWEST BITS
13044 2013114      JMS XON    /SEND CONFIRMING X-ON

```


/USR OVERLAY TO NTCFT V#90615

13045	2013101	JMS READ	/READ NEXT MSB'S
13046	701010	LASH 10	/SHIFT
13047	2421173	APMM @ START	/ADD ON AFTER SHIFT
13050	2013114	JMS XON	/SEND CONFIRMING X-ON
13051	2013101	JMS READ	/READ MSB'S
13052	701020	LASH 20	/SHIFT L6
13053	2421173	APMM @ START	/ADD ON AFTER SHIFT
13054	2013114	JMS XON	/SEND CONFIRMING X-ON
13055	2217173	MEMA @ START	/LOAD WHOLE WORD AND
13056	1421175	APMM CHKSUM	/ ADD TO CHECKSUM
13057	1437173	MPOM START	/ADVANCE POINTER
13060	1535174	MMOMZ SIZE	/DECREMENT SIZE, SKIP ON 0
13061	13042	JMP LOAD	/NEXT POINT
13062	2013101	JMS READ	/NOW RECEIVE CHECKSUM
13063	1413176	ACCM CSUM	/STORE 8 LSB'S
13064	2013114	JMS XON	/SEND
13065	2013101	JMS READ	/NEXT 8 MSB'S
13066	701010	LASH 10	/SHIFT
13067	1421176	APMM CSUM	/STORE SHIFTED
13070	2013114	JMS XON	/SEND CONFIRMING X-ON
13071	2013101	JMS READ	/MSB'S
13072	701020	LASH 20	/SHIFT
13073	1421176	APMM CSUM	/ADD ON
13074	2013114	JMS XON	/SEND CONFIRMING X-ON
13075	1217176	MEMA CSUM	/NOW COMPARE WITH COUNTED
13076	1323175	AMMAZ CHKSUM	/ CHECKSUM
13077	2013106	JMS ERROR	
13100	1013031	JMP @ RSREAD	/RETURN TO NTCFT
13101	0 READ,	0	/RS232 READ ROUTINE
13102	644720	RSINF	/RS232B READY?
13103	13102	JMP #-1	
13104	600721	RSIN	/READ 8 BITS INTO ACC
13105	1013101	JMP @ READ	
13106	0 ERROR,	0	/CHECKSUM ERROR ROUTINE
13107	216207	MEMA (207	/RING BELL
13110	644710	TTYPF	
13111	13110	JMP #-1	
13112	620711	PRTTY	
13113	1013106	JMP @ ERROR	
13114	0 XON,	0	/TRANSMIT X-ON ROUTINE
13115	644730	RSOUTF	/RS232B READY TO XMIT?
13116	13115	JMP #-1	/NO
13117	216021	MEMA (021	/X-ON
13120	620731	RSOUT	/TRANSMIT X-ON
13121	1013114	JMP @ XON	
13122	0 RSDUMP,	0	/START OF RS232B TRANSMIT
13123	1403175	ZERM CHKSUM	/ZERO CHECKSUM
13124	2013000	JMS STATUS	/FIND ADDRESSES TO USE
13125	2013160	JMS RXON	/WAIT UNTIL X-ON RECEIVED
13126	2013114	JMS XON	/SEND ACKNOWLEDGING X-ON
13127	2217173	MEMA @ START	/LOAD POINTED TO WORD
13130	1421175	APMM CHKSUM	/ADD TO CHECKSUM
13131	2013140	JMS WORFUN	/TRANSMIT WORD
13132	1437173	MPOM START	/INCREMENT POINTER

/USR OVERLAY TO NTCFT V#90615

```

13133 1535174      MMOMZ SIZE      /DECREMENT SIZE COUNT
13134   13127      JMP LOOP2        /NEXT WORD
13135 1217175      MEMA CHKSUM      /NOW TRANSMIT CHECKSUM
13136 2013140      JMS WOPUN
13137 1013122      JMP @ RSDUMP
13140           0 WOPUN, 0          /20 BIT WORD XMIT ROUTINE
13141 1413167      ACCM TEMP        /STORE WORD IN TEMP
13142 2013152      JMS OUT          /DUMP LOWEST 8 BITS
13143 1217167      MEMA TEMP        /RESTORE WORD AND
13144   716010      RISH 10         / SHIFT OVER 8
13145 2013152      JMS OUT          /DUMP NEXT 8
13146 1217167      MEMA TEMP        /RESTORE WORD AGAIN AND
13147   716020      RISH 20         / SHIFT OVER 16
13150 2013152      JMS OUT          /DUMP 4 MSB'S + JUNK
13151 1013140      JMP @ WOPUN
13152           0 OUT, 0          /RS232B OUTPUT ROUTINE
13153 644730       RSOUTF
13154   13153      JMP #-1
13155 620731       RSOUT
13156 2013160      JMS RXON         /WAIT FOR REC X-ON
13157 1013152      JMP @ OUT
13160           0 RXON, 0         /RECEIVE X-ON ROUTINE
13161 644720       RSINF           /WAIT FOR RS232B INPUT READY
13162   13161      JMP #-1
13163 600721       RSIN
13164 322021       AMMAZ (021
13165   13161      JMP #-4
13166 1013160      JMP @ RXON
13167           0 TEMP, 0        /TEMP STORAGE
13170           777 MASK1, 777
13171 777000 MASK2, 777000
13172 20000 STADDR, 20000
13173           0 START, 0
13174           0 SIZE, 0
13175           0 CHKSUM, 0
13176           0 CSUM, 0
13177           0 COUNT1, 0      /TEMP COUNTER
16341           *16341
16341 13031 POINT1, RSREAD
16343           *16343
16343 13122 POINT2, RSDUMP

```

-3-

/USR OVERLAY TO NTCFT V#90615

CHKSUM	13175	COUNT1	13177	CSUM	13176	ERROR	13106
LOAD	13042	LOOP1	13034	LOOP2	13127	MASK0	13027
MASK1	13170	MASK2	13171	OUT	13152	POINT1	16341
POINT2	16343	PSIZE	13030	READ	13101	RSDUMP	13122
RSREAD	13031	RXON	13160	SETADR	13011	SETSIZ	13023
SIZE	13174	STADDR	13172	START	13173	STATUS	13000
TEMP	13167	TSPW	600610	WOPUN	13140	XON	13114

APPENDIX 2 Complex Impedance Calculation

Figure A2.1 BASIC Program Listing for Calculating the Complex
 Impedance from the Magnitude and Phase of the
 Reflection Coefficient

```

100 REM  CALCULATE COMPLEX IMPEDANCE FROM REFLECTION COEFFICIENT
110 PRINT "Enter A amplitude --";
120 INPUT A
130 PRINT "Enter B amplitude --";
140 INPUT B
150 PRINT "Enter relative phase --";
160 INPUT P
170 LET B=1.0532*B
180 PRINT
190 PRINT "B CORRECTION = 1.0532"
200 PRINT "CORRECTED B =" ; B
210 LET G=B/A
220 PRINT "GAMMA =" ; G ; " (Reflection Coefficient)"
230 LET P1=P*3.14159265/180.
240 LET X1=G*COS(P1)
250 LET Y1=G*SIN(P1)
260 LET N1=(X1+1.)*(X1+1.)+(Y1*Y1)
270 LET A1=SQR(N1)
280 LET A2=ATN(Y1/(X1+1.))*180./3.141592654
290 LET N2=(1.-X1)*(1.-X1)+(Y1*Y1)
300 LET A3=SQR(N2)
310 LET A4=ATN(-Y1/(1.-X1))*180./3.141592654
320 LET N3=51.1*A1/A3
330 LET N4=(A2-A4)*3.141592654/180.
340 LET R1=N3*COS(N4)
350 LET R2=N3*SIN(N4)
360 PRINT "REAL PART OF IMPEDANCE --";R1;" OHMS"
370 PRINT "IMAGINARY PART OF IMPEDANCE --";R2;" J OHMS"
380 PRINT
390 PRINT "MAGNITUDE =" ; N3
400 PRINT "PHASE ANGLE =" ; (A2-A4) ; " DEG"
410 PRINT
420 PRINT
430 PRINT "TO QUIT, TYPE 0, OTHERWISE Enter A amplitude --";
440 INPUT A
450 IF A=0. THEN      470
460 GO TO      130
470 END

```

APPENDIX 3 Supplementary Schematic Documentation

Figure A3.1 10X2063-S1 8 - 270 MHz SSB Generator


Figure A3.2 10X2093-S1 Four Phase NMR Generator


Figure A3.3 10X2103-S1 8 - 270 MHz to 30 MHz Linear Converter

Figure A3.4 10X2113-S1 30 MHz Double Phase Sensitive Detector


Figure A3.5 10X2122-S1 30 MHz Output Monitor Module

Figure A3.6 10X2153-S1 30 MHz Reference Amplifier






		FOUR PHASE NMR GENERATOR	
		SCHEMATIC	
		DATE: 2-1-68	
		BY: C.C. G.	
		UNIVERSITY OF CALIFORNIA	
		10X 20935-1	







NOTES: UNLESS OTHERWISE SPECIFIED

1. CAPACITORS 0.01uF 100V CERAMIC CK05BX103M.
- *2. FEED-THRU CAPACITORS ERIE #2425-001-W0502AA (B-32THD MTD).
3. RESISTORS 1/4W 1% METAL FILM.
- **4. VARIABLE CAPACITORS 9-35pF CERAMIC ERIE MODEL #538-000.

				30MH± OUTPUT MONITOR MODULE			

"To do is to be."

- J.P. Sartre

"To be is to do."

- E. Kant

"Do be do be do."

- F. Sinatra